

MAR 22 2001

Section 7- 510(k) Summary of Safety and Effectiveness

**7.1
Statement** This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and CFR 807.92

**7.2
Submitter** Endius, Inc.
23 West Bacon Street
Plainville, MA. 02762

**7.3
Company
Contact** Susan Finneran
Director Regulatory Affairs
508-643-0983

**7.4
Device Name** **Proprietary Name:** Endius Bipolar Sheath
Common Name: Bipolar Coagulation Device
Classification Name: Electrosurgical cutting and coagulation device and accessories

**7.5
Predicate
Legally
Marketed
Devices** The Bipolar Sheath is substantially equivalent to the Everest Bipolar Probes, Everest Medical (Minneapolis, MN).

**7.6
Device
Description** The Endius Bipolar Sheath is a stainless Steel tube covered with an insulation material that is intended to fit over an automated tissue removal blade. The device is intended to be connected to the Valley Lab's Force 2 generator by using the Bipolar Sheath Adapter which is intended to decrease the maximum voltage of the Valley Lab's Generator from 800 volts to 100 volts. This will ensure that the appropriate level of energy is transmitted to the Bipolar Sheath for the maximum performance.

7.7

Device**Indications and****Intended use**

The Endius Bipolar Sheath is intended to be used in conjunction with the Endius XPS Microdebrider System to cut and coagulate soft tissue during various spinal surgical procedures.

7.8

Substantial**Equivalence**

The Endius Bipolar Sheath is substantially equivalent to the Everest Bipolar Forceps


Table of Substantial Equivalence

Device Name	Everest Bipolar Forcep	Endius Bipolar Sheath and MDS Device
510k #	K945975	The subject of this submission
Intended use	The Everest Bipolar Forceps are intended to be used to cut and coagulate soft tissue electrosurgically.	The Endius Bipolar Sheath is intended to be used in conjunction with the Endius XPS Microdebrider System to cut and coagulate soft tissue during various spinal surgical procedures.
Materials	Stainless Steel, Plastics, Adhesives, and ceramics	Stainless Steel, Parylene coating, plastics, adhesives
Sterilization/ Labeling	Single Use, Sterilized by Ethylene Oxide	Single use, Sterilized by gamma irradiation
Sizes	3mm and 5mm	4.5mm

Applicant



Date





MAR 22 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Susan Finneran
Director Regulatory Affairs/Clinical Sciences
Endius, Incorporated
23 West Bacon Street
Plainville, Massachusetts 02762

Re: K003897
Trade Name: Endius Bipolar Sheath and Accessories
Regulatory Class: II
Product Code: GEI
Dated: February 6, 2001
Received: February 7, 2001

Dear Ms. Finneran:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

Miriam C. Provost
for Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K003897

Device Name: Endius Bipolar Sheath and Accessories

Indications for Use: The Endius Bipolar Sheath is intended to be used to coagulate soft tissue during various spinal surgical procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Optional Format 3-10-98)

(Posted July 1, 1998)

Miriam C. Provost
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K003897

Prescription Use ✓
(Per 21 CFR 801.109)