

JUL 11 2002

## 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

### Acetaminophen method for ADVIA® IMS™

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of the Safe Medical Device Act of 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K020792

#### 1. Intended Use

This *in-vitro* diagnostic method is intended to measure acetaminophen in human serum and plasma using lithium heparin as the anticoagulant on ADVIA® IMS™. Acetaminophen (Tylenol, paracetamol, p-hydroxyacetanilide) is used in many formulations as an analgesic, generally with no adverse effects. Measurement of acetaminophen is used in the diagnosis and treatment of severe liver damage caused by overdose through chronic usage, accident or self-infliction.

#### 2. Predicate Device

Product Name	Reagent Part #	Calibrator Part #
Abbott/TDx	953669	953602

#### 3. Device / Method

Product Name	Reagent BAN	Calibrator BAN
Bayer ADVIA IMS	01352146	04919015

#### Imprecision (Serum)

ADVIA IMS		Abbott/TDx	
Level (mg/dL)	Total CV(%)	Level (mg/dL)	Total CV(%)
1.5	2.6	1.5	4.9
5.0	1.5	3.5	3.0
14.9	1.3	15	3.9

#### Correlation (Y=ADVIA IMS, X=comparison system)

Specimen type	Comparison System (X)	N	Regression Equation	Syx (mg/dL)	R	Sample Range (mg/dL)
Serum	Abbott/TDx	46	$Y=1.05X - 0.25$	0.31	0.999	0.3-19.3
Plasma(y), Serum(x)	ADVIA IMS	50	$Y=1.00X + 0.05$	0.32	0.998	2.3-17.3

#### Interfering Substances


Interfering Substance	Interfering Sub. Conc. (mg/dL)	Acetaminophen Conc (mg/dL)	Effect (% change)
Bilirubin (unconjugated)	25	5.7	+9
Bilirubin (conjugated)	25	5.7	+4
Hemoglobin	1000	5.5	-4
Lipids (Triglycerides)	500	5.6	-9

#### Analytical Range

Serum/Plasma: 0 to 20 mg/dL

#### 4. Conclusion

Performance of the ADVIA IMS Acetaminophen Assay on a *Bayer ADVIA*® IMS™ is equivalent to the performance of the Acetaminophen Assay on the predicate device (Abbott TDX, K840941) and is within proposed manufacturing specifications. No safety and effectiveness issues have been raised.



Kenneth T. Edds  
Manager Regulatory Affairs  
Bayer Corporation  
511 Benedict Avenue  
Tarrytown, New York 10591-5097

Date

2/28/02



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road.  
Rockville MD 20850

JUL 11 2002

Kenneth T. Edds, Ph.D.  
Regulatory Affairs Manager  
Bayer Corporation  
511 Benedict Avenue  
Tarrytown, NY 10591-5097

Re: k020792  
Trade/Device Name: Acetaminophen Assay for the ADVIA® IMS™  
Regulation Number: 21 CFR 862.3030  
Regulation Name: Acetaminophen test system  
Regulatory Class: Class II  
Product Code: LDP  
Dated: June 19, 2002  
Received: June 20, 2002

Dear Dr. Edds:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 -

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical Laboratory-Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number: K020792

Device Name: Acetaminophen Assay for the ADVIA® IMS™

**Indications for Use:**

The *Bayer ADVIA IMS* Acetaminophen method is an *in vitro* diagnostic device intended to measure acetaminophen levels in human serum or plasma (Lithium heparin). Such measurements are used in the diagnosis of acetaminophen toxicity and overdose.

Jean Cooper  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K020792

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)