

APR 15 2002

## **Summary of Safety and Effectiveness Lyphochek® Hemostasis Control**

### **1.0 Submitter**

Bio-Rad Laboratories  
9500 Jeronimo Road,  
Irvine, California 92618-2017  
Telephone: (949) 598-1200  
Fax: (949) 598-1555

### **Contact Person**

Maria Zeballos  
Regulatory Affairs Specialist  
Telephone: (949) 598-1367

### **Date of Summary Preparation**

March 15, 2002

### **2.0 Device Identification**

Product Trade Name: Lyphochek® Hemostasis Control  
  
Common Name: Control Plasma Normal  
Control Plasma Abnormal  
  
Classifications: Class II  
Product Code: GGC  
Regulation Number: 21 CFR 864.5425

### **3.0 Device to Which Substantial Equivalence is Claimed**

Coagulation S.A.C Assayed Control 1  
Helena Laboratories  
Beaumont, Texas

Docket Number: K941737

Coagulation S.A.C Assayed Control 2  
Helena Laboratories  
Beaumont, Texas

Docket Number: K941872

### **4.0 Description of Device**

Lyphochek® Hemostasis Control is prepared from human plasma with added purified biochemicals and preservatives. The control is provided in lyophilized form for increased stability.

## 5.0 Statement of Intended Use

Lyphocheck® Hemostasis Control is intended for use as a quality control plasma to monitor the precision of laboratory analytes listed in the package insert.

## 6.0 Comparison of the new device with the Predicate Device

Lyphocheck® Hemostasis Control claims substantial equivalence to the Coagulation S.A.C Assayed Control 1 and 2 currently in commercial distribution.

**Table 1. Similarities and Differences between new and predicate device.**

Characteristics	Bio-Rad Lyphocheck® Hemostasis Control (New Device)	Helena Laboratories Coagulation S.A.C.-1 Assayed Control 1 (Predicate Device - K941737)	Helena Laboratories Coagulation S.A.C.-2 Assayed Control 2 (Predicate Device - K941872).
<b>Similarities</b>			
<b>Intended Use</b>	Lyphocheck® Hemostasis Control is intended for use as a quality control plasma to monitor the precision of laboratory testing procedures for analytes listed in the package insert.	S.A.C.-1 is assayed control plasma, which may be used in many facets of testing in the coagulation laboratory.	S.A.C.-2 is assayed abnormal control plasma, which may be used in many facets of testing in the coagulation laboratory.
<b>Form</b>	Lyophilized	Lyophilized	Lyophilized
<b>Matrix</b>	Human plasma based	Human plasma based	Human plasma based
<b>Differences</b>			
<b>Storage (Unopened)</b>	2-8°C Until expiration date	2-6°C Until expiration date	2-6°C Until expiration date
<b>Reconstituted Vial Claim</b>	8 hours at 2°C to 25°C with the following exception:  Protein S will be stable for 8 hours at 2°C to 8°C.	4 hours 2-6°C	4 hours 2-6°C
<b>Analytes</b>	<u>Contains:</u> Prothrombin Time (PT), Activated Partial Thromboplastin Time (APTT), Fibrinogen, Antithrombin III (AT III) Thrombin Time (TT), Factor II, V, VII, VIII, IX, X, XI, XII, Protein S (Functional) and Protein C (Functional), Plasminogen  <u>Does not Contain:</u> Alpha-Antiplasmin, von Willebrand Factor Antigen, Ristocetin cofactor	<u>Contains:</u> Prothrombin Time (PT), Activated Partial Thromboplastin Time (APTT), Fibrinogen, Antithrombin III (AT III), Factor II, V, VII, VIII, IX, X, XI, XII, Total and Free Protein S, Protein C, Plasminogen, von Willebrand Factor Antigen, Alpha-antiplasmin, Ristocetin cofactor  <u>Does not Contain:</u> Thrombin Time (TT)	<u>Contains:</u> Prothrombin Time (PT), Activated Partial Thromboplastin Time (APTT), Fibrinogen, Antithrombin III (AT III), Factor II, V, VII, VIII, IX, X, XI, XII, Total and Free Protein S, Protein C, Plasminogen, Ristocetin cofactor, Alpha-antiplasmin, Monoclonal Free Protein S, von Willebrand Factor Antigen  <u>Does not Contain:</u> Thrombin Time (TT)

## 7.0 Summary of Performance Data

Stability studies have been performed to determine the reconstituted stability and shelf life for the Lyphocheck® Hemostasis Control. Product claims are as follows:

7.1 Reconstituted Stability: Once the control material is reconstituted, all analytes will be stable for 8 hours at 2°C to 25°C with the exception of Protein S which will be stable for 8 hours at 2°C to 8°C.

7.2 Shelf Life: Three years when stored at 2 to 8 °C.

Real time studies will be ongoing to support the shelf life of this product.  
All supporting data is retained on file at Bio-Rad Laboratories.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

**APR 15 2002**

Ms. Maria Zeballos  
Regulatory Affairs Specialist  
Bio-Rad Laboratories  
9500 Jeronimo Road  
Irvine, California 92618-2017

Re: k020878  
Trade/Device Name: Lyphochek® Hemostasis Control  
Regulation Number: 21 CFR § 864.5425  
Regulation Name: Control Plasma Abnormal  
Regulatory Class: II  
Product Code: GGC  
Dated: March 15, 2002  
Received: March 18, 2002

Dear Ms. Zeballos:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style with a large, stylized 'S' and 'G'.

Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510 (k) Number (if known): K020878

Device Name: **Lyphochek® Hemostasis Control**

Indications for Use:

**Lyphochek® Hemostasis Control is intended for use as a quality control plasma to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.**

(PLEASE DO NOT WRITE BELOW THE LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

*Concurrence of CDRH, Office of Device Evaluation (ODE)*

Prescription use ✓ or Over-the Counter use \_\_\_\_\_

*Josephine Bantista*  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K020878