

OCT 25 2002

K022706

**510(k) Summary of
Safety and Effectiveness**

Submitter:

- SPSmedical Supply Corp.
6789 West Henrietta Road
Rush, NY 14543 U.S.A.
Phone: (585)-359-0130 Fax: (585)-359-0167
- Establishment FDA Registration No.: 1319130
- Date Summary was Prepared October 21, 2002
- Gary J. Socola
Printed name of person submitting for 510(k)
- Gary J. Socola
Signature of person submitting for 510(k)
- Director of Quality Assurance
Title of person submitting for 510(k)

Device Name and Classification

Trade Name: SPSmedical SporView® Steam BI Test Pack
Classification Name: Biological Indicator
Common Name: Biological Test Pack
Device Classification: Class II, Regulation Number 880.2800
Product Code: 80FRC
Predicate Device: SteriTec Biological Indicator Test Pack (K001444)

Device Description:

The SPSmedical SporView® Steam BI Test Pack consists of a SporView® steam biological indicator placed inside a package of porous and non-porous material. The SPSmedical SporView® Steam BI Test Pack is designed to create a significant challenge to air removal and steam penetration.

Intended Use:

The SPSmedical SporView® Steam BI Test Pack is indicated for use in routine and challenge testing of steam gravity displacement cycles at 121°C/250°F for 30 minutes exposure time and prevacuum steam sterilization cycles at 132°C/270°F for 4 minutes exposure time or longer and can be used in conjunction with a integrator test pack.

Technical Characteristics:

The SPSmedical SporView® Steam BI Test Pack has the same intended use and technological characteristics as the AAMI biological indicator test pack and other commercially available test packs. The SPSmedical SporView® Steam BI Test Pack is designed to create a significant challenge to air removal and steam penetration. The SPSmedical SporView® Steam BI Test Pack adds resistance and impedes steam penetration to the SporView® steam biological indicator located within the pack. This provides a significant challenge to the steam sterilization process.

Non-Clinical Testing:

Two hundred ten (210) sterilization tests were run to compare performance standards/results of the SPSmedical SporView® Steam BI Test Pack to the AAMI biological indicator test pack. A SporView® self-contained biological indicator containing *Geobacillus stearothermophilus* spores was used within the AAMI biological indicator test. Of the two hundred and ten sterilization tests that were run; 60 comparison tests for failures were orchestrated, 90 comparison tests for pass/failures were orchestrated and 60 tests for passing results were orchestrated. The pass and failure sterilization testing of SPSmedical SporView® Steam BI Pack consistently showed results comparable to the AAMI biological indicator test pack.

Conclusion:

Supportive data has demonstrated that the SPSmedical SporView® Steam BI Pack is equivalent to the predicate device. Results of performance testing indicate that the SPSmedical SporView® Steam BI Pack provides a sufficient load challenge to monitor steam gravity displacement sterilization cycles at 121°C/250°F for 30 minutes exposure time and prevacuum steam sterilization cycles at 132°C/270°F for 4 minutes exposure time or longer. The SPSmedical SporView® Steam BI Pack is an effective and reliable, single use device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

OCT 25 2002

Mr. Gary J. Socola
Director of QA & Sterilization Projects
SPS Medical Supply Corporation
6789 West Henrietta Road
Rush, New York 14543

Re: K022706

Trade/Device Name: Sporview® Steam BI Test Pack
Regulation Number: 880.2800
Regulation Name: Sterilization Process Indicator
Regulatory Class: II
Product Code: FRC
Dated: August 6, 2002
Received: August 14, 2002

Dear Mr. Socola:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.


Page 2 -- Mr. Socola

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Timothy A. Ulatowski
Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS for USE STATEMENT

Applicant: SPSmedical Supply Corp.

510(k) Number (if known): K022706

Device Name: SporView® Steam BI Test Pack


Indications For Use:

The SPSmedical SporView® Steam BI Test Pack is indicated for use in routine and challenge testing of steam gravity displacement cycles at 121°C/250°F for 30 minutes exposure time and prevacuum steam sterilization cycles at 132°C/270°F for 4 minutes exposure time or longer.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Per 21 CFR 801.109)



(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K022706