

FEB 26 2004

## 510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: **K033479**.

**A. Contact:**

David Teicher  
Remel Inc. Ramsey Operations  
14000 Unity Street NW  
Ramsey, MN  
Phone: (763) 712-2367  
Fax: (763) 323-7858  
e-mail: [dteicher@remel.com](mailto:dteicher@remel.com)

**B. Date of Preparation:** January 19, 2004

**C. Name of Device and Classification**

Name: ProSpecT® Clostridium difficile Toxin A/B Microplate Assay  
Classification: 866.2660 Reagents, Clostridium difficile toxin, Class I, LLH

**D. Predicate Devices**

Meridian Premier Toxins A&B: 510(k) K993914  
Wampole *C. difficile* TOX A/B II: 510(k) K003306/K971182

**E. Device Description**

The ProSpecT® Clostridium difficile Toxin A/B test detects the presence of Toxin A and Toxin B in clinical stool specimens through the use of specific antibodies. Microwell strips are coated with mouse monoclonal anti-Toxin A and rabbit anti-Toxin B antibodies. A stool specimen is diluted in Sample Diluent or used directly if pre-diluted in modified Cary-Blair medium. The sample is added to a microwell allowing the toxins, if present, to bind to the immobilized antibodies. After washing to remove unbound components, a conjugate reagent containing goat anti-Toxin A-HRP and rabbit anti-Toxin B-HRP is added to each well. Unbound conjugate is removed by washing and a chromagenic substrate solution is added to detect the presence of bound toxin. A stop reagent is added and the test results are read visually or spectrophotometrically. The presence of a yellow color indicates the presence of toxin.

## F. Intended Use Statement

Remel's ProSpecT® Clostridium difficile Toxin A/B Microplate Assay is a qualitative enzyme immunoassay (EIA) for the detection of *C. difficile* Toxin A and B in human fecal specimens from patients suspected of having *Clostridium difficile* disease. The test is intended for use as an aid in diagnosis of *Clostridium difficile*-associated disease (CDAD).

## G. Device Comparison

Product Feature	ProSpecT®	Wampole	Merdian
Intended Use	Detection of <i>C. difficile</i> Toxins A and B in fecal specimens	Detection of <i>C. difficile</i> Toxins A and B in fecal specimens	Detection of <i>C. difficile</i> Toxins A and B in fecal specimens
Technology	Enzyme Immunoassay	Enzyme Immunoassay	Enzyme Immunoassay
Capture Antibodies or Molecules: Device	Mouse monoclonal anti-Toxin A and rabbit anti-Toxin B	Polyclonal goat antibody against Toxins A and B	Mouse monoclonal anti-Toxin A and polyclonal goat anti-Toxin B
Antibodies: Conjugate	Goat anti-Toxin A and rabbit anti-Toxin B	Toxin A monoclonal mouse antibody and Toxin B polyclonal goat antibody	Polyclonal goat anti-Toxin A and anti-Toxin B
Material: Device	Microwell	Microwell	Microwell
Material: Conjugate	Horseradish peroxidase conjugated to anti-toxins	Horseradish peroxidase conjugated to anti-toxins	Horseradish peroxidase conjugated to anti-toxins
Specimen Type	Fresh human stool specimens or specimens in modified Cary-Blair	Fresh human stool specimens	Fresh human stool specimens
Sample volume	200 µl	50 µl	50 µl

## H. Summary of Performance Data

The performance of the ProSpecT® Clostridium difficile Toxin A/B Microplate Assay was evaluated in a clinical trial at three sites in the USA. The kit and predicate devices were compared to the cellular cytotoxicity assay (CTA).

### Overall Performance:

#### Dual Wavelength (450/620-650 nm) Results:

Sensitivity:  $(149/165) \times 100 = 90.3\%$ ; 95% CI = 84.7% - 94.4%

Specificity:  $(576/599) \times 100 = 96.2\%$ ; 95% CI = 94.3% - 97.5%

		CTA Results	
		+	-
ProSpecT® EIA Results (Dual)	+	149	23
	-	16	576
Total		165	599

### Visual Interpretation of Test:

Sensitivity:  $(85/100) \times 100 = 85.0\%$ ; 95% CI = 76.5% - 91.4%

Specificity:  $(464/486) \times 100 = 95.5\%$ ; 95% CI = 93.2% - 97.1%

		CTA Results	
		+	-
ProSpecT® EIA Visual Results	+	85	22
	-	15	464
Total		100	486

Percent Agreement:  $(580/586) \times 100 = 99.0\%$

		ProSpecT® EIA Results OD 450/620-650	
		+	-
ProSpecT® EIA Visual Results	+	104	3
	-	3	476
Total		107	479

### Comparison to Predicate Devices:

EIA	Performance versus CTA			
	Sensitivity		Specificity	
	#	%	#	%
ProSpecT®	33/40	82.5	263/268	98.1
Predicate 1	33/40	82.5	260/268	97.0
ProSpecT®	115/124	92.7	302/320	94.3
Predicate 2	98/124	79.0	309/320	96.5

#### I. Reproducibility

Reproducibility testing was conducted at three sites on three separate days with four blinded samples. Each site tested eight replicate wells of each specimen on each day of testing (n=288). The specimens included one negative specimen and three positive specimens with varying levels of reactivity. The average intra-assay coefficient of variation (CV) for a mid-range sample was 7.7%. The average inter-assay CV for a mid-range sample was 18.9%.

#### J. Analytical Sensitivity

The ProSpecT® Clostridium difficile Toxin A/B Microplate Assay detects Toxin A at levels of  $\geq 0.20$  ng/ml and Toxin B at levels of  $\geq 0.61$  ng/ml.

#### K. Cross-Reactivity

Forty (40) microorganisms were evaluated with the ProSpecT® Clostridium difficile Toxin A/B Microplate Assay. Bacteria and yeast isolates were tested at  $\geq 10^8$  colony-forming units per ml (cfu/ml). Viral isolates were tested at concentrations of  $10^4$  TCID<sub>50</sub>/ml (tissue culture infectious dose per milliliter). No cross-reactivity was observed. There was no cross-reactivity observed with *Clostridium sordellii* ATCC® 9714. However, published literature indicates that certain strains of *C. sordellii* can produce toxins which may be cross-reactive with antibodies to *C. difficile* Toxins A and B. The following organisms were tested in the ProSpecT® Clostridium difficile Toxin A/B Microplate Assay.

<u>Organism</u>	<u>ATCC® designation</u>
<i>Aeromonas hydrophila</i>	35654
<i>Bacillus cereus</i>	11778
<i>Bacillus subtilis</i>	6633
<i>Bacteriodes fragilis</i>	25285
<i>Campylobacter coli</i>	33559
<i>Campylobacter jejuni</i>	33291
<i>Candida albicans</i>	10231

<u>Organism</u>	<u>ATCC® designation</u>
<i>Clostridium beijerinckii</i> (butyricum)	8260
<i>Clostridium difficile</i> (non-toxigenic)	700057
<i>Clostridium haemolyticum</i>	9650
<i>Clostridium histolyticum</i>	19401
<i>Clostridium novyi</i> (toxin A)	19402
<i>Clostridium perfringens</i> (Type A)	13124
<i>Clostridium septicum</i>	12464
<i>Clostridium sordellii</i>	9714
<i>Clostridium sporogenes</i>	19404
<i>Clostridium tetani</i>	19406
<i>Enterobacter aerogenes</i>	35028
<i>Enterobacter cloacae</i>	13047
<i>Enterococcus faecalis</i>	19433
<i>Escherichia coli</i>	11229
<i>Klebsiella pneumoniae</i>	13882
<i>Peptostreptococcus anaerobius</i>	27337
<i>Porphyromonas asaccharolytica</i>	25260
<i>Proteus vulgaris</i>	49132
<i>Pseudomonas aeruginosa</i>	27853
<i>Salmonella choleraesuis</i> (typhimurium)	23852
<i>Serratia liquefaciens</i>	27592
<i>Shigella dysenteriae</i>	11835
<i>Shigella flexneri</i>	12022
<i>Shigella sonnei</i>	25931
<i>Staphylococcus aureus</i>	25923
<i>Staphylococcus aureus</i> (Cowan)	12598
<i>Staphylococcus epidermidis</i>	12228
<i>Vibrio cholerae</i>	9459
<i>Vibrio parahaemolyticus</i>	17802
<i>Yersinia enterocolitica</i>	23715
Adenovirus type 40	VR-930
Adenovirus type 41	VR-931
Rotavirus (Complement fixation antigen)	

#### **L. Interfering Substances**

The following substances were tested with the ProSpecT® Clostridium difficile Toxin A/B Microplate Assay: Vancomycin (12.5 mg/ml), Metronidazole (12.5 mg/ml), blood, mucous, fecal fat, and the following over-the-counter anti-diarrheal products: Pepto-Bismol®, Imodium® A-D, and Kaopectate® (active ingredients: bismuth subsalicylate, loperamide HCl, and attapugite respectively). No interference with positive or negative specimens was observed.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

FEB 26 2004

Mr. David Teicher  
Manager Quality Systems and Regulatory Affairs  
Remel Inc.  
Ramsey Operation  
14000 Unity Street NW  
Ramsey, MN 55303

Re: k033479  
Trade/Device Name: ProSpecT<sup>®</sup> *C. difficile* Toxin A/B Microplate Assay  
Regulation Number: 21 CFR 866.2660  
Regulation Name: Microorganism differentiation and identification device  
Regulatory Class: Class I  
Product Code: LLH  
Dated: January 19, 2004  
Received: January 20, 2004

Dear Mr. Teicher:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Sally A. Hojvat, M.Sc., Ph.D.  
Director  
Division of Microbiology Devices  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K033479

Device Name: ProSpecT<sup>®</sup> *C. difficile* Toxin A/B Microplate Assay

### Indications For Use:

Remel's ProSpecT<sup>®</sup> *Clostridium difficile* Toxin A/B Microplate Assay is a qualitative enzyme immunoassay (EIA) for the detection of *C. difficile* Toxin A and B in human fecal specimens from patients suspected of having *Clostridium difficile* disease. The test is intended for use as an aid in diagnosis of *Clostridium difficile*-associated disease (CDAD).

For *In vitro* Diagnostic Use.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)



Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k)   K033479