

FEB 13 2004

**Section 3**  
**HemosIL Factor VIII Deficient Plasma - 510(k) Summary**  
**(Summary of Safety and Effectiveness)**

**Submitted by:**

Instrumentation Laboratory Company  
113 Hartwell Avenue  
Lexington, MA 02421  
Phone: 781-861-4467  
Fax: 781-861-4207

**Contact Person:**

Carol Marble, Regulatory Affairs Director  
Phone: 781-861-4467 / Fax: 781-861-4207

**Summary Prepared:**

December 23, 2003

**Name of the Device:**

HemosIL Factor VIII Deficient Plasma

**Classification Name(s):**

864.7290 Factor Deficiency Tests Class II  
81GJT Plasma, Coagulation Factor Deficient

**Identification of Predicate Device(s):**

K893525 Hemoliance Factor VIII Deficient Plasma on ELECTRA Series Analyzers  
K002400 IL Test Factor VIII Deficient Plasma\* on ACL Family of Analyzers  
\*NOTE: Reagent was 510(k) cleared as part of multiple analyzer systems, most recently the ACL Advance.

**Description of the Device/Intended use(s):**

HemosIL Factor VIII Deficient Plasma is human plasma immunodepleted of factor VIII and intended for the *in vitro* diagnostic quantitative determination of factor VIII activity in citrated plasma, based on the activated partial thromboplastin time (APTT) assay, on IL Coagulation and ELECTRA Systems.

Abnormalities of the intrinsic pathway factors are determined by performing a modified activated partial thromboplastin time (APTT) test. Patient plasma is diluted and added to a plasma deficient in factor VIII. Correction of the clotting time of the deficient plasma is proportional to the concentration (% activity) of the factor VIII in the patient plasma, interpolated from a calibration curve.

**Statement of Technological Characteristics of the Device Compared to Predicate Device:**

HemosIL Factor VIII Deficient Plasma is substantially equivalent to Hemoliance Factor VIII Deficient Plasma (on ELECTRA Series Analyzers) and IL Test Factor VIII Deficient Plasma (on ACL Family of Analyzers) in performance, intended use and safety and effectiveness.

### Section 3

## HemosIL Factor VIII Deficient Plasma - 510(k) Summary (Summary of Safety and Effectiveness)

### Summary of Performance Data:

#### Method Comparison

In method comparison studies evaluating approximately 60 citrated plasma samples (30 normal and 30 abnormal), the slopes and correlation coefficients (r) for HemosIL Factor VIII Deficient Plasma versus the predicate devices are shown below:

NOTE: APTT-SP and SynthASil were used as the APTT reagents in testing.

#### HemosIL Factor VIII Deficient Plasma vs. Predicate Hemoliance Factor VIII Deficient Plasma on ELECTRA

IL System	n	Slope	r
E1400C	59	0.9518	0.9873

#### HemosIL Factor VIII Deficient Plasma vs. Predicate IL Test Factor VIII Deficient Plasma on ACL Family

IL System	n	Slope	r
ACL 300	60	0.9391	0.9942
ACL Advance	63	1.0073	0.9906

#### Within Run Precision

Within run and between run precision was assessed over multiple runs (n=80) on different instruments using a specific lot of APTT reagent (APTT-SP and SynthASil) and both normal and abnormal samples.

Instrument	Control	Mean % Factor VIII	Within run CV%	Between Run CV%
<b>ACL 9000</b>	Normal Control	72.6	3.1	3.1
	Low Abnormal Control	32.2	2.2	3.3
<b>ACL Futura</b>	Normal Control	78.3	2.7	3.6
	Low Abnormal Control	31.5	3.4	4.0
<b>ELECTRA 1600C</b>	Normal Control	93.9	4.8	6.0
	Low Abnormal Control	37.2	4.1	4.6



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Ms. Carol Marble  
Regulatory Affairs Director  
Instrumentation Laboratory Company  
113 Hartwell Avenue  
Lexington, Massachusetts 02421

FEB 13 2004

Re: k034007  
Trade/Device Name: HemosIL Factor VIII Deficient Plasma  
Regulation Number: 21 CFR § 864.7290  
Regulation Name: Plasma, Coagulation Factor Deficiency Test  
Regulatory Class: II  
Product Code: GJT  
Dated: December 23, 2003  
Received: December 24, 2003

Dear Ms. Marble:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink, reading "Joseph L. Hackett". The signature is written in a cursive style with a large, stylized "J" and "H".

Joseph L. Hackett, Ph.D.  
Acting Director  
Division of Immunology and Hematology Devices  
Office of In Vitro Diagnostic Device Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Indications for Use Statement

510(k) Number (if known): K03 4007

**Device Name:** HemosIL Factor VIII Deficient Plasma

### Indications for Use:

HemosIL Factor VIII Deficient Plasma is human plasma immunodepleted of factor VIII and intended for the *in vitro* diagnostic quantitative determination of factor VIII activity in citrated plasma, based on the activated partial thromboplastin time (APTT) assay, on IL Coagulation and ELECTRA Systems.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k) K03 4007

Prescription Use ☒

OR

Over-The-Counter Use ☐