



JAN 12 2005

SonoSite, Inc.  
21919 30<sup>th</sup> Drive SE  
Bothell, WA 98021-3904 USA

Telephone: 1.425.951.1200  
Facsimile: 1.425.951.1201  
[www.sonosite.com](http://www.sonosite.com)

## 510(K) Summary of Safety and Effectiveness

This summary of safety and effectiveness is provided as part of this Premarket Notification in compliance with 21 CFR, Part 807, Subpart E, Section 807.92.

### 1) Submitter's name, address, telephone number, contact person:

SonoSite, Inc.  
21919 30<sup>th</sup> Drive SE  
Bothell, WA 98021-3904

**Corresponding Official:** Daina L. Graham  
Vice President, Regulatory Affairs and Quality Assurance  
**E-mail:** [Daina.Graham@sonosite.com](mailto:Daina.Graham@sonosite.com)  
**Telephone:** (425) 951-1275  
**Facsimile:** (425) 951-1201  
**Date prepared:** December 9, 2004

### 2) Name of the device, including the trade or proprietary name if applicable, the common or usual name, and the classification name, if known:

#### Common/ Usual Name

Diagnostic Ultrasound System with Accessories

#### Proprietary Name

SonoSite High-Resolution Ultrasound System (C3 Series) (*subject to change*)

#### Classification Names

Name	FR Number	Product Code
Ultrasonic Pulsed Doppler Imaging System	892.1550	90-IYN
Ultrasonic Pulsed Echo Imaging System	892.1560	90-IYO
Diagnostic Ultrasound Transducer	892.1570	90-ITX
Picture Archiving And Communications System	892.2050	90-LLZ

### 3) Identification of the predicate or legally marketed device:

SonoSite, Inc. believes that the System described in this Submission is substantially equivalent to a combination of the SonoSite TITAN™ Ultrasound System (K030949), Philips Ultrasound HDI® 5000 Ultrasound System (K011224), the GE LOGIQ 9 Diagnostic Ultrasound System (K011188) and the Siemens OMNIA X/XS Diagnostic Ultrasound System (K020353).

#### 4) Device Description:

The SonoSite High-Resolution Ultrasound System (C3 Series) is a full featured, general purpose, software controlled, diagnostic ultrasound system used to acquire and display high-resolution, real-time ultrasound data in 2D, M-Mode, Pulsed Wave (PW) Doppler, Continuous Wave (CW) Doppler, Color Power Doppler, and Velocity Color Doppler or in a combination of these modes.

The System has an electrocardiography (ECG) display feature and supports a 3-lead ECG cable assembly to collect data for M-mode and Doppler measurements. The System provides measurement capabilities for anatomical structures and fetal biometry that provide information used for clinical diagnostic purposes. The System has a PW and CW Doppler audio output feature and cine review, image zoom, labeling, biopsy, measurements and calculations, image storage and review, printing, and recording capabilities.

The system includes the ability to measure the intima-media thickness (IMT) of the carotid artery using digital ultrasound images. The IMT measurement of the carotid artery may be used adjunctively with other medical data obtained by a physician to help assess the cardiovascular health of a patient.

The system includes Digital Imaging and Communications (DICOM) capabilities as well as general computer communication capabilities to provide the acceptance, transfer, display, storage, and digital processing of ultrasound images and loops. Security support is also provided to facilitate HIPAA compliance.

The SonoSite High-Resolution Ultrasound System (C3 Series) is designed to accept curved or linear transducers of the types and frequency listed in the table below. All actions affecting the performance of the transducer are activated from the main system control panel.

<b>Frequency Range:</b>	<b>2.0 - 12.0 MHz</b>
<b>Transducer Types:</b>	Linear array Curved array Intracavitary array Phased array Static probes

The SonoSite High-Resolution Ultrasound System (C3 Series) is designed to comply with the standards listed below.

Reference no.	Title
EN 540	Clinical investigation of medical devices for human subjects (1993)
EN 980 A1	Graphical symbols for use in the labeling of medical devices (2003)
EN 1041	Information supplied by the manufacturer with medical devices (1998)
EN ISO 13485	Medical devices – Quality management systems – Requirements for regulatory purposes (2003)
EN ISO 14971	Medical devices – Application of risk management to medical devices (2000) (ISO 14971:2000) (Superseded standard: EN 1441)
EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing (1997) (ISO 10993-1:1997)
EN 60601-1	Medical electrical equipment. Part 1: General requirements for safety - IEC 601-1:1988
EN 60601-1	Medical electrical equipment. Part 1: General requirements for safety - IEC 601-

Reference no.	Title
Amendment A1	1:1988/A1:1991
EN 60601-1 Amendment A2	Medical electrical equipment. Part 1: General requirements for safety - IEC 601-1:1988/A2:1995 + corrigendum June 1995
EN 60601-1-1	Medical electrical equipment. Part 1: General requirements for safety - 1. Collateral standard: Safety requirements for medical electrical systems (2001) (IEC 601-1-1:2001)
EN 60601-1-2	Medical electrical equipment – Part 1: General requirements for safety; 2. collateral standard: electromagnetic compatibility; requirements and tests (2001) (IEC 60601-1-2:2001)
EN 60601-2-37	Medical electrical equipment -- Part 2-37: Particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment (2001) (IEC 60601-2-37:2001)
-	Medical Ultrasound Safety, American Institute of Ultrasound in Medicine (AIUM), 1994
-	Acoustic Output Measurement and Labeling Standard for Diagnostic Ultrasound Equipment, American Institute of Ultrasound in Medicine, 1993
-	Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment, American Institute of Ultrasound in Medicine, 1998
ANSI/AAMI EC53:1995	ECG Cables and Electrodes except for sections 4.4 and 4.5.9
ASTM D5276-98	Standard Test Methods for Drop Test of Loaded Containers by Free Fall
ASTM D999-96	Standard Methods for Vibration Testing of Shipping Containers
CEI/IEC 61157:1992	International Electrotechnical Commission, Requirements for the Declaration of the Acoustic Output of Medical Diagnostic Ultrasonic Equipment
D00017	Quality System Manual (SonoSite)
EN 60529	Degrees of protection provided by enclosures (IP Code) (1991)
ISO 9001: 1994	Quality Systems – Model for quality assurance in design, development, production, installation and servicing
JIS T 0601-1	Medical Electrical Equipment – Part 1: General Requirements for Safety (Japan) (1999)
JIS T 1507	Electronic Linear Scanning Ultrasonic Diagnostic Equipment (Japan) (1989)
NEMA UD 2-2004	Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment
NEMA UD 3-2004	Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment, American Institute of Ultrasound in Medicine
RTCA/DO160D:1997	Radio Technical Commission for Aeronautics, Environmental Conditions and Test Procedures for Airborne Equipment, Section 21.0 Emission of Radio Frequency Energy, Category B
Title 21 CFR Part 820	Quality System Regulation – Medical Devices: Current Good Manufacturing Practice (CGMP); Final Rule
UL 60601-1	Underwriters Laboratories, Medical Electrical Equipment-Part 1: General Requirements for Safety (2003)
ACR-NEMA DICOM version 3.0 - 2003	Digital Imaging and Communications in Medicine

**5) Intended Use:**

The intended uses of the SonoSite High-Resolution Ultrasound System (C3 Series) as defined by FDA guidance documents, are:

Fetal - OB/GYN	Trans-rectal
Abdominal	Trans-vaginal
Intra-operative (Abdominal organs and vascular)	Trans-urethral
Intra-operative (Neuro.)	Musculo-skel. (Conventional)
Laparoscopic	Musculo-skel. (Superficial)
Pediatric	Cardiac Adult
Small Organ (breast, thyroid, testicles.)	Cardiac Pediatric
Neonatal Cephalic	Trans-esophageal (card.)
Adult Cephalic	Peripheral vessel

Typical examinations performed using the SonoSite High-Resolution Ultrasound System (C3 Series) are:

**Abdominal Imaging Applications**

This system transmits ultrasound energy into the abdomen of patients using 2D, M Mode, color Doppler (Color), color power Doppler (CPD), Tissue Harmonic Imaging (THI), and pulsed wave (PW) Doppler to obtain ultrasound images. The liver, kidneys, pancreas, spleen, gallbladder, bile ducts, transplanted organs, abdominal vessels, and surrounding anatomical structures can be assessed for the presence or absence of pathology transabdominally.

**Cardiac Imaging Applications**

This system transmits ultrasound energy into the thorax of patients using 2D, M Mode, color Doppler (Color), Tissue Harmonic Imaging (THI), pulsed wave (PW) Doppler, and continuous wave (CW) Doppler to obtain ultrasound images. The heart muscle, cardiac valves, great vessels, and surrounding anatomical structures can be assessed for, overall cardiac performance, size, fluid, and the presence or absence of pathology. Quantification of size, function, and hemodynamic performance can also be assessed. The heart can be imaged transthoracically or transesophageally. Transesophageal imaging can be used to monitor the heart during surgical procedures such as cardiac bypass or valve replacements.

The patient's electrocardiogram (ECG) may be obtained and is used for accurate timing of diastolic and systolic function.

Warning: The ECG is not used to diagnose cardiac arrhythmias and is not designed for long term cardiac rhythm monitoring.

**Gynecology and Infertility Imaging Applications**

This system transmits ultrasound energy in the pelvis and lower abdomen using 2D, M Mode, color power Doppler (CPD), color Doppler (Color) Tissue Harmonic Imaging (THI), and pulsed wave (PW) Doppler to obtain ultrasound images. The uterus, ovaries, adnexa, and surrounding anatomical structures can be assessed for the presence or absence of pathology transabdominally or transvaginally.

### **Interventional and Intraoperative Imaging Applications**

This system transmits ultrasound energy into the various parts of the body using 2D, M-Mode, color Doppler (Color), color power Doppler (CPD), Tissue Harmonic Imaging (THI), and pulsed wave (PW) Doppler to obtain ultrasound images that provide guidance during interventional and intraoperative procedures. This system can be used to provide ultrasound guidance for biopsy and drainage procedures, vascular line placement, peripheral nerve blocks, ova harvesting, amniocentesis and other obstetrical procedures, and provide assistance during abdominal, vascular, and neurological intraoperative procedures.

Warning: This system is not intended for use in providing guidance for central nerve blocks, i.e., the brain and spinal cord, or for ophthalmic applications.

### **Obstetrical Imaging Applications**

This system transmits ultrasound energy into the pelvis of pregnant women using 2D, M Mode, color Doppler (Color), color power Doppler (CPD), Tissue Harmonic Imaging (THI), and pulsed wave (PW) Doppler to obtain ultrasound images. The fetal anatomy, viability, estimated fetal weight, gestational age, amniotic fluid, and surrounding anatomical structures can be assessed for the presence or absence of pathology transabdominally or transvaginally. CPD and color Doppler (Color) imaging is intended for high-risk pregnant women. High-risk pregnancy indications include, but are not limited to, multiple pregnancy, fetal hydrops, placental abnormalities, as well as maternal hypertension, diabetes, and lupus.

Warning: To prevent injury or misdiagnosis do not use this system for Percutaneous Umbilical Blood Sampling (PUBS) or *in vitro* Fertilization (IVF) The system has not been validated to be proven effective for these two uses. CPD, or Color images can be used as an adjunctive method, not as a screening tool, for the detection of structural anomalies of the fetal heart and as an adjunctive method, not as a screening tool for the diagnosis of Intrauterine Growth Retardation (IUGR).

### **Pediatric and Neonatal Imaging Applications**

This system transmits ultrasound energy into the pediatric or neonatal patients using 2D, M Mode, color Doppler (Color), color power Doppler (CPD), pulsed wave (PW) Doppler, and continuous wave (CW) Doppler to obtain ultrasound images. The pediatric abdominal, pelvic and cardiac anatomy, pediatric hips, neonatal heads, and surrounding anatomical structures can be assessed for the presence or absence of pathology.

### **Prostate Imaging Applications**

This system transmits ultrasound energy into the prostate of an adult male using 2D, M Mode, color power Doppler (CPD), and pulsed wave (PW) Doppler to obtain ultrasound images. The prostate gland can be assessed for the presence or absence of pathology.

### Superficial Imaging Applications

This system transmits ultrasound energy into various parts of the body using 2D, M Mode, color Doppler (Color), color power Doppler (CPD), and pulsed wave (PW) Doppler to obtain ultrasound images. The breast, thyroid, testicle, lymph nodes, hernias, musculoskeletal structures, soft tissue structures, and surrounding anatomical structures can be assessed for the presence or absence of pathology. This system can be used to provide ultrasound guidance for biopsy and drainage procedures, vascular line placement, and peripheral nerve blocks.

Warning: This system is not intended for use in providing guidance for central nerve blocks, i.e., the brain and spinal cord, or for ophthalmic applications.

### Vascular Imaging Applications

This system transmits ultrasound energy into the various parts of the body using 2D, M Mode, color Doppler (Color), color power Doppler (CPD), and pulsed wave (PW) Doppler to obtain ultrasound images. The carotid arteries, deep veins in the arms and legs, superficial veins in the arms and legs, great vessels in the abdomen, and various small vessels feeding organs can be assessed for the presence or absence of pathology.

## 6) Technological Characteristics:

This device operates identically to the predicate devices in that piezoelectric material in the transducer is used as an ultrasound source to transmit sound waves into the body. Sound waves are reflected back to the transducer and converted to electrical signals that are processed and displayed as 2D or M-mode images. Doppler shift caused by blood flow is displayed as Color Flow, or as spectrum analysis. The modes of this device (2D, PW Doppler, CW Doppler, velocity color Doppler, Color Power Doppler, and duplex imaging) are the same as a combination of the predicate devices identified in item 3. Transducer patient contact materials are biocompatible.

This device conforms to the *Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment* (AIUM/NEMA UD 3-2004) for an on-screen display feature that provides information on potential thermal and cavitation bioeffect mechanisms. A user education program provides additional information so users may moderate the system's acoustic output in accordance with the ALARA (as low as reasonably achievable) principle.

The device's acoustic output limits are:

All applications:

$I_{SPTA}$ (d)	720 mW/cm <sup>2</sup> (Maximum)
TIs/TIb/TIc	0.1 - 4.0 (Range)
Mechanical Index (MI)	1.9 (Maximum)
$I_{SPPA}$ (d)	0 - 700 W/cm <sup>2</sup> (Range)

The limits are the same as predicate Track 3 devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JAN 12 2005

SonoSite, Inc.  
% Mr. Mark Job  
Responsible Third Party Official  
Regulatory Technology Services LLC  
1394 25<sup>th</sup> Street NW  
BUFFALO MN 55313

Re: K043559

Trade Name: SonoSite High-Resolution Diagnostic Ultrasound System (C3 Series)  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulation Number: 21 CFR 892.1560  
Regulation Name: Ultrasonic pulsed echo imaging system  
Regulation Number: 21 CFR 892.1570  
Regulation Name: Diagnostic ultrasonic transducer  
Regulatory Class: II  
Product Code: 90 IYN, IYO, and ITX  
Dated: December 20, 2004  
Received: December 27, 2004

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the SonoSite High-Resolution Diagnostic Ultrasound System (C3 Series), as described in your premarket notification:

Transducer Model Number

2MHZ Doppler Pencil  
C8/8-5

C11/8-5  
C15/4-2  
C60/5-2  
C60e/5-2  
HFL38/12-5  
HST/10-5  
HST/12-7  
HST/12-8  
ICT/8-5  
LAP/10-5  
L25/10-5  
L38/10-5  
L38e/10-5  
L52/10-5  
L52-S/10-5  
MPT/7-4  
P/4-2  
P/7-4

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:




Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,

*for* 

Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure(s)

Table 4.3- 1 Diagnostic Ultrasound Indications for Use Form – C3 System

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		N/A						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
Fetal Imaging	Laparoscopic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
& Other	Pediatric	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Adult Cephalic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-rectal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-vaginal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
Cardiac	Cardiac Pediatric	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Trans-esophageal (card.)	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures, and picture archiving, communications and storage functionality were all previously cleared in K030949. An expanded intended use for imaging guidance for peripheral nerve block procedures was previously cleared in K033367.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number

K043559

Table 4.3- 2 Diagnostic Ultrasound Indications for Use Form - 2 MHz Doppler Pencil

<b>System:</b>	SonoSite High-Resolution Ultrasound System (C3 Series)							
<b>Transducer:</b>	Doppler Pencil Transducer 2.0 MHz Dual Element Circular Array							
<b>Intended Use:</b>	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:							
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Abdominal organs and vascular)							
	Intra-operative (Neuro.)							
Fetal Imaging & Other	Laparoscopic							
	Pediatric							
	Small Organ (breast, thyroid, testicles.)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non- Card.)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult				N			
Cardiac	Cardiac Pediatric				N			
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Prescription Use (Per 21 CFR 801.109)

*Erin A. Seymour*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K043559

Table 4.3- 3 Diagnostic Ultrasound Indications for Use Form – C8/8-5 Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		C8/8-5 8.0-5.0 MHz Curved Array Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal							
	Intra-operative (Abdominal organs and vascular)							
	Intra-operative (Neuro.)							
Fetal Imaging & Other	Laparoscopic							
	Pediatric							
	Small Organ (breast, thyroid, testicles.)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-vaginal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-urethral							
	Trans-esoph. (non- Card.)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was cleared in K030949.

Prescription Use (Per 21 CFR 801.109)

*David A. Syron*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K043559

Table 4.3- 4 Diagnostic Ultrasound Indications for Use Form - C11/8-5 Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		C11/8-5 8.0–5.0 MHz Curved Array						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
Fetal Imaging	Laparoscopic							
	Pediatric	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
& Other	Small Organ (breast, thyroid, testicles)							
	Neonatal Cephalic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non- Card.)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in 510(k) K030949. An expanded intended use for imaging guidance for peripheral nerve block procedures was previously cleared in K033367.

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K043559

Prescription Use (Per 21 CFR 801.109)

Indications for Use

Table 4.3- 5 Diagnostic Ultrasound Indications for Use Form - C15/4-2 Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		C15/4-2 4.0–2.0 MHz Curved Array						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Adult Cephalic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non- Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
Cardiac	Cardiac Pediatric	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in 510(k) K030949.

Prescription Use (Per 21 CFR 801.109)

Indications for Use

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

28 of 905 pages K043559 Section 4.3

Table 4.3- 6 Diagnostic Ultrasound Indications for Use Form - C60/5-2 Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		C60/5-2 5.0-2.0 MHz Curved Array						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in 510(k) K030949.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

*Leah A. Lyman*

K043559

Table 4.3- 7 Diagnostic Ultrasound Indications for Use Form - C60e/5-2 Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		C60e/5-2 5.0-2.0 MHz Curved Array						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in 510(k) K030949.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

*David A. Segura*

K043559



Table 4.3- 8 Diagnostic Ultrasound Indications for Use Form - HFL38/12-5 Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		HFL38/12-5 12.0-5.0 MHz Linear Array Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Includes imaging to assist in the placement of needles and catheters in vascular or other anatomical structures and an expanded intended use for imaging guidance for peripheral nerve block procedures.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K043559

Table 4.3- 9 Diagnostic Ultrasound Indications for Use Form – HST/10-5 Transducer

<b>System:</b>	SonoSite High-Resolution Ultrasound System (C3 Series)							
<b>Transducer:</b>	HST/10-5 10.0-5.0 MHz Linear Array Transducer							
<b>Intended Use:</b>	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:							
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode.

Prescription Use (Per 21 CFR 801.109)

*David A. Lynn*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K041559

**Table 4.3- 10 Diagnostic Ultrasound Indications for Use Form – HST/12-7 Transducer**

<b>System:</b>	SonoSite High-Resolution Ultrasound System (C3 Series)							
<b>Transducer:</b>	HST/12-7 12.0-7.0 MHz Linear Array Transducer							
<b>Intended Use:</b>	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:							
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic						B+M; B+PWD; B+CD	Note 1
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

KD43559

**Table 4.3- 11 Diagnostic Ultrasound Indications for Use Form – HST/12-8 Transducer**

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		HST/12-8 12.0-8.0 MHz Linear Array Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic						B+M; B+PWD; B+CD	Note 1
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

#### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode.

Prescription Use (Per 21 CFR 801.109)

*David H. Syman*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K043559

**Table 4.3- 12 Diagnostic Ultrasound Indications for Use Form – ICT/8-5 Transducer**

<b>System:</b>	SonoSite High-Resolution Ultrasound System (C3 Series)							
<b>Transducer:</b>	ICT/8-5 8.0-5.0 MHz Intracavitary Transducer							
<b>Intended Use:</b>	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:							
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal							
	Intra-operative (Abdominal organs and vascular)							
	Intra-operative (Neuro.)							
Fetal Imaging & Other	Laparoscopic							
	Pediatric							
	Small Organ (breast, thyroid, testicles.)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-vaginal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in 510(k) K030949.

Prescription Use (Per 21 CFR 801.109)

*David A. Byrnes*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number R043559

**Table 4.3- 13 Diagnostic Ultrasound Indications for Use Form – LAP/10-5  
Laparoscopic Transducer**


<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		LAP/10-5 10.0-5.0 MHz Linear Array Laparoscopic Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Abdominal organs and vascular)							
	Intra-operative (Neuro.)							
Fetal Imaging & Other	Laparoscopic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Pediatric							
	Small Organ (breast, thyroid, testicles.)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non- Card.)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

#### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode.

Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number RD 43559

**Table 4.3- 14 Diagnostic Ultrasound Indications for Use Form – L25/10-5 Transducer**

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		L25/10-5 10.0-5.0 MHz Linear Array						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal	N	N			N	B+M; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N				B+M; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N			N	B+M; B+CWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles)	N	N			N	B+M; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N			N	B+M; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N			N	B+M; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N			N	B+M; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

#### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in K030949. An expanded intended use for imaging guidance for peripheral nerve block procedures was previously cleared in K033367.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K0413559

Table 4.3- 15 Diagnostic Ultrasound Indications for Use Form - L38/10-5 Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		L38/10-5 10.0-5.0 MHz Linear Array Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in K030949. An expanded intended use for imaging guidance for peripheral nerve block procedures was previously cleared in K033367.

Prescription Use (Per 21 CFR 801.109)

*David A. Legman*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number KD48559



**Table 4.3- 16 Diagnostic Ultrasound Indications for Use Form - L38e/10-5 Transducer**

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		L38e/10-5 10.0-5.0 MHz Linear Array Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
& Other	Small Organ (breast, thyroid, testicles)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non- Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

#### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Includes imaging to assist in the placement of needles and catheters in vascular or other anatomical structures and imaging guidance for peripheral nerve block procedures.

Prescription Use (Per 21 CFR 801.109)

*David A. Segerson*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number RD 43559

**Table 4.3- 17 Diagnostic Ultrasound Indications for Use Form – L52/10-5 Transducer**

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		L52/10-5 10.0-5.0 MHz Linear Array						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

#### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. Imaging to assist in the placement of needles and catheters in vascular or other anatomical structures was previously cleared in K030949.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number

*David A. Depina*

K041359

**Table 4.3- 18 Diagnostic Ultrasound Indications for Use Form – L52-S/10-5 Transducer**

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 Series)						
<b>Transducer:</b>		L52-S/10-5 10.0-5.0 MHz Linear Array						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Abdominal	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging	Laparoscopic							
& Other	Pediatric	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non- Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Musculo-skel. (Superfic.)	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
Cardiac	Cardiac Pediatric							
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

#### Additional Comments:

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

Table 4.3- 19 Diagnostic Ultrasound Indications for Use Form – MPT/7-4 Trans-Esophageal Echocardiography Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 series)						
<b>Transducer:</b>		MPT/7-4 7.0-4.0 MHz Trans-esophageal Echocardiography Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Abdominal organs and vascular)							
	Intra-operative (Neuro.)							
	Fetal Imaging & Other							
	Laparoscopic							
	Pediatric							
	Small Organ (breast, thyroid, testicles.)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
Cardiac	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (card.)	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Other (spec.)							
	Peripheral Vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

 Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K043559

Table 4.3- 20 Diagnostic Ultrasound Indications for Use Form – P/4-2 Phased Array Transducer

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 series)						
<b>Transducer:</b>		P/4-2 4.0-2.0 MHz Phased Array Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Abdominal	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging & Other	Laparoscopic							
	Pediatric	N	N	N		N	B+M; B+PWD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Neonatal Cephalic	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Adult Cephalic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
Cardiac	Cardiac Pediatric	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. This submission includes imaging to assist in the placement of needles and catheters in vascular or other anatomical structures.

Prescription Use (Per 21 CFR 801.109)

Indications for Use

*David A. Leggett*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number HOH3559  
 43 of 905 pages

Section 4.3

**Table 4.3- 21Table 4.3- 20 – Diagnostic Ultrasound Indications for Use Form – P/7-4 Phased Array Transducer Intended Use Form**

<b>System:</b>		SonoSite High-Resolution Ultrasound System (C3 series)						
<b>Transducer:</b>		P/7-4 7.0-4.0 MHz Phased Array Transducer						
<b>Intended Use:</b>		Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
Clinical Application		Mode of Operation						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
	Fetal	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Abdominal	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Intra-operative (Abdominal organs and vascular)	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Intra-operative (Neuro.)							
Fetal Imaging & Other	Laparoscopic							
	Pediatric	N	N	N		N	B+M; B+PWD	Note 1
	Small Organ (breast, thyroid, testicles.)	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Neonatal Cephalic	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Adult Cephalic	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Convent.)	N	N	N		N	B+M; B+PWD B+CD	Note 1
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
	Cardiac Adult	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
Cardiac	Cardiac Pediatric	N	N	N	N	N	B+M; B+PWD; B+CWD; B+CD	Note 1
	Trans-esophageal (card.)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	B+M; B+PWD; B+CD	Note 1
	Other (spec.)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Note 1: Other includes Color Power Doppler, combined B and Color Power Doppler, 3-D Imaging, Tissue Harmonic Imaging, Tissue Doppler Imaging and imaging for guidance of biopsy. Color Doppler includes Velocity Color Doppler. Color Doppler can be combined with any imaging mode. This submission includes imaging to assist in the placement of needles and catheters in vascular or other anatomical structures.

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number 44 of 905 pages

Indications for Use

Section 4.3