

SEP 23 2005

**PART B: 510(k) SUMMARY**

**Submitter:** Alliance Medical Corporation  
10232 South 51<sup>st</sup> Street  
Phoenix, Arizona 85044

**Contact:** Elizabeth Renken  
Regulatory Affairs Specialist  
(480) 763-5394 (o)  
(480) 763-5320 (f)  
[erenken@alliance-medical.com](mailto:erenken@alliance-medical.com)

**Date of preparation:** February 25, 2005

**Name of device:** Trade/Proprietary Name: Reprocessed Phacoemulsification Tips  
Classification Name: Phacofragmentation System

**Predicate Device**      **510(k) Title**      **Manufacturer**  
K911808      Alcon® Series 20,000® Legacy®      Alcon® Surgical, Inc.

**Device description:** Phacoemulsification Tips are used to emulsify and excise cataract tissue in ophthalmic microsurgical procedures. When connected to the ultrasonic handpiece of a phacoemulsification system and activated, the Phacoemulsification Tip vibrates at an ultrasonic frequency that emulsifies cataract tissue. The extracted tissue is then aspirated away through the hollow tip. Irrigation of the eye with a saline solution compensates for the loss of volume in the eye when the cataract tissue is removed.

**Intended use:** Reprocessed Phacoemulsification Tips are intended to emulsify and excise cataract tissues in ophthalmic microsurgical procedures.

**Indications statement:** Reprocessed Phacoemulsification Tips are indicated for use to emulsify and excise cataract tissues in patients requiring eye surgery.

**Technological characteristics:** The design, materials, and intended use of Reprocessed Phacoemulsification Tips are identical to the predicate devices. The mechanism of action of Reprocessed Phacoemulsification Tips is identical to the predicate devices in that the same standard mechanical design, materials, and sizes are utilized. There are no changes to the claims, intended use, clinical applications, patient population, performance specifications, or method of operation. In addition, Alliance Medical Corporation's reprocessing of Phacoemulsification Tips includes removal of adherent visible soil and decontamination. Each individual Phacoemulsification Tip is tested for appropriate function of its components prior to packaging and labeling operations.

**Performance data:** Bench and laboratory testing was conducted to demonstrate performance (safety and effectiveness) of Reprocessed Phacoemulsification Tips. This included the following tests:

- Biocompatibility
- Validation of reprocessing
- Sterilization Validation
- Function test(s)
- Packaging Validation

Performance testing demonstrates that Reprocessed Phacoemulsification Tips perform as originally intended.

**Conclusion:** Alliance Medical Corporation concludes that the modified devices (Reprocessed Phacoemulsification Tips) are safe, effective, and substantially equivalent to the predicate devices as described herein.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 23 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Alliance Medical, Inc.  
c/o Ms. Elizabeth Renken  
Regulatory Specialist  
10232 South 51<sup>st</sup> Street  
Phoenix, AZ 85044

Re: K050518  
Trade/Device Name: Alliance Reprocessed Phacofragmentation Needles  
Regulation Number: 21 CFR 886.4670  
Regulation Name: Phacofragmentation System  
Regulatory Class: Class II  
Product Code: NKX  
Dated: May 4, 2005  
Received: August 29, 2005

Dear Ms. Renken:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 827-8910. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in cursive script, reading "David M. Whipple".

David M. Whipple  
Acting Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

# ATTACHEMENT I

K050518

The Alcon Phacoemulsification Tip Models  
to be

Reprocessed by Alliance Medical Corporation:

MicroTip (0.9MM OD)	8065740478
MicroTip (0.9MM OD)	30RTS
MicroTip (0.9MM OD)	45RTS
MicroTip (0.9MM OD)	30KTS
MicroTip (0.9MM OD)	45KTS
ABS MicroTip (0.9MM OD)	8065790019
ABS MicroTip (0.9MM OD)	8065790020
ABS MicroTip (0.9MM OD)	8065790021
ABS MicroTip (0.9MM OD)	8065790022
ABS MicroTip (0.9MM OD)	8065790023
Flared ABS MicroTip (0.9MM OD)	8065740836
Flared ABS MicroTip (0.9MM OD)	8065740837
Flared ABS MicroTip (0.9MM OD)	8065740838
Flared ABS MicroTip (0.9MM OD)	8065740839
Flared ABS MicroTip (0.9MM OD)	8065740840
TurboSONICS Standard U/S Tip (1.1mm OD)	8065740476
TurboSONICS Standard U/S Tip (1.1mm OD)	15RT
TurboSONICS Standard U/S Tip (1.1mm OD)	30RT
TurboSONICS Standard U/S Tip (1.1mm OD)	45RT
TurboSONICS Standard U/S Tip (1.1mm OD)	30KT
TurboSONICS Standard U/S Tip (1.1mm OD)	45KT
TurboSONICS Standard ABS Tip (1.1mm OD)	8065740791
TurboSONICS Standard ABS Tip (1.1mm OD)	8065740792
TurboSONICS Standard ABS Tip (1.1mm OD)	8065740793
TurboSONICS Standard ABS Tip (1.1mm OD)	8065740794
TurboSONICS Standard ABS Tip (1.1mm OD)	8065740795

TurboSONICS Standard ABS Flared Tip (1.1mm OD)	8065740805
TurboSONICS Standard ABS Flared Tip (1.1mm OD)	8065740806
TurboSONICS Standard ABS Flared Tip (1.1mm OD)	8065740807
TurboSONICS Standard ABS Flared Tip (1.1mm OD)	8065740808
TurboSONICS Standard ABS Flared Tip (1.1mm OD)	8065740809

## 2. Indications for Use Statement

**510(k) Number (if known):**

**Device Name:** Reprocessed Phacoemulsification Tips

**Indications for Use:** Reprocessed Phacoemulsification Tips are indicated for use to emulsify and excise cataract tissues in patients requiring eye surgery.

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓  
(per 21 CFR 801.109)

or

Over-the-Counter Use \_\_\_\_\_

MB Nicholas

(Division Sign-Off)  
Division of Ophthalmic Ear,  
Nose and Throat Devices

510(k) Number K050518