

K050945

MAR 29 2006

510(k) Summary of Safety and Effectiveness

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

Introduction

According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

Submitter name, Address, and Contact

March 20, 2006

Lin-Zhi International, Inc.
687 North Pastoria Avenue
Sunnyvale, CA 94085
Phone: (408) 732-3856
Fax: (408) 732-3849

Contact: Marie Lin, Ph.D.
President, R&D Director

Device Name and Classification

Classification Name: The Cocaine and Cocaine Metabolite test system has been placed in Class II by the Bureau of Medical Devices.
Classification Number: DIO (21 CFR 862.3250)
Panel: 91Toxicology

The "Drug Specific, Calibrators" has been placed in Class II by the Bureau of Medical Devices.
Classification No.: DLJ, 21 CFR 862.3200
Panel: 91Toxicology

The "Single (Specified) Analyte Controls" has been placed in Class I by the Bureau of Medical Devices.
Classification No.: LAS, 21 CFR 862.3280
Panel: 91Toxicology

Common Name: Cocaine Oral Fluid Homogeneous Enzyme Immunoassays

Proprietary Name: None

Predicate Device(s)

The LZI Cocaine Metabolite (Benzoylecgonine) Oral Fluid Enzyme Immunoassay is substantially equivalent to the Cocaine Metabolite Intercept® Micro-plate EIA (K001197) manufactured by OraSure Technologies Inc. (formerly known as STC Technologies, Inc) for its general intended use.

Device Description

LZI's Oral Fluid Cocaine Metabolite (Benzoylecgonine) Enzyme Immunoassay is a ready-to-use, liquid reagent, homogeneous enzyme immunoassay. The assay uses specific antibodies that can detect benzoylecgonine in oral fluid with minimal cross-reactivity to various, common prescription drugs and abused drugs. The assay is based on competition between drug labeled with glucose-6-phosphate dehydrogenase (G6PDH) enzyme and free drug from the saliva sample for a fixed amount of specific antibody. In the absence of free drug from the saliva sample the specific antibody binds to the drug labeled G6PDH enzyme causing a decrease in enzyme activity. It is therefore the drug concentration is proportional to the enzyme activity. The G6PDH enzyme activity is determined spectrophotometrically at 340 nm by measuring its ability to convert nicotinamide adenine dinucleotide (NAD) to NADH.

Intended Use

The Cocaine Metabolite (Benzoylecgonine) Enzyme Immunoassays for Drugs of Abuse in Oral Fluid is a homogeneous enzyme immunoassay system to detect benzoylecgonine in human saliva with a cutoff of 15 ng/mL when testing oral fluid specimen collected with Salivette collector (manufactured by Sarstedt) and diluted with 1 mL of buffer. The calibrators and controls of the analyte (Benzoylecgonine) are prepared with oral fluid buffer so that it can be used to verify and validate the assay. The assay is intended for use in the qualitative determination for cocaine/cocaine metabolite drugs. The assay is designed for professional use with a number of automated clinical chemistry analyzers.

The Cocaine Metabolite (Benzoylecgonine) Oral Fluid Enzyme Immunoassay is a homogeneous enzyme immunoassay system provides only a preliminary analytical test result. A more specific alternative chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are used.

Comparison to Predicate Device

The LZI Cocaine Metabolite (Benzoylecgonine) Oral Fluid Homogeneous Enzyme Immunoassay, including calibrators and controls, is substantially equivalent to OraSure's Cocaine Metabolite Intercept® Micro-plate EIA in its intended use and in for the qualitative determination of cocaine/cocaine metabolite drugs in human oral fluid.

Device Characteristics	Subject Device (LZI Cocaine Metabolite (Benzoylecgonine) Oral Fluid Homogeneous EIA)	Predicate Device (OraSure Cocaine Metabolite Intercept® Micro-plate EIA)
Intended Use	<p>The LZI Cocaine Metabolite (Benzoylecgonine) Oral Fluid Homogeneous EIA is a homogeneous enzyme immunoassay system to detect cocaine metabolite in human saliva with a cutoff of 15 ng/mL when testing oral fluid specimen collected with Salivette collector (manufactured by Sarstedt) and diluted with 1 mL of buffer. The calibrators and controls of the analyte (Benzoylecgonine) are prepared with oral fluid buffer so that it can be used to verify and validate the assay. The assay is intended for use in the qualitative determination for cocaine/cocaine metabolite drugs.</p> <p><i>The Cocaine Metabolite (Benzoylecgonine) Oral Fluid Enzyme Immunoassay is a homogeneous enzyme immunoassay system provides only a preliminary analytical test result. A more specific alternative chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are used.</i></p>	<p>The OraSure Cocaine Metabolite Intercept® Micro-plate EIA is intended for use by clinical laboratories in the qualitative determination of cocaine and cocaine metabolite in oral fluid collected with Intercept® DOA Oral Specimen Collection Device using a 5.0 ng/mL cutoff. For <i>In Vitro</i> Diagnostic Use.</p> <p><i>The OraSure Cocaine Metabolite Intercept® Micro-plate EIA provides only a preliminary analytical test result. A more specific alternative chemical method should be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgement should be applied to any drugs of abuse test result, particularly when a preliminary, positive result is observed.</i></p>
Analyte	Benzoylecgonine	Benzoylecgonine
Matrix	Saliva	Saliva
Calibrators/ Controls Level	5 levels including a negative	4 levels including a negative

LZI Cocaine Metabolite (Benzoylcegonine) Oral Fluid Assay

Feature	Oral Fluid Cocaine Metabolite (Benzoylcegonine) EIA			
		Mean.	SD	% CV
Qualitative : (n=120) mA/min	Negative	240.0	1.71	0.71%
Within Run Precision:	5 ng/mL	255.1	2.00	0.78%
	10 ng/mL	270.3	2.17	0.80%
	20 ng/mL	285.8	2.32	0.81%
	50 ng/mL	321.9	2.49	0.77%
		Mean.	SD	% CV
Total Precision:	Negative	240.0	2.10	
	5 ng/mL	255.1	2.26	0.88%
				0.89%
	10 ng/mL	270.3	2.37	0.88%
	20 ng/mL	285.8	2.51	0.88%
	50 ng/mL	321.9	2.68	0.83%
Accuracy: Clinical patients samples (n=118) vs. GC/MS	95.8 % Agreement			
Specificity:	See attached Assay package insert			

OraSure Cocaine Metabolite Micro=Plate EIA

Feature	Benzoylcegonine	% CV	
Precision			
Intra-assay N=64	0 ng/mL	3.7	
	2.5 ng/mL	3.4	
	5.0 ng/mL	4.3	
	7.5 ng/mL	7.6	
Inter-assay N=4/day, 20 days	0 ng/mL	8.0	
	2.5 ng/mL	9.0	
	5.0 ng/mL	9.6	
	7.5 ng/mL	10.5	
Accuracy: Clinical patients samples (n=220) vs. GC/MS	93.2 % Agreement		
Specificity	See OraSure product insert		

Summary

The information provided in this pre-market notification demonstrates that the LZI Oral Fluid Cocaine Metabolite (Benzoylcegonine) Homogeneous EIA is substantially equivalent to the legally marketed predicated device for its general intended use. Data and results provided in this premarket notification were collected and prepared in accordance with the NCCLS guidance. Substantial equivalence was demonstrated through comparison of intended use and physical properties to the commercially available predicate device as confirmed by gas chromatography/mass spectrometry, an independent analytical method. The information supplied in this pre-market notification

provides reasonable assurance that the LZI Oral Fluid Cocaine Metabolite (Benzoylecgonine) Homogeneous EIA is safe and effective for its stated intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

MAR 29 2006

Ms. Marie Lin, Ph.D.
President, R&D Director
Lin-Zhi International, Inc.
687 North Pastoria Ave
Sunnyvale, CA 94085-2917

Re: k050945

Trade/Device Name: Oral Fluid Metabolite (Benzoylecgonine) Homogeneous Enzyme
Immunoassay, Calibrators and Controls

Regulation Number: 21 CFR 862.3250

Regulation Name: Cocaine and cocaine metabolite test system

Regulatory Class: Class II

Product Code: DIO, DLJ, LAS

Dated: March 08, 2006

Received: March 10, 2006

Dear Ms. Lin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Alberto Gutierrez', with a stylized flourish at the end.

Alberto Gutierrez, Ph.D.

Director

Division of Chemistry and Toxicology

Office of In Vitro Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): k050945

Device Name: Oral Fluid Cocaine Metabolite (Benzoylecgonine) Homogeneous Enzyme Immunoassay, Calibrators and Controls.

Indications For Use:

The Cocaine Metabolite (Benzoylecgonine) Enzyme Immunoassays for Drugs of Abuse in Oral Fluid is a homogeneous enzyme immunoassay system to detect cocaine metabolite in human saliva with a cutoff of 15 ng/mL when testing oral fluid specimen collected with Salivette collector (manufactured by Sarstedt) and diluted with 1 mL of buffer. The calibrators and controls of the analyte (Benzoylecgonine) are prepared with oral fluid buffer so that it can be used to verify and validate the assay. The assay is intended for use in the qualitative determination for cocaine/cocaine metabolite drugs. The assay is designed for professional use with a number of automated clinical chemistry analyzers.

The Cocaine Metabolite (Benzoylecgonine) Oral Fluid Enzyme Immunoassay is a homogeneous enzyme immunoassay system provides only a preliminary analytical test result. A more specific alternative chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are used.


Prescription Use ✓
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)


Division Sign-CT

Page 1 of 1

Office of In Vitro Diagnostic Device
Evaluation and Safety

k050945