

MAY - 2 2005

**3.1 Summary of Safety and Effectiveness****Non-Confidential Summary of Safety and Effectiveness**

Page 1 of 2

March 2, 2005

Medical Industries America, Inc.  
dba Evo Medical Solutions  
2636 – 289<sup>th</sup> Place  
Adel, IA 50003-8021

Tel: (515) 993-5001  
Fax: (515) 993-4172

**Official Contact:** Keith Theisen, Director of Engineering

**Proprietary or Trade Name:** RemRest heated humidifier

**Common/Usual Name:** Heated humidifier

**Classification Name:** Humidifier, Respiratory Gas (Direct Patient Interface)

**Predicate Devices:** Respironics – RemStar – K010263  
Fisher & Paykel – HC 150 – K003973

**Device Description:**

The RemRest Heated Humidifier is a simple heated humidifier to be used with a CPAP device to provide heated and humidified air to the patient. It operates by a conductive heating surface warming a water reservoir and air passes over the water. The device has the following basic specifications:

Size	13.1" L x 8.4" W x 4.7" H
Weight	2.5 lbs. dry
Water capacity	400 ml (1 2/3 cups)
Power	100-240 VAC 50/60Hz or 1.0 A max.
Heater setting	1-5 (85°F to 149°F / 29.4°C to 65°C)
Humidity Range	30 mg H <sub>2</sub> O/l @ 25 lpm
Pressure drop	0.5 cm H <sub>2</sub> O @ 40 lpm

**Indications:**

**Indications for Use --** The evo Heated Humidifier is intended to warm and add moisture to the breathing gases for administration to a patient. It is used for patients requiring Continuous Positive Airway Pressure (CPAP) therapy for the treatment of Obstructive Sleep Apnea. The addition of heated humidification to this therapy relieves the drying and irritating effects on the patient airways, which may arise from the use of a CPAP system. Use only with RemRest 900 series CPAP devices.

**Patient Population --** Patients utilizing CPAP devices for OSA and desire to have heated humidification

**Environment of Use --** Home, Hospital, and Institutional settings

**Non-Confidential Summary of Safety and Effectiveness**

Page 2 of 2  
March 2, 2005

**Comparison to Predicate Devices:**

	<b>RemRest Heated Humidifier</b>	<b>Predicates</b>
<b>Attributes</b>		
<b>Indications for use</b>	To warm and add moisture to the breathing gases for administration to a patient. It is used for patients requiring Continuous Positive Airway Pressure (CPAP) therapy for the treatment of Obstructive Sleep Apnea. The addition of heated humidification to this therapy relieves the drying and irritating effects on the patient airways, which may arise from the use of a CPAP system.	F&P – HC 150 – K003973  Respironics – RemStar with humidifier - K010263  Same
<b>Environments of use</b>	Home, Hospital, Institutional	Same
<b>Patient Population</b>	Patients utilizing CPAP devices for OSA and desire to have heated humidification	Same
<b>Contraindications</b>	None	Same
<b>Technology</b>		
Method of heating	Heated plate	Yes
Temperature controlled	Yes	Yes
Connects between CPAP and patient	Yes	Yes
Thermal cutoff	Yes	Yes

**Differences Between Other Legally Marketed Predicate Devices**

There are no significant differences between the proposed device, RemRest Heated Humidifier, and the identified predicates.

## Third Party Review Quality Assessment

### Section 1 – Submission Information

510(k) No.: K050990 Third Party Organization: Intetec / Intetec  
Third Party's Primary Reviewer(s): Will Devine  
ODE/OIVD Division: DAGAD Branch/Team: ARDB

### Section 2 – 510(k) Decision

Third party recommendation: SE ☒ NSE ☐ Other (specify): \_\_\_\_\_  
ODE/OIVD final decision: SE ☒ NSE ☐ Other (specify): \_\_\_\_\_

### Section 3 – Assessment of Third Party Review

Review Element	Rating (check one)		
	Adequate	Minor Issue(s)	Major Issue(s)
a. Determination of device eligibility for third party review	<input checked="" type="checkbox"/>		
b. Extent of pre-submission consultation with ODE/OIVD division	<input checked="" type="checkbox"/>		
c. Organization and format of review documentation	<input checked="" type="checkbox"/>		
d. Determination of 510(k) administrative completeness (screening review)	<input checked="" type="checkbox"/>		
e. Summary of device characteristics, intended use, and performance (including accessories, if applicable) and reason for 510(k) submission	<input checked="" type="checkbox"/>		
f. Comparison to legally marketed devices—identification and analysis of key similarities and differences	<input checked="" type="checkbox"/>		
g. Rationale for conclusions and recommendation	<input checked="" type="checkbox"/>		
h. Use of guidance documents and standards	<input checked="" type="checkbox"/>		
i. Resolution of 510(k) deficiencies and FDA requests for additional information	<input checked="" type="checkbox"/>		
j. Scope of reviewer expertise and use of consulting reviewers	<input checked="" type="checkbox"/>		
k. Other (specify):			

Comments (explanation of ratings/issues): OK.

### Section 4 – ODE/OIVD Assessor Information

Assessed by: [Signature] Date: 4/26/05 Tel. No.: 827 4479

**Routing:** Division--Clip completed assessment (this page only) to inside front cover of 510(k).  
DMC--Forward this page only to Eric Rechen, POS/ODE, Rm. 120J, Corp. Blvd. (HFZ-402).



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY - 2 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Medical Industries America, Incorporated  
C/O Mr. Ned E. Devine  
Responsible Third Party Official  
Intertek Testing Services NA, Incorporated  
70 Codman Hill Road,  
Boxborough, Massachusetts 01719

Re: K050990  
Trade/Device Name: RemRest Heated Humidifier  
Regulation Number: 868.5450  
Regulation Name: Respiratory Gas Humidifier  
Regulatory Class: II  
Product Code: BTT  
Dated: April 18, 2005  
Received: April 19, 2005

Dear Mr. Devine:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

3.3 Indications for Use

Page 1 of 1

510(k) Number: K050990 (To be assigned)

Device Name: RemRest Heated Humidifier

Indications for Use: The Evo Heated Humidifier is intended to warm and add moisture to the breathing gases for administration to a patient. It is used for patients requiring Continuous Positive Airway Pressure (CPAP) therapy for the treatment of Obstructive Sleep Apnea. The addition of heated humidification to this therapy relieves the drying and irritating effects on the patient airways, which may arise from the use of a CPAP system. Use only with RemRest 900 series CPAP devices.


It is portable and is intended to be used in the home, hospital, and institutional settings.

Prescription Use XX  
(Per CFR 801.109)

or

Over-the-counter use     

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
\_\_\_\_\_  
(Division Sign)  
Division of ~~Medical Devices~~ ~~Medical Devices~~ ~~Medical Devices~~  
Infection Control, Dental Devices

510(k) Number: K050990