

K051038

-X-

Y<sub>2</sub>

AUG 12 2005

15 March 2005

## 510K Summary

### Model 506 Patient Monitor

**Contact:** Alex Kaplan  
Director of QA & RA  
Criticare Systems, Inc.  
20925 Crossroads Circle  
Waukesha, WI 53186 USA  
262-798-8282 Voice  
262-798-8290 FAX

**Trade Name:** 506 Patient Monitor

**Common Name:** Vital Signs Monitor

**Classification Name:** Monitor, Physiological, Patient (74 MWT)

**Substantial Equivalence is claimed to :** 507EL Vital Signs Monitoring System (K022435).

#### Device Description:

The 506 monitor measures and displays real time physiological data of the patient, including a pulse bar plethysmogram and numerical data. The 506 can be used to monitor one or more of the following parameters: Noninvasive BP (NIBP), SpO<sub>2</sub>, and Temperature. For all these vital parameters, the 506 will be capable of limit alarms and alerts, printing of strip chart recordings and storing trends for retrospective review.

#### Intended Use:

This system is intended to monitor physiological parameters of patients within any healthcare environment. The user, responsible to interpret the monitored data made available, will be a professional health care provider. Physiological data, system alarms and patient data analysis will be available to the care provider from the monitor.

#### Comparison with predicate device:

Criticare Systems Inc. has developed and distributed physiological monitoring devices worldwide since its inception in 1984. The 506 monitor utilizes existing core technologies from the predicate 507EL monitor for patient monitoring of NIBP, SpO<sub>2</sub>,

-12-

2/2

and Temperature. Temperature measurement is provided by Kendall's Filac FasTemp module which has been integrated into the 506 monitor. This temperature measurement technology allows predictive or continuous temperature measurement on adult or pediatric patients in oral, axillar, or rectal placements. SpO2 can optionally be provided by Nellcor's SpO2 module. This SpO2 technology has been listed with the capability to be used in an environment where resistance to clinical motion is required. The patient data collected by the 506 monitor is displayed for the user on an alphanumeric display and LEDs equivalent to the predicate device. The 506 monitor utilizes alphanumeric display technology in combination with LED numeric displays. Membrane key panels provide a user interface equivalent to the predicate device. The packaging design of the 506 monitor is molded plastic and allows for it to be either a stationary monitor or to be used during patient translocation within the healthcare facility, as did the predicate 506EL.

#### **Determination of Substantial Equivalence:**

The 506 monitor performance for each monitoring modality has been confirmed to be equivalent to the predicate device. Additionally, the 506 complies with applicable safety and performance standards (detailed below) for each monitoring modality and verification of compliance has been completed. The patient monitoring technologies present in the 506 monitor have been in clinical use for at least six years in the 507EL monitor and its predicates. CSI's field experience with these modalities in the predicate devices has been satisfactory. This combination of equivalence testing, applicable objective standards compliance and field experience substantiates a high level of confidence in the safety and efficacy of the 506 monitor.

Therefore, the 506 monitor is substantially equivalent to the predicate devices.

#### **Compliance to standards and regulations:**

The 506 Vital Signs Monitor complies with the following national and international standards:

##### **Safety**

UL 60601-1 Medical Electrical Safety  
EN 60601-1-2 EMC Compliance  
IEC 60601-2-49 Multiparameter Monitor Safety  
ISO 10993-5,10-11 Biocompatibility

##### **Performance**

IEC 60601-2-30 NIBP Safety  
EN1060-1 NIBP Performance  
EN 1060-3 NIBP Performance {including EN 475 Alarm Performance}  
AAMI SP-10 NIBP Performance  
EN 865 Oximetry Performance (Equivalent to ASTM F 1415)



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

AUG 12 2005

Criticare Systems, Inc.  
c/o Mr. Alex Kaplan  
Director of QA and RA  
20925 Crossroads Circle  
Waukesha, WI 53186

Re: K051038  
Trade Name: 506 Patient Monitor  
Regulation Number: 21 CFR 870.2300  
Regulation Name: Cardiac Monitor (including cardiometer and rate alarm)  
Regulatory Class: Class II (two)  
Product Code: DRT  
Dated: August 5, 2005  
Received: August 9, 2005

Dear Mr. Kaplan:


We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4648. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "B. Zuckerman for".

Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications For Use

510(k) Number (if known): K051038

Device Name: 506 Patient Monitor

### Indications for Use:

The 506 Patient monitor interprets and displays real time physiological data of the patient including numerical data. The 506 is configured to monitor Noninvasive BP (NIBP),  $S_pO_2$ , and Temperature. For each patient vital parameter, the 506 will be capable of providing limit alarms and alerts, printing of strip chart recordings and storing data trends for retrospective review. The 506 monitor utilizes existing core technologies from the predicate 507EL monitor for patient monitoring of NIBP,  $SpO_2$ , and Temperature. Temperature measurement is provided by Kendall's Filac FasTemp module, which has been integrated into the 506 monitor. This temperature measurement technology allows predictive or continuous temperature measurement on adult or pediatric patients in oral, axillar, or rectal placements.  $SpO_2$  can optionally be provided by Nellcor's  $SpO_2$  module.


This system is intended to monitor physiological parameters of patients within any healthcare environment. The user, responsible to interpret the monitored data made available, will be a professional health care provider. Physiological data, system alarms and patient data analysis will be available to the care provider from the monitor.

Prescription Use ☒ \_\_\_\_\_  
(Part 21 CFR 801.109)

AND/OR

Over-The-Counter Use \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

  
(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number K051038