

**NOV 18 2005**

**510K SUMMARY**

**This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92**

**The assigned 510(k) number is: K051211**

**COMPANY/CONTACT PERSON**

Seradyn, Inc  
7998 Georgetown Road, Suite 1000  
Indianapolis, IN 46268

Establishment registration No: 1836010

Jack Rogers  
Manager of Regulatory Affairs  
Telephone: (317) 610-3823  
Fax: (317) 610-0018

**DATE PREPARED**

November 18, 2005

**DEVICE NAME**

Trade Name: QMS® Zonisamide  
Common Name: Homogeneous Particle-Enhanced Turbidimetric Immunoassay  
Device Classification: 21 CFR 862.3350 Diphenylhydantoin Test System; Class II

**INTENDED USE**

The QMS® Zonisamide assay is intended for the quantitative determination of zonisamide in human serum or plasma on automated clinical chemistry analyzers.

Zonisamide concentrations can be used as an aid in management of patients treated with zonisamide.

The QMS® Zonisamide Calibrator set is intended for use in calibration of the QMS Zonisamide assay.

The QMS® Zonisamide Control set is intended for use in quality control of the QMS Zonisamide assay.

**LEGALLY MARKETING DEVICE TO WHICH EQUIVALENCY IS CLAIMED**

Innofluor® Phenytoin (K955562)

**DESCRIPTION OF DEVICE**

The QMS® Zonisamide assay system is a homogeneous assay utilizing particle agglutination technology and is based on the competitive binding principle.

The assay consists of reagents R1: anti-Zonisamide rabbit polyclonal antibody and R2: Zonisamide-coated microparticles. A six-level set of QMS® Zonisamide Calibrators (A through F) is used to calibrate the assay. A three-level set of QMS® Zonisamide Controls (1 through 3) is used for quality control of the assay.

## COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

	<b>Device</b> QMS® Zonisamide	<b>Predicate</b> Innofluor® Phenytoin
Intended Use	The QMS Zonisamide assay is intended for the quantitative determination of zonisamide in human serum or plasma on automated clinical chemistry analyzers.	The Innofluor Phenytoin assay is intended for the quantitative determination of total phenytoin.
Indications for Use	Zonisamide concentrations can be used as an aid in management of patients treated with zonisamide.	Therapeutic drug monitoring.
Methodology	Homogeneous particle-enhanced turbidimetric immunoassay (particle agglutination)	Fluorescence Polarization Immunoassay (FPIA) technology.
Reagent Components	Two (2) reagent system: <ul style="list-style-type: none"> <li>• Anti-Zonisamide Antibody Reagent (R1) in buffers containing stabilizers with sodium azide</li> <li>• Zonisamide-coated Microparticle Reagent (R2) in buffer containing stabilizers with sodium azide</li> </ul>	Three (3) reagent system: <ul style="list-style-type: none"> <li>• A Phenytoin Antiserum (Sheep) in buffer with protein stabilizer and Sodium azide.</li> <li>• T Phenytoin Fluorescein Tracer in buffer with protein stabilizer, surfactant and Sodium azide</li> <li>• B Pre-Treatment Buffer with surfactant and Sodium azide</li> </ul>
Calibration	QMS Zonisamide Calibrators – six levels	Innofluor Phenytoin Calibrators – six levels
Quality Control	QMS Zonisamide Controls – three levels	N/A

## SUMMARY OF CLINICAL TESTING

### Accuracy and Linearity

Accuracy and linearity were determined by a study based on the NCCLS guideline *EP6: Evaluation of the Linearity of Quantitative Measurement*.

### Sensitivity

The Analytical Sensitivity or Least Detectable Dose (LDD) of the assay was determined to be 1.0 µg/mL.

The Functional Sensitivity or Limit of Quantitation (LOQ) of the assay was determined to be 3.0 µg/mL.

### Assay Range

Based on the Accuracy, Linearity, and Sensitivity (LDD and LOQ) data, the package insert claim for the reportable range for the assay is 3.0 to 50.0 µg/mL.

### Method Comparison

Correlation studies were conducted using NCCLS Guideline *EP9: Method Comparison and Bias Estimation Using Patient Samples* as a guideline to compare accuracy of recovery of Zonisamide in serum and plasma assayed by the QMS® Zonisamide assay to an HPLC reference method.

**Precision**

A precision study was performed using the National Committee for Clinical Laboratory Standards (NCCLS) guideline *EP5: Evaluation of Precision Performance of Clinical Chemistry Devices*.

**Specificity**

There are two known major metabolites of zonisamide: N-acetyl zonisamide (NAZ) and sulfamoyl acetyl phenol (SMAP). Testing results indicated that there is no significant cross-reactivity for either of these metabolites.

**Interferences**

Interference studies were conducted using NCCLS Guideline *EP7: Interference Testing in Clinical Chemistry*. The results of the study indicated that of 26 drugs tested, none showed cross-reactivity in the QMS zonisamide assay system.

**CONCLUSION**

As summarized above, the QMS<sup>®</sup> Zonisamide assay is substantially equivalent to the Innofluor<sup>®</sup> Phenytoin assay. Substantial equivalence has been demonstrated through performance testing to verify that the device functions as intended and that design specifications have been satisfied.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

NOV 18 2005

Mr. Jack Rogers  
Manager of Regulatory Affairs  
Seradyn, Inc.  
7998 Georgetown Road, Suite 1000  
Indianapolis, IN 46268

Re: k051211  
Trade/Device Name: QMS® Zonisamide Assay  
QMS® Zonisamide Calibrator Set  
QMS® Zonisamide Control Set  
Regulation Number: 21 CFR 862.3350  
Regulation Name: Diphenylhydantoin test system  
Regulatory Class: Class II  
Product Code: NWM, LAS, DLJ  
Dated: November 9, 2005  
Received: November 10, 2005

Dear Mr. Rogers:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Alberto Gutierrez', with a stylized flourish at the end.

Alberto Gutierrez, Ph.D.

Director

Division of Chemistry and Toxicology

Office of In Vitro Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K051211

Device Name: QMS® Zonisamide

### Indications for Use:

The QMS® Zonisamide assay is intended for the quantitative determination of zonisamide in human serum or plasma on automated clinical chemistry analyzers.

Zonisamide concentrations can be used as an aid in management of patients treated with zonisamide.

The QMS® Zonisamide Calibrator set is intended for use in calibration of the QMS Zonisamide assay.

The QMS® Zonisamide Control set is intended for use in quality control of the QMS Zonisamide assay.

Prescription Use   X   AND/OR Over-The-Counter Use             
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Signatory

Office of In Vitro Diagnostic  
Device Evaluation and Safety

510(k) K051211

Page 1 of   1