

K051388  
JUN 1 - 2005

## 510(k) Summary

The following 510(k) summary has been prepared pursuant to requirements specified in 21CFR1807.92(a).

### 807.92(a)(1)

#### **Submitter Information**

Carri Graham, Official Correspondent  
The Anson Group  
7992 Castleway Drive  
Indianapolis, Indiana 46250  
Phone: (317) 849-1916 x103  
Facsimile: (317) 577-9070

Contact Person: Carri Graham

Date: May 16, 2005

### 807.92(a)(2)

Trade Name: (6150) MyLab70 Ultrasound Imaging System

Common Name: Ultrasound Imaging System

Classification Name(s): Ultrasonic pulse doppler imaging system 892.1550  
Ultrasonic pulsed echo imaging system 892.1560

Classification Number: 90IYN; 90IYO

### 807.92(a)(3)

#### **Predicate Device(s)**

Esaote, S.p.A.	7250 Ultrasound Imaging System	K982444
Esaote, S.p.A.	7250 Ultrasound Imaging System	K994369
Esaote, S.p.A.	7350 Ultrasound Imaging System	K050326
Esaote, S.p.A.	Technos Ultrasound Imaging System	K014168
Esaote, S.p.A.	Technos Ultrasound Imaging System	K023255

**807.92 (a)(4)**

**Device Description**

The MyLab 70 is a mainframe ultrasound system used to perform diagnostic general ultrasound studies. Its primary modes of operation are: B-Mode, M-Mode, Doppler and Color Flow Mapping and, on lower frequency probes, Tissue Enhancement Imaging (TEI). The MyLab70 is equipped with a CRT Color Display. The full alphanumeric keyboard allows complete on-screen data entry of patient information and on-screen annotations.

The MyLab70 can drive phased (PA), convex (CA), linear array (LA) and Doppler probes.

The MyLab70 is equipped with a DVD-RW disk drive that can be used for image storage. Data can also be stored directly to a Personal Computer via a LAN port. Optional accessory devices available for the MyLab70 include a S-VHS video recorder and a monochrome or color page printer. The MyLab70 is equipped with an isolation transformer to adequately insulate the system's peripherals.

**807.92(a)(5)**

**Intended Use(s)**

Esaote's MyLab70 is a mainframe ultrasound system used to perform diagnostic general ultrasound studies including Cardiac, Transesophageal, Peripheral Vascular, Neonatal Cephalic, Small organ, Musculoskeletal (Conventional and Superficial), Abdominal, Fetal, Transvaginal, Transrectal, Adult Cephalic, Pediatric and Other: Urologic.

807.92(a)(6)

### Technological Characteristics

	<b>MyLab 70 this submission</b>	<b>Technos (K014168 &amp; K023255)</b>	<b>7350 MyLab 50 (K050326)</b>
Electrical Safety	IEC60601-1	IEC60601-1	IEC60601-1
Ultrasound Safety	Track 3 (Acoustic Output Display)	Track 3 (Acoustic Output Display)	Track 3 (Acoustic Output Display)
<b>Indication for Use</b>			
• OB/Fetal	YES	YES	YES
• Abdominal	YES	YES	YES
• Pediatric	YES	YES	YES
• Small organ	YES	YES	YES
• Neonatal Cephalic	YES	YES	YES
• Adult Cephalic	YES	YES	YES
• Cardiac	YES	YES	YES
• Transesophageal	YES	YES	YES
• Transrectal	YES	YES	YES
• Transvaginal	YES	YES	YES
• Peripheral Vascular	YES	YES	YES
• Musculoskeletal (conventional & superficial)	YES	YES	YES
<b>Probe Technology</b>			
• Phased Array	YES	YES	YES
• Linear Array	YES	YES	YES
• Convex Array	YES	YES	YES
• Doppler Probes	YES	YES	YES
• Bi-Scan	YES	NO	NO
<b>Modes of operation</b>			
2D, M-Mode, PW, CW, CFM, Amplitude Doppler (PD), TEI	YES	YES	YES
CnTI	YES	YES	No
TVM	YES	YES	YES
VPAN	YES	YES	NO
Imaging Frequencies	1 ÷ 16 MHz	1.5 ÷ 16 MHz	2 ÷ 10 MHz
CFM/Doppler Frequencies	2 ÷ 12 MHz	2 ÷ 12 MHz	2 ÷ 8 MHz
Tissue Velocity Mapping feature	YES	YES	YES
<b>Biopsy Guidance</b>			
• Biopsy Intended Uses	General Purpose, Transrectal, Transvaginal	General Purpose, Transrectal, Transvaginal	General Purpose, Transrectal, Transvaginal
Display type	CRT	CRT	CRT
Display Standard	SVGA	SVGA	SVGA
Digital Archival Capabilities	YES	YES	YES

	<b>MyLab 70 this submission</b>	<b>Technos (K014168 &amp; K023255)</b>	<b>7350 MyLab 50 (K050326)</b>
DICOM Classes:	Media Storage, Storage SCU	Media Storage, Storage SCU	Media Storage, Storage SCU
VCR / Page Printer	YES	YES	YES
M&A Capabilities	Cardiac, Vascular, OB and general purpose measurements	Cardiac, Vascular, OB and general purpose measurements	Cardiac, Vascular, OB and general purpose measurements
Weight	110 kg	140 kg	90 kg
Dimensions	60(w) x 160(h) x 90(d) cm	60(w) x 160(h) x 105(d) cm	60(w) x 155(h) x 90(d) cm



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUN 1 - 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Esaote, S.p.A.  
% Ms. Carri Graham  
Consultant  
The Anson Group  
7992 Castleway Drive  
INDIANAPOLIS IN 46250

Re: K051308

Trade Name: Model 6150 (MyLab70) Ultrasound Imaging System  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulatory Number: 21 CFR 892.1560  
Regulatory Name: Ultrasonic pulsed echo imaging system  
Regulation Number: 21 CFR 892.1570  
Regulation Name: Diagnostic ultrasonic transducer  
Regulatory Class: II  
Product Code: IYN, IYO, and ITX  
Dated: May 16, 2005  
Received: May 19, 2005

Dear Ms. Graham:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the Model 6150 (MyLab70) Ultrasound Imaging System, as described in your premarket notification:

Transducer Model Number

CA123  
CA421

CA430  
CA621

EC123  
LA424  
LA522  
LA523  
LA532  
PA121  
PA122

PA230  
TEE022  
TRT23  
2CW  
5CW  
BS230  
PA023

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850


This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled,

"Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,

*for* 

Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure(s)

Mod. 6150

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Abdominal		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Small Organ (specify)		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Neonatal Cephalic		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Adult Cephalic		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Cardiac		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Transesophageal		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Transrectal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Transvaginal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Laparoscopic										
Musculo-skeletal Conventional		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Musculo-skeletal Superficial		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Other (Urological)		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments:

Small Organs (thyroid, testicles, penis and breast);

Peripheral Vascular to include Vein Mapping & Sclerotherapy

Note (1): Combinations: any combination of the following modes: B+M+PW+CW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode  
Bi-Scan

*Ruth Phillips*  
(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number K051308

Prescription Use ✓



CA123

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Small Organ (specify)		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Neonatal Cephalic		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Adult Cephalic										
Cardiac		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments:

Small Organs (thyroid, testicles, penis and breast);

Peripheral Vascular to include Vein Mapping & Sclerotherapy

Note (1): Combinations: any combination of the following modes: B+M+PW+ CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Rachel Phillips*

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K051308

Prescription Use ☒

CA421

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Abdominal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Urological)		N	N	N		N	N		N (see Note 1)	N (see Note 2)

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW +CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Rachael Phillips*

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number 2051308

Prescription Use

## CA430

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Abdominal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW +CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Rae Phillips*

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*K057308*

Prescription Use ☒

## CA621

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Abdominal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Urological)		N	N	N		N	N		N (see Note 1)	N (see Note 2)

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Rita Phyllis*

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*K057302*

Prescription Use ☒

## EC123

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Transvaginal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Urological)		N	N	N		N	N		N (see Note 1)	N (see Note 2)

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW +CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode


*Rachel Phillips*

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

FDX(k) Number

1051308

Prescription Use 

**LA424**[illegible]

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

**Small Organs** (thyroid, testicles, penis and breast);

Peripheral Vascular to include Vein Mapping & Sclerotherapy

Note (1): Combinations: any combination of the following modes: B+M+PW +CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Rita Phillips  
(Division Sign Off)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number 2051308

Prescription Use

LA522

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N		N	N		N (see Note 1)	
Small Organ (specify)		N	N	N		N	N		N (see Note 1)	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N		N	N		N (see Note 1)	
Laparoscopic										
Musculo-skeletal Conventional		N	N	N		N	N		N (see Note 1)	
Musculo-skeletal Superficial		N	N	N		N	N		N (see Note 1)	
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments:

Small Organs (thyroid, testicles, penis and breast);

Peripheral Vascular to include Vein Mapping & Sclerotherapy

Note (1): Combinations: any combination of the following modes: B+M+PW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

*Racal*

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

1051308

Prescription Use

## LA523

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Small Organ (specify)		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Laparoscopic										
Musculo-skeletal Conventional		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Musculo-skeletal Superficial		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

**Additional Comments:**

Small Organs (thyroid, testicles, penis and breast);

Peripheral Vascular to include Vein Mapping & Sclerotherapy

Note (1): Combinations: any combination of the following modes: B+M+PW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Rae Phillips*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number 4051302

✓  
 Description Use



LA532

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Small Organ (specify)		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments:

Small Organs (thyroid, testicles, penis and breast);

Peripheral Vascular to include Vein Mapping & Sclerotherapy

Note (1): Combinations: any combination of the following modes: B+M+PW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Rosa P. Kelly*

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*K05,308*

PA121

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N	N		N	N		N (see Note 1)	N (see Note 2)
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Ra-A Phillips*

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K051300

PA122

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N	N	N	N		N (see Note 1)	
Small Organ (specify)										
Neonatal Cephalic		N	N	N	N	N	N		N (see Note 1)	
Adult Cephalic										
Cardiac		N	N	N	N	N	N		N (see Note 1)	
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N	N	N	N		N (see Note 1)	
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW+CW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

*Rosa Phillips*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K051308

*Prescription File*

PA230

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Cardiac		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW+CW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode

*Rachael Kelly*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number       K051308

TEE022

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N	N	N	N	N		N (see Note 1)	
Transesophageal		N	N	N	N	N	N		N (see Note 1)	
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW+CW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.


*Rachael Phillips*  
 Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Cardiovascular Devices  
 Division Number K051308


TRT23

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		N	N	N		N	N		N (see Note 1)	
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Urological)		N	N	N		N	N		N (see Note 1)	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW +CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Thoracic Devices  
 Device Number K051308

Prescription Use 

2 CW

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac					N					
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

*Rita Phyllis*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number 1052308

*[Signature]*

5 CW

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular					N					
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Urological)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

*Rachel Phillips*

Deputy Chief

Reproductive, Abdominal,  
Neurological Devices

File Number

*K051308*

*Prescription Only*



## BS230

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Abdominal		N	N	N	N	N	N		N (see Note 1)	N (see Note 2)
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW+CW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

Note (2): TEI (Tissue Enhanced Imaging) mode  
Bi-Scan

*Racal P...*

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*1051308*

PA023

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD (CW)	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric		N	N	N	N	N	N		N (see Note 1)	
Small Organ (specify)										
Neonatal Cephalic		N	N	N	N	N	N		N (see Note 1)	
Adult Cephalic										
Cardiac		N	N	N	N	N	N		N (see Note 1)	
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N	N	N	N	N		N (see Note 1)	
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Note (1): Combinations: any combination of the following modes: B+M+PW+CW+CFM+PD, where only one mode is live; B+CFM, M+CFM, B+PW, B+M, B+CFM+PW where all single modes are live.

*Ruth Pelly*

(Division Sign-Off)

Division of Reproductive, Abdominal,

Radiological Devices

File Number

*1057308*

*Division Chief* ✓