

K051686
p1/3

Abbreviated 510(k) Notification		getemed
CardioMem® CM 3000-12		Project ID: 0404H1
510(k) - Summary		Section 16-0001-Rev B

510(k) - Summary

Submitted By: getemed AG
Oderstr. 59
14513 Teltow
Germany
Tel.: +49 3328 – 3942-0
Fax: +49 3328 – 3942-99

Contact: Dr. Bert Schadow
Regulatory Affairs Manager

Manufacturing Facility: getemed Medizin- und Informationstechnik AG
Oderstr. 59
14513 Teltow
Germany

Date of Preparation: 2005-05-25

Trade Name: CardioMem® CM 3000-12 Holter recorder

Common Name: Holter recorder

Classification Name: Electrocardiograph, ambulatory (without Analysis)

Product Classification: 21 CFR 870.2800, Class II

Product Code: MWJ

Legally Marketed Devices: DR180+ Holter recorder (K004007, NorthEast Monitoring, Inc.)
H12+ Holter recorder (K021373, Mortara Instrument, Inc.)
CardioID+ (RZ 153+) Holter recorder (K022540, Rozinn Electronics, Inc.)

JUL 14 2005

Reason for Submission

Premarket notification (Abbreviated 510(k)) for CardioMem® CM 3000-12, a New Device, seeking authority to market the device under Section 510(k) as a device that is substantially equivalent to the DR180+ Holter recorder (K004007, NorthEast Monitoring, Inc.), the H12+ Holter recorder (K021373, Mortara Instrument, Inc.) and CardioID+ (RZ 153+) Holter recorder (K022540, Rozinn Electronics, Inc.).

Device Description

The CM 3000-12 is a Holter recorder designed to be used in conjunction with the evaluation software CardioDay® (not included in this 510(k)). This recorder is not capable of any diagnosis nor can it provide any interpretation of the data. The CM 3000-12 acquires, digitizes and stores data to be analyzed by CardioDay®. The CM 3000-12 utilizes a 10-lead electrode hookup and placement to provide CardioDay® with 12 channels of full disclosure for Holter analysis. The cardiac data provided by CardioDay® is used by trained medical personnel to assist in the diagnosis of patients with various rhythm patterns. The CM 3000-12 Holter recorder stores 12 ECG channels continuously for up to 48-hours. A keypad is available to set up system configuration, to enter patient's ID and name, to check lead quality during hook-up, and to start the recording. During the recording, the keypad can be used to enter event markers. The CM 3000-12 has a LCD screen to allow ECG display during the hook-up, lead

Abbreviated 510(k) Notification		getemed
CardioMem® CM 3000-12		Project ID: 0404H1
510(k) - Summary		Section 16-0001-Rev B

quality check, system configuration and various messages for the hook-up technician. The CM 3000-12 uses one or two AA batteries, and a removable memory card for data storage.

Intended use

The CardioMem® CM 3000-12 digital Holter recorder is intended to continuously record up to 48 hours of ECG data on a digital flash memory card. The CardioMem® CM 3000-12 performs no cardiac analysis by itself and is intended to be used with the analysis evaluation software CardioDay®. The recorded data are downloaded to a PC for analysis and following evaluation by a trained physician or health care professional.

Federal law restricts CardioMem® CM 3000-12 to use on order of a physician.

This device is available only upon the order of a physician or other licensed medical professional.


Indications for use

The CardioMem® CM 3000-12 is a Holter recorder which is indicated for patients who may benefit from a long-term continuous electrocardiographic (ECG) recording, including, but not limited to, those with complaints of palpitations, syncope, chest pain, shortness of breath, or those that need to be monitored to judge their current cardiac functionality.

Comparison of Technology Characteristics

The CardioMem® CM 3000-12 Holter recorder, the DR180+ Holter recorder, the H12+ Holter recorder, and the CardiolD+ (RZ153+) Holter recorder have the following technology specifications:

Specification	Legally Marketed Device DR180+ Holter recorder	Legally Marketed Device H12+ Holter recorder	Legally Marketed Device CardiolD+ (RZ153+)	New Device CardioMem CM 3000-12
Online data monitoring & alarm	No	No	No	No
Patient hookup	10 ECG electrodes	10 ECG electrodes	7 ECG electrodes	10 ECG electrodes
Number of ECG channels derived	12	12	3	12
ECG lead names	RA, LA, LL, RL, V1, V2, V3, V4, V5, V6	RA, LA, LL, RL, V1, V2, V3, V4, V5, V6	Holter configuration per AAMI EC38	RA, LA, LL, RL, V1, V2, V3, V4, V5, V6
A to D sample rate	720 samples/sec.	180 samples/sec.	1024 samples/sec.	1024 samples/sec.
A to D resolution	12 bit	20 bit	12 bit	12 bit
Pacemaker detection	Yes	Yes	Yes	No
Open-Lead detection	not specified	Yes	Yes	Yes
Memory type	CompactFlash™ Memory Card	CompactFlash™ Memory Card	CompactFlash™ Memory Card	CompactFlash™ Memory Card
Data transfer method	Via removable memory card	Via removable memory card	Via removable memory card	Via removable memory card
Memory card data format	Standard file system	Standard file system	Standard file system	Standard file system
Liquid crystal display (LCD)	Yes	Yes	Yes	Yes
Display purpose	Display ECG, check lead quality, input patient ID, display messages	Display ECG, check lead quality, input patient ID, display messages	Display ECG, check lead quality, input patient ID, display messages	Display ECG, check lead quality, input patient ID, display messages

Abbreviated 510(k) Notification		
CardioMem® CM 3000-12		Project ID: 0404H1
510(k) - Summary		Section 16-0001-Rev B

Specification	Legally Marketed Device DR180+ Holter recorder	Legally Marketed Device H12+ Holter recorder	Legally Marketed Device CardioID+ (RZ153+)	New Device CardioMem® CM 3000-12
Keyboard	Protected touch keys (membrane)	Protected touch keys (membrane)	Protected touch keys (membrane)	Protected touch keys (membrane)
System configuration	Check lead quality, input patient ID, recording mode, start recording	Check lead quality, input patient ID, start recording	Check lead quality, input patient ID and name, start recording	Check lead quality, input patient ID and name, start recording
System configuration method	Per keyboard and LCD display	Per keyboard and LCD display	Per keyboard and LCD display	Per keyboard and LCD display
Marker button	Yes	Yes	Yes	Yes
Size	125 x 70 x 25 mm	64 x 91 x 25 mm	108 x 79 x 22 mm	108 * 86 * 22 mm
Weight	142 g (without battery)	125 g (without battery)	<160 g (without battery)	<160 g (without battery)
Belt clip / Carrier pouch	not specified	Carrier bag (Pouch)	Carrier bag (Pouch)	Carrier bag (Pouch)
Battery	1 or 2 x 1.5 V AA alkaline	1 x 1.5 V AA alkaline	1 or 2 x 1.5 V AA alkaline	1 or 2 x 1.5 V AA alkaline
Battery check prior to recording	Yes	Yes	Yes	Yes
External patient cable	Yes	Yes	Yes	Yes
Record identification procedure	Yes	Yes	Yes	Yes
ECG channel preview	Yes	Yes	Yes	Yes

Conclusion

The CardioMem® CM 3000-12 Holter recorder, the DR180+ Holter recorder, the H12+ Holter recorder, and the CardioID+ (RZ153+) Holter recorder are intended to store continuously ECG Data for a period of at least 24 hours. The CM 3000-12, the DR180+, the H12+, and the RZ153+ do not analyze the data at the time of the recording. The recorded data of the mentioned recorders are downloaded to a PC for analysis and following evaluation by a trained physician or health care professional.

Quality system regulation 21 CFR 820 (outlined by the FDA) is a basis for the development of the CardioMem® CM 3000-12 Holter recorder. This recorder is safe and effective for the application for which it is intended and has been tested (environmental and safety tests, including EMC (Electromagnetic Compatibility) tests, according to the IEC 60601-1, IEC 60601-1-2, IEC 60601-2-25, IEC 60601-2-47 and ANSI/AAMI EC 38:1998 standards) to confirm the safety and efficacy of the recorder.

CardioMem® CM 3000-12 is found to be **substantially equivalent** to the DR180+ Holter recorder (K004007), the H12+ Holter recorder (K021373) and the CardioID+ (RZ153+) Holter recorder (K022540).



JUL 14 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Getemed Medizin-und Informationstechnik AG
c/o Ms. Carolann Kotula
MDI Consultants Inc.
55 Northern Blvd., Suite 200
Great Neck, NY 11021

Re: K051686
CardioMem® CM3000-12
Regulation Number: 21 CFR 870.2800
Regulation Name: Medical magnetic tape recorder
Regulatory Class: Class II (two)
Product Code: MWJ
Dated: June 20, 2005
Received: June 27, 2005

Dear Ms. Kotula:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0295. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in cursive script, appearing to read "B. Zuckerman for".

Bram D. Zuckerman, M.D.

Director

Division of Cardiovascular Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

Device Name: CardioMem® CM 3000-12

Indications For Use:

The CardioMem® CM 3000-12 is a Holter recorder which is indicated for patients who may benefit from a long-term continuous electrocardiographic (ECG) recording, including, but not limited to, those with complaints of palpitations, syncope, chest pain, shortness of breath, or those that need to be monitored to judge their current cardiac functionality.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Cardiovascular Devices

510(k) Number K051686

Page 1 of 1