

MAY - 5 2006

4/3

## Appendix 1: 510(k) Summary of Safety and Effectiveness

---

<b>Statement</b>	Information supporting claims of substantial equivalence, as defined under the Federal Food, Drug and Cosmetic Act, respecting safety and effectiveness is summarized below. For the convenience of the Reviewer, this summary is formatted in accordance with the Agency's final rule "...510(k) Summaries and 510(k) Statements..." (21 CFR §807) and can be used to provide a substantial equivalence summary to anyone requesting it from the Agency.
------------------	---

---

<b>Device description</b>	The Navigant™ Workstation with Niobe® Magnetic Navigation System [NWS05] is an interventional workstation for the navigation of appropriately equipped, magnetically adapted, devices (e.g., catheters or guidewires) through tissue to designated target sites. The system uses computer-controlled permanent magnets for orienting the tip of a magnetic device.
---------------------------	--

The system employs magnetic fields to *orient* or *steer* the tip of a magnetic device.

The NWS05 is a modification to the Navigant™ Workstation/Niobe® Magnetic Navigation System (K032937). The changes introduce new software design and mode of operation, new hardware, and a new range of motion to improve imaging, but maintain the existing technology for orientation of magnetically-adapted devices and clinical utility.

The NWS05 requires a digital fluoroscopy system to function properly.

---

<b>Intended use</b>	The Navigant™ Workstation with Niobe® Magnetic Navigation System is intended to navigate a compatible magnetic device through tissue to designated target sites in the right and left heart and coronary vasculature by orienting the device tip in a desired direction. The Navigant™ feature provides an enhanced navigation interface for the physician to control the MNS.
---------------------	--

---

<b>Substantial equivalence</b>	The NWS05 is substantially equivalent to the Navigant™ Workstation/Niobe® Magnetic Navigation System (K032937).
--------------------------------	---

---

*Continued on next page*

## Appendix 1: 510(k) Summary of Safety and Effectiveness, Continued

### Technological characteristics

The NWS05 employs magnetic fields to orient the distal tip of a magnetically-adapted device (e.g., catheter or guidewire).

### Device comparisons – steering control

The following is a comparison of the key features of The Navigant™ Workstation with Niobe® Magnetic Navigation System [NWS05] vs. the predicate device, the Navigant™ Workstation/Niobe® Magnetic Navigation System (NWS2; K032937).

Device Characteristics	New Device - NWS05	Predicate Device - NWS2
Intended use	To navigate a compatible magnetic device through tissue to designated target sites in the right and left heart and coronary vasculature by orienting the device tip in a desired direction.	To navigate a compatible magnetic device through tissue to designated target sites in the right and left heart and coronary vasculature by orienting the device tip in a desired direction.
Direct contact with patient tissue	No	No
Remote physician control of steerable device distal orientation	Yes	Yes
Computer control of steerable device distal orientation	Yes	Yes
Conducted under fluoroscopic visualization	Yes	Yes
Guided magnetic device employed	Specially designed magnetic catheters/guidewires	Specially designed magnetic catheters/guidewires
Steering control	Via magnetic fields, from a control room or at patient table side	Via magnetic fields, from a control room or at patient table side

Continued on next page

## Appendix 1: 510(k) Summary of Safety and Effectiveness, Continued

### Device comparisons – steering control (continued)

Device Characteristics	New Device - NWS05	Predicate Device - NWS2
System command	Physician-directed computer command	Physician-directed computer command
Magnetic field source	Two permanent magnets – positioned mechanically	Two permanent magnets – positioned mechanically
Operating field strength	Up to 0.10 T	Up to 0.10 T

**Physical testing** Performance testing has demonstrated substantial equivalence of the new device to the predicate device.

**Preclinical animal and clinical performance data** The Navigant™ Workstation with Niobe® Magnetic Navigation System [NWS05] is a modification of the predicate NWS2. Animal and clinical data are not necessary to support the modifications. Application data (animal and clinical) for magnetic navigation were provided and/or referenced in K032937.

**Contact** Peter A. Takes, Ph.D., RAC  
Director, Clinical & Regulatory Affairs  
Stereotaxis, Inc.  
4041 Forest Park Avenue  
St. Louis, Missouri 63108  
Ph. 314-615-6964  
Fax 314-615-6912  
takesp@stereotaxis.com

**Date** March 11, 2005



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAY - 5 2006

Stereotaxis, Inc.  
c/o Mr. Uwe Degenhardt  
TUV Product Service  
1775 Old Highway 8  
New Brighton, MN 55112-1891

Re: K051760  
Trade Name: Navigant™ Workstation with Niobe® Magnetic Navigation System,  
Version NWS05.  
Regulation Number: 21 CFR 870.1290  
Regulation Name: Steerable Catheter Control System  
Regulatory Class: Class II (two)  
Product Code: 74 DXX  
Dated: April 19, 2006  
Received: April 20, 2006

Dear Mr. Degenhardt:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Mr. Uwe Degenhardt

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Bram D. Zuckerman, M.D.

Director

Division of Cardiovascular Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K051760

Device Name: Navigant™ Workstation with Niobe® Magnetic Navigation System

Indications For Use: The Navigant™ Workstation with Niobe® Magnetic Navigation System is intended to navigate a compatible magnetic device through tissue to designated target sites in the right and left heart and coronary vasculature by orienting the device tip in a desired direction. The Navigant™ feature provides an enhanced navigation interface for the physician to control the MNS.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

B. Bimmuma  
(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number K051760

Page 1 of   1