

SEP 27 2005

K051865
1 of 6

**SUMMARY OF SAFETY AND EFFECTIVENESS INFORMATION
PERTAINING TO SUBSTANTIAL EQUIVALENCE**

A. Device Name

Proprietary Name

TERUMO® SurGuard2™ Safety Needle or similar proprietary name

Classification Name

Hypodermic Single Lumen Needle (880.5570) with antistick

Product Code: 80 FMI / 80 MEG

Classification: Class II

Common Name

Hypodermic needle with safety sheath or needle with needle protection device

B. Intended Use

The TERUMO® SurGuard2™ SAFETY NEEDLE device is intended for use in the aspiration and injection of fluids for medical purposes. The Terumo Safety Needle is compatible for use with standard luer slip and luer lock syringes.

Additionally, after withdrawal of the needle from the body, the attached needle safety shield can be manually activated to cover the needle immediately after use to minimize risk of accidental needlestick.

C. Device Description

The SurGuard2™ Safety Needle consists of a hypodermic needle with a hinged safety sheath attached to the connector hub with or without an attached hypodermic syringe. The safety sheath contains a locking mechanism which is activated when the sheath is manually pressed over the needle immediately after use and just prior to disposal to minimize the possibility of sharps injury. The safety sheath is activated with a one-handed operation by pressing the sheath against a firm surface. The locking mechanism is located at a designated position within the body of the short or long sheath appropriate for the needle size it is to contain. The needle gauge sizes fall within the 18 to 30gauge range and the needle lengths are 3/8" to 1 1/2". The hinge feature allows the user to set the sheath to the desired position for use. For user convenience, when the needle is in the "bevel up" position the sheath is located to the right. The SurGuard2™ Safety Needle will be individually packaged and sterilized by electron beam as a safety needle only or as a safety needle with attached Terumo syringe. The syringe sizes are 1, 3, 5, and 10cc/ml.

The specification of the needle and safety feature portions are identical to the SurGuard2™ Safety Needle cleared under K031453 and K040531, manufactured by Terumo Medical Corporation, Elkton, MD, USA. The attached 3, 5, and 10cc/ml syringes are identical to those syringe cleared under K023271, manufactured by Terumo Philippines Corporation. Information provided in this 510k will show that the 1cc/ml syringe manufactured by Terumo Philippines Corporation is substantially equivalent to the 1cc/ml manufactured by Terumo Medical Corporation, USA, cleared under K771205.

D. Substantial Equivalence

The SurGuard2™ Safety Needle manufactured by Terumo Philippines Corporation is substantially equivalent to:

1. K031453 TERUMO® SURGUARD2™ SAFETY NEEDLE (18 to 25g) with or without syringe manufactured by Terumo Medical Corporation, Elkton, Maryland.
2. K040531 TERUMO® SURGUARD2™ SAFETY NEEDLE (26 to 30g) with or without syringe manufactured by Terumo Medical Corporation, Elkton, Maryland.
3. K023271 TERUMO® SYRINGE manufactured by Terumo (Philippines) Corporation, Binan, Laguna, Philippines
4. K771205 TERUMO® HYPODERMIC SYRINGE (1cc/ml) manufactured by Terumo Medical Corporation, Elkton, Maryland

All of these cleared devices serve as predicates for the device which is subject of this 510k.

E. Principle of Operation and Technology

The Terumo SurGuard2™ Safety Needle device with and without syringe manufactured by Terumo (Philippines) Corporation and Terumo Medical Corporation, USA (K040531 and K031453) and all referenced predicate devices are operated manually.

F. Materials

The type of materials used for the hypodermic needle, safety feature and syringe of the SurGuard2™ Safety Needle device are identical to the materials used for the cleared SurGuard2™ Safety Needle cleared under K031453 and K040531. Any differences in materials between the SurGuard2™ Safety Needle manufactured by Terumo (Philippines) Corporation and the SurGuard2™ predicate manufactured by Terumo Medical Corporation, USA, raise no new issues of safety and effectiveness.

G. Specifications

Product Descriptions
18gauge x 1" (25mm) safety needle
18gauge x 1 1/2" (38mm) safety needle
19gauge x 1" (25mm) safety needle
19gauge x 1 1/2" (38mm) safety needle
20gauge x 1" (25mm) safety needle
20gauge x 1 1/2" (38mm) safety needle
21gauge x 1" (25mm) safety needle
21gauge x 1 1/2" (38mm) safety needle
22gauge x 1" (25mm) safety needle
22gauge x 1 1/2" (38mm) safety needle
23gauge x 1" (25mm) safety needle
23gauge x 1 1/2" (38mm) safety needle
25gauge x 5/8" (16mm) safety needle
25gauge x 1" (25mm) safety needle
25gauge x 1 1/2" (38mm) safety needle
26gauge x 1/2" (13mm) safety needle
27gauge x 1/2" (13mm) safety needle

30 gauge x 1/2" (13mm) safety needle
1cc/mL syringe with 25gauge x 5/8" (16mm) safety needle
1cc/mL syringe with 26gauge x 3/8" (9mm) safety needle
1cc/mL syringe with 27gauge x 1/2" (13mm) safety needle
3cc/mL syringe with 20gauge x 1" (25mm) safety needle
3cc/mL syringe with 20gauge x 1 1/2" (38mm) safety needle
3cc/mL syringe with 21gauge x 1" (25mm) safety needle
3cc/mL syringe with 21gauge x 1 1/2" (38mm) safety needle
3cc/mL syringe with 22gauge x 1" (25mm) safety needle
3cc/mL syringe with 22gauge x 1 1/2" (38mm) safety needle
3cc/mL syringe with 23gauge x 1" (25mm) safety needle
3cc/mL syringe with 25gauge x 5/8" (16mm) safety needle
3cc/mL syringe with 25gauge x 1" (25mm) safety needle
5cc/mL syringe with 20gauge x 1" (25mm) safety needle
5cc/mL syringe with 20gauge x 1 1/2" (38mm) safety needle
5cc/mL syringe with 21gauge x 1 1/2" (38mm) safety needle
10cc/mL syringe with 20gauge x 1" (25mm) safety needle
10cc/mL syringe with 20gauge x 1 1/2" (38mm) safety needle
10cc/mL syringe with 21gauge x 1" (25mm) safety needle

H. Performance

The following tests were performed on the SurGuard2™ Safety Needle:

- Activation Force
- Deactivation Force
- Puncture Resistance
- Sheath Removal Force
- Collar Removal Force
- Sheath Radial Force
- Protector Fit
- Adhesive Hold
- Simulated Use Study

The following tests were performed on the Terumo Philippines 1cc/ml syringe:

- Leakage (Aspiration and Injection)
- Plunger Gasket Fit
- Nominal Graduation Capacity

- Deadspace
- Conical Fitting

None of the data raises any new issues of safety and effectiveness. Additionally, a risk analysis was conducted and there were no new issues of safety and effectiveness.

The performance of the SurGuard2™ Safety Needle submitted in this 510k is substantially equivalent in intended use, design, technology/principles of operation, materials and performance to the legally marketed predicate devices.

I. Additional Safety Information

Manufacturing controls include visual, functional, and sterility tests.

The sterility of the device is assured using a sterilization method validated in accordance with ANSI/AAMI/ISO 11137-1994 Medical Devices – Validation and Routine Control of Radiation Sterilization. The SurGuard2™ Safety Needle is sterilized to provide a Sterility Assurance Level (SAL) of 10^{-6} .

The Terumo SurGuard2™ Safety Needle is classified as Externally Communicating Device, Blood Path Indirect, Limited Duration of Contact (< 24 hr). The device's blood contacting materials were tested in accordance with the tests recommended in the FDA General Program Memorandum #G95-1 (5/1/95): Use of International Standard ISO-10993, "Biological Evaluation of Medical Devices Part-1: Evaluation and Testing". Results of the testing demonstrate that the blood contacting materials are biocompatible.

J. Conclusion

The SurGuard2™ Safety Needle manufactured by Terumo Philippines Corporation is substantially equivalent to:

1. K031453 TERUMO® SURGUARD2™ SAFETY NEEDLE (18 to 25g) with or without syringe manufactured by Terumo Medical Corporation, Elkton, Maryland.
2. K040531 TERUMO® SURGUARD2™ SAFETY NEEDLE (26 to 30g) with or without syringe manufactured by Terumo Medical Corporation, Elkton, Maryland.
3. K023271 TERUMO® SYRINGE manufactured by Terumo (Philippines) Corporation, Binan, Laguna, Philippines
4. K771205 TERUMO® HYPODERMIC SYRINGE (1cc/ml) manufactured by Terumo Medical Corporation, Elkton, Maryland

Differences between the devices do not raise any significant issues of safety or effectiveness.

Terumo's statement of substantial equivalence is done solely to comply with the requirements of the Federal Food, Drug and Cosmetic Act and is not intended whatsoever to be the basis for a patent infringement action.

Date Prepared: July 8, 2005

Prepared By: Barbara Smith
Sr. Regulatory Affairs Specialist
Terumo Medical Corporation
Phone: 410-392-7241
Fax: 410-398-6079

Prepared For: Terumo (Philippines) Corporation
#124 East Main Avenue
Laguna Technopark
Binan, Laguna 4026
Philippines
Phone: 011-63-49-541-2111



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 27 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Barbara Smith
Senior Regulatory Affairs Specialist
Terumo Medical Corporation
Regulators Affairs Department
125 Blue Ball Road
Elkton, Maryland 21921

Re: K051865

Trade/Device Name: Terumo Surguard 2 Safety Needle or Similar
Regulation Number: 21 CFR 880.5570
Regulation Name: Hypodermic Single Lumen Needle
Regulatory Class: II
Product Code: MEG
Dated: July 6, 2005
Received: July 11, 2005

Dear Ms. Smith:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Chiu Lin', with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center of Devices and

Radiological Health

Enclosure

K051865
1 of 1

Indications For Use

510(k) Number (if known): _____

Device Name: _____ TERUMO® SurGuard2™ Safety Needle

Indications For Use:

The TERUMO® SurGuard 2™ SAFETY NEEDLE device is intended for use in the aspiration and injection of fluids for medical purposes. The Terumo Safety Needle is compatible for use with standard luer slip and luer lock syringes.

Additionally, after withdraw of the needle from the body, the attached needle safety shield can be manually activated to cover the needle immediately after use to minimize risk of accidental needlestick.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Anthony D. [Signature]
(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K051865

Page 1 of _____