


**BERKELEY ADVANCED BIOMATERIALS, INC.**

901 Grayson Street, Suite 101, Berkeley, CA 94710, USA

Tel: (510) 883 0500; Fax: (510) 883 0511

Email: [info@hydroxyapatite.com](mailto:info@hydroxyapatite.com)<http://www.hydroxyapatite.com>

K 051914

NOV - 3 2005

**510(K) Summary Statement**  
***GenerOs Bone Void Filler***

In accordance with the Food and Drug Administration Rule to implement provisions of the Safe Medical Devices Act of 1990 and in conformance with 21 CFR 807, this information serves as a Summary of Safety and Effectiveness for the use of GenerOs™ 80 and GenerOs 60 Bone Void Filler.

Submitted By:	Berkeley Advanced Biomaterials, Inc.
Date:	10 July 2005
Contact Person:	François Génin, Ph.D.
Position:	Chief Executive Officer
Contact Information:	Phone: 510-883-0500; Fax: 510-883-0511
Proprietary Name:	GenerOs™ 80, GenerOs™ 60
Regulation Name:	Resorbable Calcium Salt Bone Void Filler Device
Regulation Number:	888.3045
Classification:	Class II
Device Product Code/	Orthopedics/87/MQV
Panel Code:	

**DEVICE INFORMATION****A. INTENDED USES/INDICATIONS**

GenerOs (80 or 60) is an osteoconductive bone substitute shaped as granules or blocks that are intended to be used to fill voids and gaps that are not intrinsic to the stability of the bone structure. These gaps or voids may be located in the extremities, spine, pelvis, or cranium. The granules or blocks may be pressed into the void or into the surgical site by hand. The granules or blocks provide void filling material that acts as a temporary support medium. The granules or blocks are not intended to provide structural support during the healing process. The implant is radio-opaque. It is biocompatible and resorbs in the body as bone ingrowth occurs.

**B. DEVICE DESCRIPTION**

GenerOs (60 or 80) is a sterile osteoconductive bone void filler. This synthetic bone graft comes in the shape of granules or blocks. GenerOs (60 or 80) is supplied sterile for single patient use only. It is biocompatible and resorbs in the human body as bone ingrowth occurs when applied according to its indications for use. The implant is bioresorbable and radio-opaque.

**C. SUBSTANTIAL EQUIVALENCE INFORMATION**

The intended use, materials and design features of GenerOs 80 and GenerOs 60 are substantially equivalent to the predicate devices previously cleared for market. The safety

and effectiveness of GenerOs are adequately supported by the substantial equivalence information provided within the Premarket Notification.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

NOV - 3 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

François Génin, Ph.D.  
Chief Executive Officer  
Berkeley Advanced Biomaterials, Inc.  
901 Grayson St., Suite 101  
Berkeley, California 94710

Re: K051914 /S1  
Trade/Device Name: GenerOs™ 80, GenerOs™ 60  
Regulation Number: 21 CFR 888.3045  
Regulation Name: Resorbable calcium salt bone void filler device  
Regulatory Class: II  
Product Code: MQV  
Dated: September 19, 2005  
Received: October 6, 2005

Dear Dr. Génin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Mark N. Melkerson', is written over the typed name.

Mark N. Melkerson  
Acting Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K051914

Device Name: GenerOs Bone Void Filler

Indications for Use:

GenerOs™ is an osteoconductive bone substitute shaped as granules or blocks that are intended to be used to fill voids and gaps that are not intrinsic to the stability of the bone structure. These gaps or voids may be located in the extremities, spine, pelvis, or cranium.


The granules or blocks may be pressed into the void or into the surgical site by hand. The GenerOs™ granules or blocks provide void filling material that acts as a temporary support medium. The granules or blocks are not intended to provide structural support during the healing process. The implant is radio-opaque. GenerOs™ is biocompatible and resorbs in the body as bone ingrowth occurs.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

  
Concurrence of CDRH, Office of Device Evaluation (ODE)

**(Division Sign-Off)**  
**Division of General, Restorative,**  
**and Neurological Devices**

510(k) Number K051914