

DEC 23 2005

**SUMMARY OF SAFETY AND EFFECTIVENESS****Cardinal Health, Alaris Products®****VENTED VIAL ACCESS DEVICE****SUBMITTER INFORMATION**

- A. Company Name: Cardinal Health, Alaris Products
- B. Company Address: 10221 Wateridge Circle  
San Diego, CA 92121-2733
- C. Company Phone: (858) 458-7830  
Company Fax: (858) 458-6114
- D. Contact Person: Stacy L. Lewis  
Sr. Regulatory Affairs Specialist  
Cardinal Health, Alaris Products
- E. Date Summary Prepared: September 30, 2005

**DEVICE IDENTIFICATION**

- A. Generic Device Name: IV Fluid Transfer Pin
- B. Trade/Proprietary Name: SmartSite® Vented Vial Access Device
- C. Classification: Class II
- D. Product Code: LHI

**DEVICE DESCRIPTION**

The SmartSite® Vented Vial Access Device (Vented VAD) is a stand-alone, single use, multiple-dose device with a SmartSite® Needle Free Valve and a 0.2 micron hydrophobic air-venting filter for use with rubber-stopper vials. The device is intended for use when reconstituting lyophilized drugs and/or when dispensing medications, particularly hazardous drugs such as chemotherapeutics.

**SUMMARY OF SAFETY AND EFFECTIVENESS****Cardinal Health, Alaris Products®****SmartSite® Vented Vial Access Device****Page 2 of 2****SUBSTANTIAL EQUIVALENCE**

The SmartSite® Vented Vial Access Device is of comparable type and is substantially equivalent to the following predicate devices:

<b>Predicate Device</b>	<b>Manufacturer</b>	<b>510(k) No.</b>	<b>Date Cleared</b>
Chemo Dispensing Pin	B. Braun	K983794	3/23/99
SmartSite® Access Pin	Cardinal Health, Alaris Products*	K970485	4/21/97

\* Originally submitted by IVAC Holdings. Cardinal Health, Alaris Products was previously known as Alaris Medical Systems, which was the result of the merger of IVAC and IMED Corporation.

**INTENDED USE**

The SmartSite® Vented Vial Access Device is intended for use by healthcare professionals in a wide variety of healthcare environments including hospitals, healthcare facilities, and pharmacies when reconstituting lyophilized drugs and/or dispensing of medication, particularly hazardous drugs such as chemotherapeutics.

**TECHNOLOGICAL CHARACTERISTICS**

A comparison of the technological characteristics of the SmartSite® Vented Vial Access Device and the predicate devices has been performed. The results of this comparison demonstrate that the SmartSite® Vented Vial Access Device is equivalent to the marketed predicate devices in technological characteristics.

**PERFORMANCE DATA**

The performance data indicate that the SmartSite® Vented Vial Access Device meets specified requirements and is substantially equivalent to the predicate devices.



**DEPARTMENT OF HEALTH & HUMAN SERVICES**

Public Health Service

DEC 23 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Stacy L. Lewis  
Senior Regulatory Affairs Specialist  
Cardinal Health, Alaris® Products  
10221 Wateridge Circle  
San Diego, California 92121-2772

Re: K052790  
Trade/Device Name: SmartSite Vented Vial Access Device, Model TBD  
Regulation Number: 880.5440  
Regulation Name: Intravascular Administration Set  
Regulatory Class: II  
Product Code: LHI  
Dated: September 30, 2005  
Received: October 3, 2005

Dear Ms. Lewis:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

**INDICATIONS FOR USE**

510(k) Number: \_\_\_\_\_ (To Be Assigned By FDA)

Device Trade Name:

**Indications For Use:**

The SmartSite® Vented Vial Access Device is intended for use by healthcare professionals in a wide variety of healthcare environments including hospitals, healthcare facilities, and pharmacies for reconstitution or dispensing of medication. The SmartSite® Vented Vial Access Device is indicated for use with rubber-stopper medication vials for reconstitution or dispensing of medications, including chemotherapy agents.

Prescription Use   X   OR Over-The-Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109)

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

 for ADW 12/23/05

K052790