

NOV - 1 2005

## 510K SUMMARY

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92

The assigned 510(k) number is: K052815

### COMPANY/CONTACT PERSON

Seradyn, Inc  
7998 Georgetown Road, Suite 1000  
Indianapolis, IN 46268

Establishment registration No: 1836010

Jack Rogers  
Manager of Regulatory Affairs  
Telephone: (317) 610-3823  
Fax: (317) 610-0018

### DATE PREPARED

October 25, 2005

### DEVICE NAME

Trade Name: QMS® Amikacin  
Common Name: Homogeneous Particle-Enhanced Turbidimetric Immunoassay  
Device Classification: 21 CFR 862.3035; Amikacin Test System; Class II

### INTENDED USE

The QMS® Amikacin assay is intended for the quantitative determination of amikacin in human serum or plasma on automated clinical chemistry analyzers.

The results obtained are used in the diagnosis and treatment of amikacin overdose and in monitoring levels of amikacin to ensure appropriate therapy.

### LEGALLY MARKETING DEVICE TO WHICH EQUIVALENCY IS CLAIMED

Abbott TDx/TDxFLx Amikacin (K802669)

### DESCRIPTION OF DEVICE

The QMS® Amikacin assay system is a homogeneous assay utilizing particle agglutination technology and is based on the competitive binding principle.

In particle agglutination assays, the degree of agglutination is inversely proportional to the quantity of free drug in the reaction well. Hence, if no drug is present in the sample, the antibodies in the QMS® Amikacin Antibody Reagent (R1) will bind only to the bound drug on the particle which will cause it to agglutinate and will result in higher absorbance. If increased amount of competing drug is present in the sample, this will result in decreased binding of bound drug by the antibody, resulting in a relative decrease in particle agglutination. This in turn results in lower absorbance.

The precise relationship between particle agglutination and concentration of the unlabeled drug in the sample is established by measuring the absorbance values of calibrators with known concentration of the drug. The absorbance of unknown samples can be interpolated from the absorbance values of the calibration curve and the concentration of the drug present in the sample can be calculated.

The assay consists of reagents R1: anti-amikacin monoclonal antibody and R2: amikacin-coated microparticles. A six-level set of QMS® Amikacin Calibrators (A through F) is used to calibrate the assay.

## COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

	<b>Device</b> Seradyn QMS® Amikacin	<b>Predicate</b> Abbott TDx/TDxFLx Amikacin
Intended Use	The QMS Amikacin assay is for the quantitative determination of amikacin in human serum or plasma on automated clinical chemistry analyzers	The TDx/TDxFLx Amikacin assay is a reagent system for the quantitative measurement of amikacin, an aminoglycoside antibiotic drug, in serum or plasma.
Indications for Use	The results obtained are used in the diagnosis and treatment of amikacin overdose and in monitoring levels of amikacin to ensure appropriate therapy.	The measurements obtained are used in the diagnosis and treatment of amikacin overdose and in monitoring levels of amikacin to ensure appropriate therapy.
Methodology	Homogeneous particle-enhanced turbidimetric immunoassay (particle agglutination)	Fluorescence Polarization Immunoassay (FPIA) technology.
Reagent Components	Two (2) reagent system: <ul style="list-style-type: none"> <li>• Anti-amikacin Antibody Reagent (R1) in buffers containing stabilizers with sodium azide</li> <li>• Amikacin-coated Microparticle Reagent (R2) in buffer containing stabilizers with sodium azide</li> </ul>	Three (3) reagent system: <ul style="list-style-type: none"> <li>• Pretreatment Solution (P) Surfactant in buffer containing protein stabilizer and sodium azide.</li> <li>• S Amikacin Antiserum (Sheep) in buffer with protein stabilizer and Sodium azide.</li> <li>• T Amikacin Fluorescein Tracer in buffer with protein stabilizer, surfactant and Sodium azide</li> </ul>
Calibration	QMS Amikacin Calibrators – six levels	Amikacin Calibrators – six levels

## SUMMARY OF CLINICAL TESTING

### Accuracy

Accuracy by Recovery was determined by spiking USP traceable amikacin into human serum negative for the drug to achieve concentrations of 18.4 and 9.2µg/mL. The samples were analyzed in duplicate with the QMS Amikacin assay.

THEORETICAL CONC. (µG/ML)	Rep 1	Rep 2	Mean Recovered Conc.	SD	CV	% Recovery Acceptance Criteria: 100±10%
9.2	8.92	8.72	8.82	0.14	1.59%	95.87
18.4	16.79	17.13	16.96	0.17	1.00%	92.17
Mean Percent Recovery						94.02

### Linearity

Linearity by Dilution was determined by a study based on the NCCLS guideline *EP6: Evaluation of the Linearity of Quantitative Measurement*.

A linear regression analysis plot of USP Amikacin against recovered amikacin resulted in a line with a correlation coefficient ( $R^2$ ) of 0.9998, demonstrating that the assay is linear.

THEORETICAL CONC. ( $\mu\text{g/mL}$ )	Rep 1	Rep 2	Mean Recovered Conc.	SD	CV	% Recovery
1.5	1.77	1.57	1.67	0.10	5.99	111.33%
6.5	6.44	6.51	6.48	0.04	0.54	99.69%
15.0	14.84	14.49	14.67	0.18	1.19	97.80%
27.5	26.69	25.95	26.32	0.37	1.41	95.71%
42.5	41.24	41.63	41.44	0.20	0.47	97.51%
Mean Percent Recovery						100.41%

### Sensitivity

The Analytical Sensitivity or Least Detectable Dose (LDD) of the assay is defined as the concentration at which the lowest concentration is distinguishable from zero with 95% confidence.

The average LDD is 0.54  $\mu\text{g/mL}$ , supporting a claim of 0.8  $\mu\text{g/mL}$ .

### Assay Range

Based on the Accuracy, Linearity, and Sensitivity (LDD) data, the package insert claim for the reportable range for the assay will be 1.5 to 50  $\mu\text{g/mL}$ .

### Method Comparison

A study was conducted according to NCCLS Guideline *EP9: Method Comparison and Bias Estimation Using Patient Samples* to compare accuracy of recovery of amikacin in serum assayed by the QMS<sup>®</sup> Amikacin assay to the Abbott TDx/TDxFLx<sup>®</sup> Amikacin assay.

Mean values for the TDx reference method were plotted against those for the QMS on Hitachi 717. The results, using Passing – Bablok parameters, are:

N = 56

Slope = 1.00

y-intercept = 0.25

R = 0.996

$R^2$  = 0.992

Results show excellent correlation between the two assays.

## Precision

A precision study was performed using the National Committee for Clinical Laboratory Standards (NCCLS) guideline *EP5: Evaluation of Precision Performance of Clinical Chemistry Devices*.

	N	Mean µg/mL	Within Run		Between Day		Total	
			SD	CV (%)	SD	CV (%)	SD	CV (%)
Low Control	80	4.09	0.22	5.37	0.19	4.77	0.41	9.94
Mid Control	80	12.00	0.21	1.79	0.08	0.70	0.74	6.22
High Control	80	24.37	0.47	1.93	0.40	1.65	1.54	6.32

Acceptance Criteria: < 10% total CV

## Specificity

There are no metabolites of amikacin.

## Interferences

Interference studies were conducted using NCCLS Guideline *EP7: Interference Testing in Clinical Chemistry*.

### 1) Endogenous Substances

Interfering Substance	Interferent Concentration	N	Target (No Interferent) µg/mL	Mean Recovery µg/mL	% Recovery  Acceptance Criteria: 100±10%
Bilirubin	15mg/dL	2	21.65	20.87	96.40
Hemoglobin	10g/L	2	17.32	16.18	93.42
Triglyceride	1691 mg/dL	3	24.03	23.14	96.30
Total Protein	12 g/dL	3	24.03	23.07	96.00

### 2) HAMA

	Rep 1 µg/mL	Rep 2 µg/mL	Mean Recovery µg/mL	SD	CV	% Recovery  Acceptance Criteria: 100±10%
HAMA Type-1	20.68	20.13	20.41	0.27	1.35	100.5
Control	20.29	20.33	20.31	0.02	0.10	100.0
HAMA Type-2	16.97	16.99	16.98	0.01	0.05	98.04
Control	16.99	17.65	17.32	0.33	1.91	100.0

### 3) Common Co-Administered Drugs

Cross-reactant Drug	Conc. Tested µg/mL	Percent Cross- Reactivity
5-Fluorocytosine	30	-0.39
Amphotericin	100	1.33
Ampicillin	50	ND
Carbenicillin	2500	ND
Cephalexin	320	ND
Cephalosporin C	1000	ND
Cephalothin	1000	ND
Chloramphenicol	250	0.55
Clindamycin	2000	ND
Erythromycin	500	ND
Ethacrynic acid	400	ND
Furosemide	100	1.00
Fusidic acid	1000	ND
Gentamicin	100	ND
Kanamycin A	400	ND
Kanamycin B	400	ND
Lincomycin	2000	ND
Methicillin	200	0.41
Methotrexate	500	ND
Methylprednisolone	200	0.638
Neomycin	1000	ND
Netilmycin	125	ND
Oxytetracycline	2000	ND
Penicillin V	100	1.38
Prednisolone	12	2.36
Rifampicin	500	ND
Spectinomycin	100	ND
Streptomycin	400	ND
Sulfadiazine	1000	ND
Sulfamethoxazole	400	ND
Tetracycline	2000	ND
Tobramycin	100	0.32
Trimethoprim	200	ND
Vancomycin	400	ND

\*ND Not Detected

### 4) Anticoagulants

Studies were conducted to determine the performance characteristics of the assay for both serum and plasma samples containing amikacin.

The results indicate that there is no significant difference between the recovery of amikacin in serum or plasma. The collection tubes evaluated show no adverse effects on the recovery of amikacin, within the experimental error for the spiking study.

A claim for assay application to both serum and plasma samples is thus supported.

### On-Board Stability

#### 1) Calibration Curve stability

Calibration curve stability of a period of 28 days is supported by the data.

#### 2) Reagent On-Board Stability

A 40 day on-board reagent stability claim is supported by the data.

### CONCLUSION

As summarized above, the QMS® Amikacin assay is substantially equivalent to the Abbott TDx®/TDxFLx® Amikacin assay. Substantial equivalence has been demonstrated through performance testing to verify that the device functions as intended and that design specifications have been satisfied.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

NOV - 1 2005

Mr. Jack Rogers, MS, RAC  
Manger of Regulatory Affairs  
Seradyn, Inc.  
7998 Georgetown Road  
Suite 1000  
Indianapolis, IN 46268-5620

Re: k052815  
Trade/Device Name: QMS® Amikacin  
Regulation Number: 21 CFR 862.3035  
Regulation Name: Amikacin test system  
Regulatory Class: Class II  
Product Code: KLQ  
Dated: October 3, 2005  
Received: October 4, 2005

Dear Mr. Rogers

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 --

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Alberto Gutierrez", with a stylized flourish at the end.

Alberto Gutierrez, Ph.D.

Director

Division of Chemistry and Toxicology

Office of In Vitro Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K052815

Device Name: QMS® Amikacin

### Indications for Use:

The QMS® Amikacin assay is intended for the quantitative determination of amikacin in human serum or plasma on automated clinical chemistry analyzers.

The results obtained are used in the diagnosis and treatment of amikacin overdose and in monitoring levels of amikacin to ensure appropriate therapy.


Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

Page 1 of 1

510(k) K052815