

FEB 9 2006

K 060083

## 510(k) Summary of Safety and Effectiveness

**510(k) Submitter:** Streck  
7002 South 109<sup>th</sup> Street  
Omaha, NE 68128

**Official Correspondent:** Carol Thompson, Quality Assurance Manager  
(402)-537-5213

**Date Prepared:** January 10, 2006

**Name of Device:**

Trade Name: nRBC-Chex for LH  
Common Name: Assayed Hematology Control  
Classification Name: Nucleated Cell and Red Blood Cell Control (864.8625)

**Predicate Device:** Para 12 Plus Retics (K000945) Manufactured by Streck

**Description:**

nRBC-Chex for LH is stabilized suspension of human and animal blood, in a solution containing biological salts and anti-microbial preservatives. The product is packaged in plastic vials containing 4ml. The closures are polypropylene screw caps with polyethylene liners. There are two different levels. Level 1 has a low count and Level 2 has a higher count. The vials will be packaged in a six (6) or twelve (12) welled vacuum formed "clam-shell" container with the package insert / assay sheet. The product must be stored at 2 - 10°C.

**Intended Use:**

nRBC-Chex for LH is an assayed whole blood control designed to evaluate the accuracy and precision of Beckman Coulter® LH 750 in its measurement of the nucleated red blood cell parameter.

**Comparison to Predicate Device:**

Para 12 Plus Retics and nRBC-Chex for LH are both multi-parameter hematology control materials. They both contain RBC, WBC, and NRBC components. Unlike Para 12 Plus Retics, nRBC-Chex for LH does not include reticulocyte or platelet components. The suspending diluents are similar. Both Para 12 Plus Retics and nRBC-Chex for LH have closed vial stability performance claims of 75 days.

**Discussion of Tests and Test Results:**

Four types of studies were conducted to establish performance of nRBC-Chex for LH. The four tests conducted were Closed Vial Stability, Open Vial Stability, Run to Run Reproducibility, and Site to Site recovery of values. All testing showed that nRBC-Chex for LH is consistently reproducible, substantially equivalent to the predicate product and stable for the shelf life claimed.

**Conclusions Drawn From Tests:**

nRBC-Chex for LH is an effective quality control material for evaluating the accuracy and precision of the Beckman Coulter LH 750 in its measurement of the nucleated red blood cell parameter. It meets the claim of a 75 day closed vial, and a 14 day open vial stability and consistent run-to-run performance. Reproducibility studies and Closed Vial Stability results confirm lot-to-lot consistency in the manufacture of nRBC-Chex for LH. Customers can be assured of a reliable quality control material that meets their expectations.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Ms. Kerrie Oetter  
Quality Assurance Coordinator  
Streck Laboratories, Inc.  
7002 South 109<sup>th</sup> Street  
La Vista, NE 68128

FEB 9 2006

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Re: k060083

Trade/Device Name: nRBC-Chex for LH  
Regulation Number: 21 CFR 864.8625  
Regulation Name: Hematology quality control mixture  
Regulatory Class: Class II  
Product Code: GJR, GGL  
Dated: January 10, 2006  
Received: January 11, 2006

Dear Ms. Oetter:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 –

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,

A handwritten signature in cursive script, reading "Robert L. Becker, Jr.", written in dark ink.

Robert L. Becker, Jr., MD, Ph.D

Director

Division of Immunology and Hematology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K060083

Device Name: nRBC-Chex for LH

Indications For Use:

nRBC-Chex for LH is an assayed whole blood control designed to evaluate the accuracy and precision of the Beckman Coulter® LH 750 in its measurement of the nucleated red blood cell parameter.

Prescription Use   X    
(Per 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use         
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

*Josephine Bautista*  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

Page   1   of   1  

510(k) K060083