MAR 2 1 2006

510(k) Summary

In accordance with the Safe Medical Devices Act of 1990 and in compliance with 21CFR 807, the following serves as the 510(k) Summary information upon which the substantial equivalence determination is based.

Contact Information

Submitter:

BioTex, Inc.

8058 El Rio St.

Houston, TX 77054

Phone:

713,741.0111

Contact Person:

Matthew Fox

Date Prepared:

03/14/2006

Device Names

Trade/Proprietary Name:

PhoTex₁₅ Diode Laser Series: 980, 810, 940

Common Name:

Diode Laser Series

Classification Name:

Powered surgical laser instrument

Product Code:

GEX

Reg.Class:

П

Reg. Number:

878.4810

Predicate Device

SLT Thermalite Diode Laser Series: 980, 810, 940 (K952661)

Description of Device

The PhoTex₁₅ Diode Laser Series are diode lasers emitting radiation in either a continuous-wave (CW) or pulsed-mode in the infrared range at one of the following wavelengths: 980nm, 810nm, and 940nm. The PhoTex₁₅ Diode Laser Series provides a means for cutting, coagulation, and vaporization of tissue using a compatible fiber optic delivery accessory. The laser is compatible with any fiber optic delivery accessory terminated with a standard SMA905 connector whose core fiber diameter is 400 micron or larger with a numerical aperture of at least .37.

Indications for Use

The PhoTex15 Diode Laser Series is indicated for use in surgical applications requiring the ablation, vaporization, excision, incision, and coagulation of soft tissue in areas of surgery including: gastroenterology, general surgery, plastic surgery, genitourinary (urology), gynecology (GYN), neurosurgery, otolaryngology (ENT) head and neck, orthopedics, ophthalmology, pulmonology, and thoracic surgery.

Comparison to Predicate Device

The PhoTex₁₅ Diode Laser Series has been shown to be substantially equivalent to the predicate devices, the SLT Thermalite Diode Laser Series and PhotoMedex LaserPro Diode Laser Systems. Based on the technological features, device performance, and indications for use, BioTex, Inc. believes that no significant differences exist between the PhoTex₁₅ Diode Laser Series and the predicate device. Differences were determined to be minor and are each within the specifications listed by the predicate device and does not raise any concerns regarding the overall safety and effectiveness of the device.

Non-clinical Performance Tests:

Engineering studies have demonstrated the substantial equivalence of the PhoTex₁₅ Diode Laser Series to the SLT Thermalite Diode Laser Series (K 952661).

The studies concluded that the lasers are in compliance with FDA standards 21CFR1040.10 and 21CFR1040.11. In all instances, the lasers functioned as intended and performed in a manner similar to the predicate device when used in accordance with the labeled directions for use and specified indications.

Conclusion

BioTex has demonstrated the PhoTex₁₅ Diode Laser Series is substantially equivalent to the predicate device based on design, test results, and indications for use.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

MAR 2 1 2006

BioTex, Incorporated c/o Mr. Matthew Fox Development Engineer 8058 El Rio Street Houston, Texas 77054

Re: K060304

Trade/Device Name: PhoTex15 Diode Laser Series: 980, 810, 940

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser surgical instrument for use in general and plastic surgery and

in dermatology

Regulatory Class: II Product Code: GEX Dated: January 6, 2006 Received: February 8, 2006

Dear Mr. Fox:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments. Or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Mark N. Melkerson

Acting Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and

Radiological Health

Enclosure

Indications for Use

K060304

510(k) Number (if known):

| Device Name: | PhoTex15 Diode Laser Series: 980, 810, 940 |
|---|---|
| Indications For Use: | |
| | |
| | |
| | |
| ablation, vaporization, excising including: gastroenterology, | Series is indicated for use in surgical applications requiring the fon, incision, and coagulation of soft tissue in areas of surgery general surgery, plastic surgery, genitourinary (urology), rgery, otolaryngology (ENT) head and neck, orthopedics, and thoracic surgery. |
| | |
| Prescription Use X (Part 21 CFR 801 Subpart D) | AND/OR Over-The-Counter Use(21 CFR 801 Subpart C) |
| (PLEASE DO NOT WRITE NEEDED) | BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF |
| Concurrence of CDRH, Office of Device Evaluation (ODE) (Division Sign-Off) Division of General, Restorative, and Neurological Devices | |
| and redrological Devices | |
| 510(k) Number K060304 K060304 Response Letter Page 3 of 6 | |
| Tradeout response Letter Fage 5 01 0 | |