



K061033

DEC 12 2006

589 Davies Drive  
York, PA 17402  
phone 717 840-9335  
toll-free 800 221-1344  
fax 717 840-9347  
www.miltex.com

### 510(k) Summary [21 CFR §807.92]

Prepared: April 13, 2006

Device Trade Name: Miltex Rigid Sterilization Container System.

Device Common Name: Rigid Sterilization Container.

Classification Name: Sterilization wrap containers, trays, cassettes, and other accessories.

Class of Device: Class II device, product code KCT

Predicate Device: SteriTite® Rigid Sterilization Container System with MediTray  
Products- Case Medical, Incorporated- K023614

Official Contact: Charles Weaver, Regulatory Affairs Specialist

#### Device Description:

The Miltex Rigid Sterilization Container System consists of a family of rigid, re-usable, sealed containers that provides an effective sterilization packaging method for medical devices. Container bottoms and lids, within a given size, are interchangeable. The system is composed of the following components:

- Container bottoms (both perforated and non-perforated versions)
- Container baskets,
- Container lids (perforated only), and
- Container color-coding "labels."

The container system is designed for sterilant penetration through perforations in the lid and container bottom models that are perforated.

#### Intended Use:

The Miltex Rigid Sterilization Container System is intended for use in hospitals and health care facilities to contain other medical devices that are to be sterilized and to allow sterilization of the enclosed medical devices, including surfaces and lumens, using ethylene oxide gas sterilizers. Metal lumens 3-mm in diameter and up to 200-mm in length were validated. The container is intended for use with non-porous materials such as stainless steel surgical instruments.

Sterilized devices may be stored and transported within the container. The container is intended to maintain sterility of the medical devices within it until used.

Containers should not be stacked when used within ethylene oxide sterilizers.

Technological Characteristics:

A comparison of the technology characteristics of the Miltex Rigid Sterilization Containers to the predicate device's.

<b>Properties</b>	<b>Miltex System</b>	<b>SteriTite System</b>
Intended use to contain instruments being sterilized in ethylene oxide sterilizers	Yes	Yes
Intended to be re-used	Yes	Yes
Closed System	Yes	Yes
Sealed	Yes	Yes
<b>Design</b>		
Incorporates a filter system to permit entry of sterilant agent	Yes	Yes
Incorporates a filter system to prevent microbial migration during transport.	Yes	Yes
<b>Materials</b>		
Container	Aluminum alloy, Stainless Steel, & Silicone	Aluminum alloy, Stainless Steel, & Silicone

Performance Data:

A comparison of the non-clinical performance of the Miltex Rigid Sterilization Containers to the predicate device's.

Properties	Miltex System	SteriTite System
<b>Performance Standards</b>		
Conformance to appropriate AAMI standards	Yes, conforms to AAMI ST 77 Draft- <i>Containment Devices for Reusable Medical Device Sterilization</i>	Yes, conforms to AAMI ST 33- <i>Guidelines for the Selection and Use of Reusable Rigid Sterilization Container Systems for ETO Sterilization and Steam Sterilization in Health Care Facilities</i>
<b>Validation Testing</b>		
Ethylene Oxide Sterile Efficacy Testing	Yes	Yes
Aeration Time and EO Residuals	8-hours/ below maximum acceptable levels	12-hours/ below maximum acceptable levels
Sterility Maintenance	30-day real time w/ weekly handling events.	90-day real time w/ weekly handling events. 30-day real time w/ daily handling events.
Load	Up to 16-lbs. (small) Up to 20-lbs. (med.) Up to 25-lbs. (large)	Up to 22-lbs.
<b>Test Organisms/ Inoculated Product</b>		
Inoculated Lumens	Yes-- 3-mm I.D. x 200-mm, metal	Yes-- 2.2-mm I.D. x 457-mm, metal
Inoculated Stainless Steel Medical Devices	Yes (hinge area of medical pliers and knurled instruments)	Yes (blades)

Conclusion:

The Miltex Rigid Sterilization Containers is substantially equivalent to the SteriTite Container (K023614).



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Charles Weaver  
Regulatory Affairs Specialist  
Miltex, Incorporated  
589 Davies Drive  
York, Pennsylvania 17402

DEC 12 2006

Re: K061033

Trade/Device Name: Miltex Rigid Sterilization Container System  
Regulation Number: 880.6850  
Regulation Name: Sterilization Wrap  
Regulatory Class: II  
Product Code: KCT  
Dated: November 20, 2006  
Received: November 21, 2006

Dear Mr. Weaver:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Mr. Weaver

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin", with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

### Indications for Use

510(k) Number (if known): K061033

Device Name: Miltex Rigid Sterilization Container System

#### Indications for Use:

The Miltex Rigid Sterilization Container System is intended for use in hospitals and health care facilities to contain other medical devices that are to be sterilized and to allow sterilization of the enclosed medical devices, including surfaces and lumens, using ethylene oxide gas sterilizers. Metal lumens 3-mm in diameter and length up to 200-mm in length were validated. The container intended for use with non-porous materials such as stainless steel surgical instruments.

Sterilized devices may be stored and transported within the container. The container is intended to maintain sterility of the medical devices within it until used.

Containers should not be stacked when used within ethylene oxide sterilizers

The device models that are the subject of this pre-market notification are listed on page 2 of this Indications for Use statement.

Prescription Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use XX  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 2

Page 4-1, rev 2

*Shirley A. Murphy* (S)  
K 061033

Part Number	Description	UM
3-5100-10	STERI CONT 1/2 4 NON-PERFORATED BOTTOM, SILVER	EA
3-5100-13	STERI CONT 1/2 5 NON-PERFORATED BOTTOM, SILVER	EA
3-5100-15	STERI CONT 1/2 6 NON-PERFORATED BOTTOM, SILVER	EA
3-5110-10	STERI CONTAINER 1/2 4 PERFORATED BOTTOM, SILVER	EA
3-5110-13	STERI CONTAINER 1/2 5 PERFORATED BOTTOM, SILVER	EA
3-5110-15	STERI CONTAINER 1/2 6 PERFORATED BOTTOM, SILVER	EA
3-5172-05	STERI CONTAINER WIRE BASKET 1/2 2-INCHES	EA
3-5172-07	STERI CONTAINER WIRE BASKET 1/2 3-INCHES	EA
3-5172-10	STERI CONTAINER WIRE BASKET 1/2 4-INCHES	EA
3-5174-05	STERI CONTAINER PERFORATED BASKET 1/2 2-INCHES	EA
3-5174-07	STERI CONTAINER PERFORATED BASKET 1/2 3-INCHES	EA
3-5174-09	STERI CONTAINER PERFORATED BASKET 1/2 4-INCHES	EA
3-5300-10	STERI CONTAINER 3/4 4 NON-PERFORATED BOTTOM, SILVER	EA
3-5300-13	STERI CONTAINER 3/4 5 NON-PERFORATED BOTTOM, SILVER	EA
3-5300-15	STERI CONTAINER 3/4 6 NON-PERFORATED BOTTOM, SILVER	EA
3-5310-10	STERI CONTAINER 3/4 4 PERFORATED BOTTOM, SILVER	EA
3-5310-13	STERI CONTAINER 3/4 5 PERFORATED BOTTOM, SILVER	EA
3-5310-15	STERI CONTAINER 3/4 6 PERFORATED BOTTOM, SILVER	EA
3-5372-05	STERI CONTAINER WIRE BASKET 3/4 2-INCHES	EA
3-5372-07	STERI CONTAINER WIRE BASKET 3/4 3-INCHES	EA
3-5372-10	STERI CONTAINER WIRE BASKET 3/4 4-INCHES	EA
3-5374-05	STERI CONTAINER PERFORATED BASKET 3/4 2-INCHES	EA
3-5374-07	STERI CONTAINER PERFORATED BASKET 3/4 3-INCHES	EA
3-5374-09	STERI CONTAINER PERFORATED BASKET 3/4 4-INCHES	EA
3-5500-10	STERI CONTAINER 1/1 4 NON-PERFORATED BOTTOM, SILVER	EA
3-5500-13	STERI CONTAINER 1/1 5 NON-PERFORATED BOTTOM, SILVER	EA
3-5500-15	STERI CONTAINER 1/1 6 NON-PERFORATED BOTTOM, SILVER	EA
3-5510-10	STERI CONTAINER 1/1 4 PERFORATED BOTTOM, SILVER	EA
3-5510-13	STERI CONTAINER 1/1 5 PERFORATED BOTTOM, SILVER	EA
3-5510-15	STERI CONTAINER 1/1 6 PERFORATED BOTTOM, SILVER	EA
3-5572-05	STERI CONTAINER WIRE BASKET 1/1 2-INCHES	EA
3-5572-07	STERI CONTAINER WIRE BASKET 1/1 3-INCHES	EA
3-5572-10	STERI CONTAINER WIRE BASKET 1/1 4-INCHES	EA
3-5574-05	STERI CONTAINER PERFORATED BASKET 1/1 2-INCHES	EA
3-5574-07	STERI CONTAINER PERFORATED BASKET 1/1 3-INCHES	EA
3-5574-09	STERI CONTAINER PERFORATED BASKET 1/1 4-INCHES	EA
3-5940-01	STERI CONTAINER LABEL, SILVER	EA
3-5940-02	STERI CONTAINER LABEL, BLUE	EA
3-5940-03	STERI CONTAINER LABEL, RED	EA
3-5940-04	STERI CONTAINER LABEL, GOLD	EA
3-5940-05	STERI CONTAINER LABEL, GREEN	EA
3-5940-06	STERI CONTAINER LABEL, BLACK	EA
3-6161-00	STERI CONTAINER LID 1/2 PERFORATED, SILVER	EA
3-6361-00	STERI CONTAINER LID 3/4 PERFORATED, SILVER	EA
3-6561-00	STERI CONTAINER LID 1/1 PERFORATED, SILVER	EA