

K061703

AUG 02 2006

**510(k) Summary for the  
Dimension Vista™ System Chemistry 2 Calibrator  
(CHEM 2 CAL – KC120)**

**A. 510(k) Number:**

**B. Analytes:** Phosphorus (PHOS), Salicylate (SAL), and Triglycerides (TRIG).

**C. Type of Test:** Calibrator Material

**D. Applicant:** Dade Behring Inc., P.O. Box 6101, Newark, DE 19714-6101  
Victor M. Carrio, Regulatory Affairs and Compliance Manager  
Office: (302) 631-0376 Fax: (302) 631-6299

**E. Proprietary and Established Names:**

Dimension Vista™ System Chemistry 2 Calibrator  
(CHEM 2 CAL – KC120)

**F. Regulatory Information:**

1. Regulation section: 21 CFR § 862-1150 – Calibrator
2. Classification: Class II
3. Product Code: JIX – Calibrator, Multi-Analyte Mixture
4. Panel: Clinical Chemistry

**G. Intended Use:** The CHEM 2 CAL is an *in vitro* diagnostic product for the calibration of Phosphorus (PHOS), Salicylate (SAL), and Triglycerides (TRIG) methods on the Dimension Vista™ System.

**H. Device Description:**

CHEM 2 CAL is a liquid, multi-analyte, bovine serum albumin based product containing phosphorus, salicylate and glycerol. The kit consists of six vials, three vials of Calibrator A, and three vials of Calibrator B which are ready for use (no preparation is required). The volume per vial is 2.0 mL.

# **I. Substantial Equivalence Information:**

Item	Device		Predicates	
	Dimension Vista™ System Chemistry 2 Calibrator	Chemistry II Calibrator K861700	Salicylate Calibrator K904307	
Intended Use	The CHEM 2 CAL is an <i>in vitro</i> diagnostic product for the calibration of Phosphorus (PHOS), Salicylate (SAL), and Triglycerides (TRIG) methods on the Dimension Vista™ System.	CHEM II Calibrator is an <i>in vitro</i> diagnostic product to be used to calibrate the Dimension® clinical chemistry system for the magnesium (MG), phosphorus (PHOS) and triglycerides (TRIG) methods.	The Salicylate Calibrator is an <i>in vitro</i> diagnostic product to be used to calibrate the Dimension® clinical chemistry system for the Salicylate (SAL) method.	
Analytes	Phosphorus (PHOS), Salicylate (SAL), and Triglycerides (TRIG).	Phosphorus (PHOS) and Triglycerides (TRIG).	Salicylate (SAL).	
Form	Liquid.	Liquid.	Liquid.	
Traceability	PHOS – NIST SRM 2186I. SAL – Sodium Salicylate ACS grade. TRIG - Glycerol Anhydrous ACS grade.	PHOS - NIST SRM 2186I. TRIG - Glycerol Anhydrous ACS grade.	SAL - Sodium Salicylate ACS grade.	
Matrix	Bovine serum albumin based product.	Pure magnesium dissolved in a dilute solution of HCL, reagent grade potassium dihydrogen phosphate and reagent grade glycerol.	Aqueous solution of ACS grade sodium salicylate.	
Number of Levels	Two levels.	Three levels.	Three levels.	

#### J. Standard/Guidance Document Referenced:

1. Guidance: Guidance for Industry - Abbreviated 510(k) Submissions for In Vitro Diagnostic Calibrators; Final, 02/22/1999.  
Guidance for Industry and FDA Staff - Use of Symbols on Labels and in Labeling of In Vitro Diagnostic Devices Intended for Professional Use, 11/30/2004.
2. Standards: CEN 13640 Stability testing of In-Vitro Diagnostic Devices.  
ISO 14971:2000 Medical devices - Application of risk management to medical devices.

#### K. Performance Characteristics:

1. Stability: Target shelf life for the Dimension Vista™ Chemistry 2 Calibrator is 12 months. Calibrator shelf life is determined by comparing results of the product stored at 4°C with control stored at -20°C. The method is calibrated from this stored material. The 4°C material values are recovered versus the calibration. Recovery versus time is monitored and percent change over time is determined. Percent change should be less than or equal to 5%. Shelf-life stability (expiration) dating assignment at commercialization reflects the real-time data on file at Dade Behring, Inc.  
A vial punctured by the instrument and stored on board is stable for 7 days.  
An open vial not on instrument, but recapped and stored in a refrigerator is stable for 30 days.

For testing, vials are opened /punctured on day zero. A quantity sufficient for multiple calibrations is removed and the vials are recapped and stored at 2 – 8 °C. Opened/punctured vials are tested on days 1, 8, 15, 22, 31, and 32 versus freshly opened vials.

2. Traceability: The assigned values of the Chemistry 2 Calibrator are traceable to:

Constituent	Traceability
PHOS	NIST <sup>a</sup> SRM 2186I
SAL	Sodium Salicylate ACS <sup>b</sup> grade
TRIG	Glycerol Anhydrous (ACS grade)

- a. NIST-SRM: National Institute of Standards and Technology- Standard Reference Material.
- b. ACS: American Chemical Society.

### 3. Bottle Value Assignment:

The new calibrator is made by adding weighed in aqueous solutions of potassium dihydrogen phosphate, sodium salicylate, and glycerol reference materials to a stock solution at target concentrations. The concentration is verified using an instrument calibrated with primary standards. The stock solutions are added to base matrix in appropriate concentrations and verified using an instrument calibrated with primary standards. The final bottle assignment for each level of the commercial lot is tested  $N = 45$  replicates, with multiple reagent lots on multiple instruments.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Mr. Victor M. Carrio  
RA/QS Compliance Manager  
Dade Behring, Inc.  
500 GBC Drive  
PO Box 6101, Mail Stop 514  
Newark, DE 19714-6101

AUG 02 2006

Re: k061703  
Trade/Device Name: CHEM 2 CAL  
Regulation Number: 21 CFR§ 862.1150  
Regulation Name: Calibrator  
Regulatory Class: Class II  
Product Code: JIX  
Dated: June 15, 2006  
Received: June 16, 2006

Dear Mr. Carrio:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in cursive script, appearing to read "Alberto Gutierrez".

Alberto Gutierrez, Ph.D.

Director

Division of Chemistry and Toxicology

Office of In Vitro Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): k061703

Device Name: CHEM 2 CAL

Indications For Use:

The CHEM 2 CAL is an in vitro diagnostic product for the calibration of Phosphorus (PHOS), Salicylate (SAL), and Triglycerides (TRIG) methods on the Dimension Vista™ System.

Prescription Use  X   
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)



k061703