

WU'S TECH CO., LTD.

K062790

NO. 225, YUAN-PIER ST., HSIN CHU CITY, CHINA (TAIWAN)

TEL: 886-3-5382105 FAX: 886-3-5382191

Homepage: [www.wustech.com.tw](http://www.wustech.com.tw) Email: [wustis@ms45.hinet.net](mailto:wustis@ms45.hinet.net)

---

**" 510(k) SUMMARY "**

OCT 11 2006

Submitter's Name: **WU'S TECH CO., LTD.**

NO. 225, YUAN-PIER ST., HSIN CHU CITY, CHINA (TAIWAN)

Date summary prepared:

September 12, 2006

Device Name:

Proprietary Name: WU'S 3-WHEELED NEO SCOOTER WT-T3E

Common or Usual Name: Electrical Scooter

Classification Name: Motor Three-Wheeled Vehicle, Class II,  
21 CFR 890.3800

Indications for Use:

The device is intended for medical purposes to provide mobility to persons restricted to a seated position.

Description of the device:

The WU'S 3-WHEELED NEO SCOOTER WT-T3E is an indoor / outdoor electric scooter that is battery operated. It has a base with three-wheeled with a seat, armrests, and a front basket. The movement of the scooter is controlled by the rider who uses hand controls located at the top of the steering column. The device can be disassembled for transport and is provided with an onboard battery charger.

Performance Testing:

EMC Report ANSI / RESNA WC/Vol.2-1998, CISPR 11: 1990, EN61000-4-2: 1995, IEC61000-4-3: 1995 (Electrically powered wheelchairs, scooters, and their chargers – requirements and test methods)

Legally marketed device for substantial equivalence comparison:

WU'S 3-WHEELED NEO SCOOTER WT-T3D (K032488)

# WU'S TECH CO., LTD.

NO. 225, YUAN-PIER ST., HSIN CHU CITY, CHINA (TAIWAN)

TEL: 886-3-5382105 FAX: 886-3-5382191

Homepage: [www.wustech.com.tw](http://www.wustech.com.tw) Email: [wustis@ms45.hinet.net](mailto:wustis@ms45.hinet.net)

---

## Summary for substantial equivalence comparison:

The intended uses, weight limit, cruising range, maximum speed, back upholstery, safety climbing angle, and warranty period between the new device WT-T3E and the predicate device WT-T3D are all the same. Especially the electronic systems between two devices are all passed by the UL certificated, for instance the electronic controller, batteries and recharge. Besides, the back upholstery is the same material, and also passed the resistance ignition test by SGS. Thus the same safety level for the two devices is assured.

The major difference existing for new device is the overall dimension, the size of tires, and the weight are differences between the two devices. The overall appearance differences are not safety aspect. Thus the new device is substantially equivalent to the predicate devices in this aspect.

# WU'S TECH CO., LTD.

NO. 225, YUAN-PIER ST., HSIN CHU CITY, CHINA (TAIWAN)

TEL: 886-5-5382105 FAX: 886-5-5382191

Homepage: [www.wustech.com.tw](http://www.wustech.com.tw) Email: [wustis@ms45.hinet.net](mailto:wustis@ms45.hinet.net)

---

## C.2 COMPARISON SUMMARY

*( We place the related information for the predicate device in the following pages. )*

The intended uses, weight limit, cruising range, maximum speed, back upholstery, safety climbing angle, and warranty period between the new device WT-T3E and the predicate device WT-T3D are all the same. Especially the electronic systems between two devices are all passed by the UL certificated, for instance the electronic controller, batteries and recharge.. Besides, the back upholstery is the same material, and also passed the resistance ignition test by SGS. Thus the same safety level for the two devices is assured.

**The major difference existing for new device is the overall dimension, the size of tires, and the weight are differences between the two devices. The overall appearance differences are not safety aspect. Thus the new device is substantially equivalent to the predicate devices in this aspect.**

*Based on the above the information and the analysis, we know that the subject device and the predicate devices have the same intended use, the same technological aspects and only minor dimensions or data differences exist. We believe that FDA can decide the subject device and the predicate device are substantially equivalent.*



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

OCT 11 2006

Wu's Tech Co., Ltd.  
% Dr. Jen Ke-Min  
ROC Chinese-European Industrial Research Society  
No. 225, Yuan-Pier Street,  
Hsin Chu City,  
China (Taiwan)

Re: K062790

Trade/Device Name: WU's 3-Wheeled Neo Scooter, WT-T3E  
Regulation Number: 21 CFR 890.3800  
Regulation Name: Motorized three-wheeled vehicle  
Regulatory Class: Class II  
Product Code: INI  
Dated: September 12, 2006  
Received: September 18, 2006

Dear Dr. Ke-Min:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Dr. Jen Ke-Min

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120 . Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", written in a cursive style.

Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510 (K) Number ( If Known ):   K  

Device Name: WU'S 3-WHEELED NEO SCOOTER, WT-T3E

Indications for Use:

The device is intended for medical purposes to provide mobility to persons restricted to a sitting position.

Prescription Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use   √  

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Barbara Bruchmann  
(Division Sign-Off)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Division of General, Restorative,  
and Neurological Devices

510(k) Number K062790

Page   1   of   1