

K063695

---

## SECTION II. 510(K) SUMMARY

APR 27 2007

### **A. Device Name**

|                     |                       |
|---------------------|-----------------------|
| Proprietary Name    | Runthrough NS         |
| Classification Name | Wire, Guide, Catheter |
| Common Name         | Guide Wire            |

### **B. Intended Use**

The Runthrough NS are used to facilitate placement of balloon dilatation catheters for percutaneous transluminal coronary angioplasty (PTCA) and/or percutaneous transluminal angioplasty (PTA).

Note: This is the same intended use as the predicate device – Radifocus® Glidewire for Coronary Use with Platinum or Gold Coil, K961445.

### **C. Device Description**

The Runthrough NS is a coil-type guide wire. The main components of the wire include a core wire, a tip coil, a tip coil marker, and surface coating. The core wire is constructed of a Ni/Ti alloy wire and a stainless steel wire joined together. The tip coil marker, a Pt/Ni alloy, is radiopaque. The tip coil has lubricous coating (silicone coating) and/or hydrophilic coating on the surface depending on product code. The shaft is surface-coated with silicone and PTFE. Tip flexibility can be selected among three types, Extra Floppy, Floppy, and Intermediate, flexibility decreasing in the named order. The wire is also available in a hyper-coating type which is more lubricous than the standard wire. The device may be accompanied by a torque device, inserter, and stylet accessories.

### **D. Principle Of Operation / Technology**

The Runthrough NS is operated manually or by a manual process.

00014

***E. Design / Materials***

Differences in materials between the modified device and the predicate device the Radifocus® Glidewire for Coronary Use with Platinum or Gold Coil cleared under K961445 raise no new issues of safety and effectiveness.

***F. Specifications***

| Feature            | Modified Runthrough NS             | Unmodified Device<br>Radifocus® Glidewire<br>for Coronary Use with<br>Platinum or Gold Coil<br>cleared under K961445 |
|--------------------|------------------------------------|--|
| Available diameter | 0.014" ( 0.36mm )                  | 0.014" ( 0.36mm )  |
| Available length   | 180cm, 300cm                       | 180cm, 300cm   |
| Tip marker length  | 30mm                               | 20mm   |
| Accessory Devices  | Torque Device, Inserter,<br>Stylet | Torque Device, Inserter  |
| Shelf life         | 36 months                          | 24 months  |

### ***G. Performance***

The following verification tests were performed to demonstrate the substantial equivalence of the modified device (Runthrough NS) to the unmodified device (Radifocus® Glidewire for Coronary Use with Platinum or Gold Coil, K961445).

- Dimensional Inspection
- Tip sliding resistance test
- Tensile strength of tip
- Tip butting load
- Tensile strength of shaft
- Wire stability in holder loop
- Bend strength
- Seal strength of packaging
- Reshapability
- Torque Strength
- Torqueability
- Tip Flexibility
- Coating Adherence/Integrity

None of the data raises any new issues of safety and effectiveness. Additionally, a risk analysis was conducted and there were no new issues of safety and effectiveness.

Therefore the performance of the modified Runthrough NS is substantially equivalent to the performance of the predicate device the Radifocus® Glidewire for Coronary Use with Platinum or Gold Coil, K961445.

### ***H. Additional Safety Information***

Manufacturing controls include visual, functional, dimensional and sterility tests.

Blood contacting materials were tested in accordance with the tests recommended in the FDA General Program Memorandum #G95-1 (5/1/95): Use of International Standard ISO-10993, “Biological Evaluation of Medical Devices Part-1: Evaluation and Testing”.

The guide wire is classified as Externally Communicating Devices, Circulating Blood, Limited Contact ( $\leq 24$  hrs). Results of the testing demonstrate that the blood contacting materials are biocompatible.

Sterilization conditions have been validated in accordance with ANSI / AAMI / ISO 11135-1994, *Medical Devices – Validation and routine control of ethylene oxide sterilization* and EN 550. The device is sterilized to a SAL of  $10^{-6}$ .

**H. Substantial Equivalence**

The modified Runthrough NS is substantially equivalent in intended use, design, technology/principles of operation, materials and performance to the predicate device the Radifocus® Glidewire for Coronary Use with Platinum or Gold Coil, K961445. Differences between the two devices do not raise any significant issues of safety or effectiveness.

**I. Submitter Information**

Prepared By: Mr. Mark Unterreiner  
Sr. Regulatory Affairs Specialist

Prepared For: Terumo Medical Corporation  
950 Elkton Blvd.  
Elkton, MD 21921  
Phone: (410) 392-7213  
Fax: (410) 398-6079  
Email: mark.unterreiner@terumomedical.com

Date Prepared: December 12, 2006



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

APR 27 2007

Terumo Medical Corporation  
c/o Mr. Mark Unterreiner  
Sr. Regulatory Affairs Specialist  
950 Elkton Blvd.  
Elkton, MD 21921

Re: K063695  
Runthrough NS  
Regulation Number: 21 CFR 870.1330  
Regulation Name: Catheter guidewire  
Regulatory Class: II (two)  
Product Code: DQX  
Dated: April 4, 2007  
Received: April 5, 2007

Dear Mr. Unterreiner:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

*Bram D. Zuckerman*

*BZ*

Bram D. Zuckerman, M.D.

Director

Division of Cardiovascular Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K063695

Device Name: Runthrough NS

Indications For Use:

The Runthrough NS are used to facilitate placement of balloon dilatation catheters for percutaneous transluminal coronary angioplasty (PTCA) and/or percutaneous transluminal angioplasty (PTA).

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use             
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Diana R. Vachner  
Division Sign-Off)  
Division of Cardiovascular Devices

510(k) Number K063695

00005