

SEP 11 2007

**510(k) Summary for
Dimension Vista® Protein 1 Calibrator and
Dimension Vista® Protein 1 Control L, M and H**

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is:

K071980

1. Manufacturer's Name, Address, Telephone, and Contact Person, Date of Preparation:

Manufacturer: Dade Behring Marburg GmbH
Emil-von-Behring Str. 76
35041 Marburg, Germany

Contact Information: Dade Behring Inc.
P.O. Box 6101
Newark, Delaware 19714-6101
Attn: A. Kathleen Ennis
Tel: 302-631-9352
Fax: 302-631-6299

Preparation date: June 8, 2007

2. Device Name: Dimension Vista® Protein 1 Calibrator
Dimension Vista® Protein 1 Control L
Dimension Vista® Protein 1 Control M
Dimension Vista® Protein 1 Control H

Classification: Class II; Class I
Product Code: JIX; JJY
Panel: Clinical Chemistry (75)

3. Identification of the Legally Marketed Device:

Dimension Vista® Protein 1 Calibrator	K063663
Dimension Vista® Protein 1 Control L	K063663
Dimension Vista® Protein 1 Control M	K063663
Dimension Vista® Protein 1 Control H	K063663

4. Device Descriptions:

Dimension Vista® Protein 1 Calibrator

Protein 1 Calibrator is a multi-analyte, liquid human serum based product containing:

α ₁ -acid glycoprotein	immunoglobulin G
α ₁ -antitrypsin	immunoglobulin G Subclass 1
β ₂ -microglobulin	immunoglobulin G Subclass 2
C3 complement	immunoglobulin G Subclass 3
C4 complement	immunoglobulin G Subclass 4
ceruloplasmin	immunoglobulin M
haptoglobin	prealbumin
hemopexin	retinol binding protein
homocysteine	soluble transferrin receptor
immunoglobulin A	transferrin
immunoglobulin E	

Dimension Vista® Protein 1 Control L

Protein 1 Control L is a multi-analyte, low level liquid human serum based product containing :

α ₁ -acid glycoprotein	immunoglobulin G
α ₁ -antitrypsin	immunoglobulin G Subclass 1
C3 complement	immunoglobulin G Subclass 2
C4 complement	immunoglobulin G Subclass 3
ceruloplasmin	immunoglobulin G Subclass 4
haptoglobin	immunoglobulin M
hemopexin	prealbumin
homocysteine	retinol binding protein
immunoglobulin A	soluble transferrin receptor
immunoglobulin E	transferrin

Dimension Vista® Protein 1 Control M and H

Protein 1 Control M and H are multi-analyte, mid and high level respectively, liquid human serum based products containing :

α ₁ -acid glycoprotein	immunoglobulin G
α ₁ -antitrypsin	immunoglobulin G Subclass 1
β ₂ -microglobulin	immunoglobulin G Subclass 2
C3 complement	immunoglobulin G Subclass 3
C4 complement	immunoglobulin G Subclass 4
ceruloplasmin	immunoglobulin M
haptoglobin	prealbumin
hemopexin	retinol binding protein
homocysteine	soluble transferrin receptor
immunoglobulin A	transferrin
immunoglobulin E	

5. Device Intended Uses:

Dimension Vista® Protein 1 Calibrator

PROT1 CAL is an *in vitro* diagnostic product for the calibration of the Dimension Vista® System for:

α ₁ -Acid Glycoprotein (A1AG)	Immunoglobulin G (IGG)
α ₁ -Antitrypsin (A1AT)	Immunoglobulin G Subclass 1 (IGG1)
β ₂ -Microglobulin (B2MIC)	Immunoglobulin G Subclass 2 (IGG2)
C3 Complement (C3)	Immunoglobulin G Subclass 3 (IGG3)
C4 Complement (C4)	Immunoglobulin G Subclass 4 (IGG4)
Ceruloplasmin (CER)	Immunoglobulin M (IGM)
Haptoglobin (HAPT)	Prealbumin (PREALB)
Hemopexin (HPX)	Retinol Binding Protein (RBP)
Homocysteine (HCYS)	soluble Transferrin Receptor (STFR)
Immunoglobulin A (IGA)	Transferrin (TRF)
Immunoglobulin E (IGE)	

Dimension Vista® Protein 1 Control L

PROT1 CON L is an assayed, low level, intralaboratory quality control for assessment of precision and analytical bias on the Dimension Vista® System in the quantitative determination of:

α ₁ -Acid Glycoprotein (A1AG)	Immunoglobulin G (IGG)
α ₁ -Antitrypsin (A1AT)	Immunoglobulin G Subclass 1 (IGG1)
C3 Complement (C3)	Immunoglobulin G Subclass 2 (IGG2)
C4 Complement (C4)	Immunoglobulin G Subclass 3 (IGG3)
Ceruloplasmin (CER)	Immunoglobulin G Subclass 4 (IGG4)
Haptoglobin (HAPT)	Immunoglobulin M (IGM)
Hemopexin (HPX)	Prealbumin (PREALB)
Homocysteine (HCYS)	Retinol Binding Protein (RBP)
Immunoglobulin A (IGA)	soluble Transferrin Receptor (STFR)
Immunoglobulin E (IGE)	Transferrin (TRF)

Dimension Vista® Protein 1 Control M and H

PROT1 CON M and PROT1 CON H are assayed, mid-level and high level respectively, intralaboratory quality controls for assessment of precision and analytical bias on the Dimension Vista® System in the quantitative determination of:

α ₁ -Acid Glycoprotein (A1AG)	Immunoglobulin G (IGG)
α ₁ -Antitrypsin (A1AT)	Immunoglobulin G Subclass 1 (IGG1)
β ₂ -Microglobulin (B2MIC)	Immunoglobulin G Subclass 2 (IGG2)
C3 Complement (C3)	Immunoglobulin G Subclass 3 (IGG3)
C4 Complement (C4)	Immunoglobulin G Subclass 4 (IGG4)
Ceruloplasmin (CER)	Immunoglobulin M (IGM)
Haptoglobin (HAPT)	Prealbumin (PREALB)
Hemopexin (HPX)	Retinol Binding Protein (RBP)
Homocysteine (HCYS)	soluble Transferrin Receptor (STFR)
Immunoglobulin A (IGA)	Transferrin (TRF)
Immunoglobulin E (IGE)	

6. Medical device to which equivalence is claimed and comparison information:

The Dimension Vista® Protein 1 Calibrator and Dimension Vista® Protein 1 Control L, M and H modified to include hemopexin (HPX), immunoglobulin G subclass 1 (IGG1), immunoglobulin G subclass 2 (IGG2), immunoglobulin G subclass 3 (IGG3), immunoglobulin G subclass 4 (IGG4), and retinol binding protein (RBP) are substantially equivalent in intended use to the current Dimension Vista® Protein 1 Calibrator and Dimension Vista® Protein 1 Control L, M and H (K063663). The modified Dimension Vista® Protein 1 Calibrator and Protein 1 Control L, M, and H like the current products are intended to be used for the calibration of human protein assays and for use as assayed intralaboratory quality controls respectively, on the Dimension® Vista System.

7. Conclusion

The modified Dimension Vista® Protein 1 Calibrator and Dimension Vista® Protein 1 Control L, M and H, are substantially equivalent to the legally marketed devices based on the information described above.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Dade Behring, Inc.
c/o Ms AnnaMarie Kathleen Ennis
Regulatory Affairs and Compliance Specialist
Glasgow Site
P.O. Box 6101
Newark, DE 19714-6101

SEP 11 2007

Re: k071980

Trade/Device Name: Dimension Vista® Protein 1 Calibrator
Dimension Vista® Protein 1 Control H
Dimension Vista® Protein 1 Control M
Dimension Vista® Protein 1 Control L

Regulation Number: 21 CFR 862.1150

Regulation Name: Calibrator

Regulatory Class: Class II

Product Code: JIX, JJY

Dated: July 18, 2007

Received: July 18, 2007

Dear Ms. Ennis:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Robert L. Becker, Jr.", written in a cursive style.

Robert L. Becker, Jr., M.D., Ph.D.

Director

Division of Immunology and Hematology Devices

Office of In Vitro Diagnostic Device Evaluation and Safety

Center for Devices and Radiological Health

Enclosure

Indications Statement

Device Name: Dimension Vista® Protein 1 Calibrator

Indications for Use: K071980

Dimension Vista® Protein 1 Calibrator

PROT1 CAL is an *in vitro* diagnostic product for the calibration of the Dimension Vista® System for:

α_1 -Acid Glycoprotein (A1AG)	Immunoglobulin G (IGG)
α_1 -Antitrypsin (A1AT)	Immunoglobulin G Subclass 1 (IGG1)
β_2 -Microglobulin (B2MIC)	Immunoglobulin G Subclass 2 (IGG2)
C3 Complement (C3)	Immunoglobulin G Subclass 3 (IGG3)
C4 Complement (C4)	Immunoglobulin G Subclass 4 (IGG4)
Ceruloplasmin (CER)	Immunoglobulin M (IGM)
Haptoglobin (HAPT)	Prealbumin (PREALB)
Hemopexin (HPX)	Retinol Binding Protein (RBP)
Homocysteine (HCYS)	soluble Transferrin Receptor (STFR)
Immunoglobulin A (IGA)	Transferrin (TRF)
Immunoglobulin E (IGE)	

Prescription Use X
(Per 21 CFR 801 Subpart D)

Over-The-Counter-Use
(21 CFR 801)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Maria M Chan
Division Sign-Off

Office of In Vitro Diagnostic
Device Evaluation and Safety

510(k) K071980

CONFIDENTIAL

Indications Statement

Device Name: Dimension Vista® Protein 1 Control L

Indications for Use:

K071980

Dimension Vista® Protein 1 Control L

PROT1 CON L is an assayed, low level, intralaboratory quality control for assessment of precision and analytical bias on the Dimension Vista® System in the quantitative determination of:

α ₁ -Acid Glycoprotein (A1AG)	Immunoglobulin G (IGG)
α ₁ -Antitrypsin (A1AT)	Immunoglobulin G Subclass 1 (IGG1)
C3 Complement (C3)	Immunoglobulin G Subclass 2 (IGG2)
C4 Complement (C4)	Immunoglobulin G Subclass 3 (IGG3)
Ceruloplasmin (CER)	Immunoglobulin G Subclass 4 (IGG4)
Haptoglobin (HAPT)	Immunoglobulin M (IGM)
Hemopexin (HPX)	Prealbumin (PREALB)
Homocysteine (HCYS)	Retinol Binding Protein (RBP)
Immunoglobulin A (IGA)	soluble Transferrin Receptor (STFR)
Immunoglobulin E (IGE)	Transferrin (TRF)

Prescription Use X
(Per 21 CFR 801 Subpart D)

Over-The-Counter-Use
(21 CFR 801)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Maria M Chan
Division Sign-Off

Office of In Vitro Diagnostic
Device Evaluation and Safety

510(k) 071980

Indications Statement

Device Name: Dimension Vista® Protein 1 Control M

Indications for Use: K071980

Dimension Vista® Protein 1 Control M

PROT1 CON M is an assayed, mid-level, intralaboratory quality control for assessment of precision and analytical bias on the Dimension Vista® System in the quantitative determination of:

α_1 -Acid Glycoprotein (A1AG)	Immunoglobulin G (IGG)
α_1 -Antitrypsin (A1AT)	Immunoglobulin G Subclass 1 (IGG1)
β_2 -Microglobulin (B2MIC)	Immunoglobulin G Subclass 2 (IGG2)
C3 Complement (C3)	Immunoglobulin G Subclass 3 (IGG3)
C4 Complement (C4)	Immunoglobulin G Subclass 4 (IGG4)
Ceruloplasmin (CER)	Immunoglobulin M (IGM)
Haptoglobin (HAPT)	Prealbumin (PREALB)
Hemopexin (HPX)	Retinol Binding Protein (RBP)
Homocysteine (HCYS)	soluble Transferrin Receptor (STFR)
Immunoglobulin A (IGA)	Transferrin (TRF)
Immunoglobulin E (IGE)	

Prescription Use X
(Per 21 CFR 801 Subpart D)

Over-The-Counter-Use
(21 CFR 801)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Maria M. Chan
Division Sign-Off

Office of In Vitro Diagnostic
Device Evaluation and Safety

510(k) K071980

Indications Statement

Device Name: Dimension Vista® Protein 1 Control H

Indications for Use: K071980

Dimension Vista® Protein 1 Control H

PROT1 CON H is an assayed, high level, intralaboratory quality control for assessment of precision and analytical bias on the Dimension Vista® System in the quantitative determination of:

α_1 -Acid Glycoprotein (A1AG)	Immunoglobulin G (IGG)
α_1 -Antitrypsin (A1AT)	Immunoglobulin G Subclass 1 (IGG1)
β_2 -Microglobulin (B2MIC)	Immunoglobulin G Subclass 2 (IGG2)
C3 Complement (C3)	Immunoglobulin G Subclass 3 (IGG3)
C4 Complement (C4)	Immunoglobulin G Subclass 4 (IGG4)
Ceruloplasmin (CER)	Immunoglobulin M (IGM)
Haptoglobin (HAPT)	Prealbumin (PREALB)
Hemopexin (HPX)	Retinol Binding Protein (RBP)
Homocysteine (HCYS)	soluble Transferrin Receptor (STFR)
Immunoglobulin A (IGA)	Transferrin (TRF)
Immunoglobulin E (IGE)	

Prescription Use X
(Per 21 CFR 801 Subpart D)

Over-The-Counter-Use
(21 CFR 801)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Maria M Chan
Division Sign-Off

Office of In Vitro Diagnostic
Device Evaluation and Safety

510(k) K071980