

K072510

STERIS®



**510(k) Summary  
For  
Verify® V-PRO Chemical Indicator  
-Versions 1 and 2-**

DEC 11 2007

STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060  
Phone: (440) 354-2600  
Fax No: (440) 639-4459

Contact: Jack Scoville.  
Fellow  
Regulatory Affairs  
Telephone: (440) 392-7330  
Fax No: (440) 357-9198

Submission Date: September 05, 2007

**STERIS ABBREVIATED 510(k) PREMARKET NOTIFICATION  
VERIFY® V-PRO CHEMICAL INDICATOR – VERSIONS 1 AND 2-**

---

**1. Device Name**

Trade Name: Verify® V-PRO Chemical Indicator.

Models: Version 1: Verify® V-PRO Chemical Indicator.  
Version 2: Verify® V-PRO Chemical Indicator Adhesive Label.

Common Name: Chemical Indicator.

Classification Name: Physical/chemical sterilization process indicator (21 CFR 880.2800 (b), Product Code JOJ).

**2. Predicate Device**

Chemical Indicator Component of the STERRAD® CycleSure™ Biological Indicator (K994055)

**3. Device Description**

The Verify® V-PRO Chemical Indicator is used in each processing cycle to indicate exposure to an Amsco V-PRO 1 Low Temperature Sterilizer cycle. When exposed, the indicator exhibits a visible color change from magenta to yellow.

The Verify® V-PRO Chemical Indicator is provided as two formats:

- Version 1: Verify® V-PRO Chemical Indicator
- Version 2: Verify® V-PRO Chemical Indicator Adhesive Label

The Version 1: Verify® V-PRO Chemical Indicator is a Class 1 process indicator in accordance with ISO 11140-1:2005 which consists of the chemical indicator applied to an inert polymeric substrate; the indicator spot is laminated with a transparent laminate.

The Version 2: Verify® V-PRO Chemical Indicator Adhesive Label is a Class 1 process indicator in accordance with ISO 11140-1:2005 which consists of the chemical indicator applied to a spun bonded polyolefin substrate with an adhesive supplied on a backing paper.

**4. Intended Use:**

The Verify® V-PRO Chemical Indicator (Version 1) and the Verify® V-PRO Chemical Indicator Adhesive Label (Version 2) are Class 1 vaporized hydrogen peroxide sterilization process indicators. They are designed to distinguish between

processed and unprocessed units when placed within (Version 1) or affixed to (Version 2) sterilization wraps, trays or pouches to indicate, through a visible change from magenta to yellow, when the device (Version 1) or pack (Version 2) has been exposed to a V-PRO 1 sterilization process. This product is designed for use exclusively in the Amsco V-PRO 1 Low Temperature Sterilization system at 50°C using Vaprox™ HC Sterilant.

**5. Description of Safety and Substantial Equivalence**

The proposed and predicate devices are single use process indicators for use in monitoring Vaporized Hydrogen Peroxide sterilization cycles. The differences between the proposed Verify® V-PRO Chemical Indicator and the predicate device are limited to differences in design, material, and parameters of the sterilization cycles these indicators are designed to monitor. These differences do not raise any new issues of safety and efficacy.

**6. Performance Testing**

Performance testing was conducted to verify that the proposed Verify® V-PRO Chemical Indicator meets the requirements for Class 1 vaporized hydrogen peroxide sterilization indicators as defined in ANSI/AAMI ISO 11140-1:2005. Additional testing was completed to simulate typical in-use applications and testing was also performed to investigate the effects of exposure to UV, visible light and aggressive chemicals to the performance of the Verify® V-PRO Chemical Indicator.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC 11 2007

Mr. Jack Scoville  
Fellow, Regulatory Affairs  
STERIS, Corporation  
5960 Heisley Road  
Mentor, Ohio 44060-1834

Re: K072510  
Trade/Device Name: Verify<sup>®</sup> V-PRO Chemical Indicator (Versions 1 and 2)  
Regulation Number: 880.2800  
Regulation Name: Sterilization Process Indicator  
Regulatory Class: II  
Product Code: JOJ  
Dated: November 20, 2007  
Received: November 21, 2007

Dear Mr. Scoville:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin, Ph.D.", written in a cursive style.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

**STERIS REQUEST FOR ADDITIONAL INFORMATION**  
**K072510/S001 Verify® V-PRO Chemical Indicator (Versions 1 and 2)**

---

**Indications for Use**

510(k) Number (if known): K072510

Device Name: Verify® V-PRO Chemical Indicator (Versions 1 and 2).

Indications For Use:

The Verify® V-PRO Chemical Indicator (Version 1) and the Verify® V-PRO Chemical Indicator Adhesive Label (Version 2) are Class 1 vaporized hydrogen peroxide sterilization process indicators. They are designed to distinguish between processed and unprocessed units when placed within (Version 1) or affixed to (Version 2) sterilization wraps, trays or pouches to indicate, through a visible change from magenta to yellow, when the device (Version 1) or pack (Version 2) has been exposed to a V-PRO 1 sterilization process. This product is designed for use exclusively in the Amsco V-PRO 1 Low Temperature Sterilization system at 50 °C using Vaprox™ HC Sterilant.

Prescription Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use  X   
(21 CFR 801 Subpart C)

**(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER  
PAGE IF NEEDED)**

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Shirley R. Murphy MD*  
*Medical Director*  
*K 072510*