

JAN 10 2008

## 510(k) Summary

## 510(k) Summary

**Trade Name:** SPY<sup>®</sup> Fluorescent Imaging System

**Device Model Number:** SP2001

**Common Name:** Fluorescent Angiographic System

**Classification:** 21 CFR 892.1600

**Product Code:** 90 IZI

**Classification:** Class II

**Manufacturer:** Novadaq Technologies Inc.  
2585 Skymark Avenue  
Suite 306  
Mississauga, Ontario  
Canada  
L4W 4L5  
905.629.3822 ext. 240

**Contact Name:** Allison Manners  
Vice President, Regulatory and Clinical Affairs

**Date 510(k) Summary Prepared:** December 15, 2007

### Legally Marketed Predicate Devices:

The Novadaq<sup>®</sup> SPY Fluorescent Imaging System (SPY System with the SP2000 Imaging Device) had received initial FDA 510(k) clearance for market for use during coronary artery bypass graft (CABG) surgery in January 2005 (K042961). Subsequent 510(k) clearance was obtained in May 2006 (K060867) for a labeling change, followed by a clearance for use in plastic, micro- and reconstructive surgery in January 2007 (K063345), and a clearance for use of an alternative brand of fluorescent ICG agent in coronary artery bypass surgery in May 2007 (K071037), as well as in the plastic, micro-, and reconstructive surgery in September 2007 (K072222). In addition, the SPY Fluorescent Imaging System (SPY System with the SP2000 Imaging Device) recently received FDA 510(k) clearance for market for its use in cardiovascular surgery (K071619).

The Leica FL800 had received FDA 510(k) clearance for market in September 2006 (K061871). The Leica FL800 is intended for use to allow neurosurgeons to view blood flow.

## Device Description:

The SPY Fluorescent Imaging System is currently cleared for use:

- For intra-operative visual assessment of blood flow in vessels and related tissue perfusion during cardiovascular surgical procedures.
- As an imaging system used in capturing and viewing fluorescent images for the visual assessment of blood flow as an adjunctive method for the evaluation of tissue perfusion, and related tissue-transfer circulation in tissue and free flaps used in plastic, micro- and reconstructive surgical procedures.

The Novadaq® Technologies SPY® Fluorescent Imaging System used in plastic, micro, and reconstructive surgery indication consists of 2 components:

- the SP2001 Imaging Device
- the SPY Paq®

The SPY Paqs are available in 2 configurations dependent on the intended indication for use<sup>1</sup>. Each SPY Paq contains sufficient number of ICG, custom sterile drapes called Novadrape®, and diluent for 6 imaging procedures. Each configuration of SPY Paq has a unique part number assigned to it, and different Instructions for Use exist for the two types of SPY Paqs. The different Instructions for Use also have unique part numbers for ease of assembly of the Paqs and guidance for the end user.

### The SP2000 Imaging Device

The SP2000 Imaging Device consists of an imaging head containing a charge coupled device (CCD) camera, a laser light source, motion sensor, and distance sensor attached via an articulating arm to a mobile cart. The mobile cart contains a flat panel display, a computer with keyboard and optical mouse, an electronics enclosure, and a printer.

The SPY® System provides the surgeon with the capability to view, record, and replay fluorescent images of blood flow in vessels and organs. A laser light source is used to illuminate the area of interest. In order to obtain the images, ICG is injected intravenously through the central or peripheral venous line, bypass pump, cardioplegia line, or down a coronary graft. While the ICG is

---

<sup>1</sup> Novadaq provides ICG as it is sold by the manufacturer and does not adulterate the integrity of the original packaging or labeling. IC-Green™ (Akorn, Inc.) is packaged in a kit that contains 6 x 25 mg vials of ICG and 6 diluents used to dissolve the ICG. A single IC-Green kit constitutes a 6 procedure SPY Paq for the organ transplant indication.

passing through the vessels, the absorption of laser light causes excitation of the ICG dye, followed by the emission of infrared energy. A CCD camera of the SP2000 Imaging Device captures the infrared emission, resulting in a fluorescent image of blood flow and related tissue perfusion. These images are used to evaluate the integrity of native and grafted vasculature and blood flow in the organs.

There have been no significant changes or modifications made to the SP2000 Imaging Device since the original 510(k) clearance in January 2005 premarket notification 510(k) K042961, the 510(k) premarket notification K060867 submitted for a label change for this device, the 510(k) clearance K063345 for use in the plastic, micro- and reconstructive surgery, and the clearances K071037 and K072222 for use of an alternative brand of fluorescent dye – ICG PULSION® in May 2007 and September 2007, respectively.

### The SP2001 Imaging Device

The SP2001 Imaging Device represents a modification of the SP2000 Imaging Device in the following ways:

- The maximum recording time for image sequences captured has been extended to 60 seconds from the 34 seconds of the SP2000 Imaging Device. This modification has been carried out to permit the users to observe the entire cycle of blood flow in the area of interest under imaging; namely the arterial inflow, the perfusion, and the venous outflow. Accomplishing this change required a replacement of watch-dog component on the controller board inside the electronics enclosure with a different one, as well as a software change. The watch-dog hardware component controls the duration that the laser is turned on
- An ability to move the SP2001 camera head in horizontal plane during image sequence acquisition has been implemented. This modification has been carried out to permit the users to image the area larger than the 7.6 cm x 7.6 cm which is imaged by a stationary camera of the SP2000 Imaging Device. This feature is very desirable, given that an organ can be large or extended. This modification required disabling of the accelerometer on the distribution board inside the camera head.
- A variable Laser Power Attenuator has been added to the device to permit laser illumination at power levels lower or equal to the laser power level of the SP2000 Imaging Device. This modification has been carried out to permit the SP2001 operators to lower laser intensity, in case the signal from the area of interest would be saturating the sensitivity of the camera. This new feature was accomplished by connecting a control box to the controller board inside the electronics enclosure.
- New HELIOS™ software has been developed to support the SPY System in its organ transplant surgery indication. The structure and workflow of the HELIOS software is based on the cleared DaqPac software. The

HELIOS software concentrated on providing the user with the most suitable interface and procedure reporting and it did not alter the way that the device is used.

The modifications of the SP2000 Imaging Device were introduced to provide the end users the increased functionality for using the SPY System in the organ transplant surgery and plastic, micro, and reconstructive surgery (traditional 510(k) currently under review, no 510(k) number yet assigned).

It is demonstrated in Section 10 – Device Description that the modifications do not introduce any additional or new concerns regarding the safety and use of the device.

SPY® Systems intended for use solely in CABG and cardiovascular surgery utilize the SP2000 Imaging Device. SPY Systems intended for use in the organ transplant surgery utilize the SP2001 Imaging Device.

This traditional premarket 510(k) notification is being made to expand the Indications for Use for the SPY® Fluorescent Imaging System in organ transplant procedures; in addition to the already cleared Indications for Use in CABG and cardiovascular surgery, as well as in plastic, micro-, and reconstructive surgery. These cleared Indications for Use are not being amended in any way with this submission.

#### **Proposed Intended Use of the SPY System:**

The SPY® Fluorescent Imaging System is intended to intra-operatively enable surgeons to visually assess blood flow and related tissue perfusion during organ transplant procedures.

#### **Testing:**

Animal studies, human experience and *in vitro* testing were conducted to support the safe and effective use of the SPY System in its original premarket notification 510(k) application (K042961).

The information contained within this Traditional premarket notification 510(k) demonstrates the utility of the SPY System in organ transplant surgery in addition to the previously cleared indications.

#### ***In Vitro* Testing:**

Testing of the SPY System was completed in conformance with the following standards. The SPY System successfully met all of the requirements for these standards.

1. Electrical per IEC 60601-1 and UL2601-1

2. Electromagnetic Compatibility per IEC 60601-1-2
3. Light Emitting Laser Products per 21 CFR 1040
4. Safe Use of Lasers in Health Care Facilities per ANSI Z136.3
5. American National Standard for Safe Use of Lasers per ANSI Z136.1

### ***In Vivo Testing:***

The SPY® System is commercially available in the United States of America, Japan, Europe, and Canada. To date, the SPY System has been used in over 7000 vascular procedures in humans and there have been no reports of adverse acute or long-term cellular, renal or hepatic effects. The data from intra-operative imaging in CABG, cardiovascular, as well as in plastic, micro-, reconstructive surgery demonstrated the clinical utility of the device in producing high quality and resolution images of the entire vascular bed of the area of interest.

Results from the use of the SPY System has been the subject of 15 peer reviewed journal articles, 13 related to its use in cardiac surgery and 2 related to its use in transplantation kidney and liver surgeries. Please refer to the bibliography in Section 19 - Clinical for a listing of all relevant journal articles.

The literature reports that when used in CABG surgery, the SPY System was able to non-invasively, quickly and safely identify 17 conduits in 311 patients that required revision during the surgical procedures. In all cases the lack of patency was visualized clearly by the SPY System using doses of ICG well below the maximum dose approved for human use. It allowed the surgeon to revise the graft, decreasing subsequent myocardial infarctions and the morbidity and mortality associated with poor graft patency. Cardiac, renal and hepatic functions were monitored during use of the SPY System and there were no reported adverse effects.

To support the original traditional 510(k) premarket notification for the indication of intra-operative visual assessment of blood flow in vessels and related tissue perfusion during cardiovascular surgical procedures, the system was used in six pig studies, please refer to Section 18 – Performance Testing: Animals of this submission for full details. These studies were primarily focused on the effects of the SPY System on the heart. They demonstrated that:

- 1) It was possible to acquire high quality images in a simple and reproducible manner using small doses of ICG well below the concentrations approved for human use.
- 2) It was possible to perform multiple imaging sequences with no detrimental effects on heart function, coronary flow or peripheral pressure.
- 3) It was possible to acquire images with no increase in myocardial tissue temperature.

- 4) It was possible to visualize all of the coronary beds with high quality images even when the heart was in a vertical position for visualizing posterior arteries.

Therefore, the *in vivo* human clinical experience with the SPY<sup>®</sup> System used for fluorescent imaging during organ transplant procedures (31 cases in the United States and 18 reported in the literature from Japan<sup>2,3</sup>) shows that:

1. The exposure for the SP2001 Imaging Device at the imaging distance of 30 cm is 31.2 mW/cm<sup>2</sup> which is far below the maximum permissible exposure (MPE) of 326 mW/cm<sup>2</sup> established by ANSI for exposure to the skin. The exposure remains below the MPE value as long as the camera head of SP2001 is not advanced closer than 5 cm from the region of interest imaged. The extension of the recording time to 60 seconds does not affect the MPE value.
2. There were no intra-operatively observed renal or hepatic effects of using the SPY System.
3. The SPY System enabled surgeons to assess blood flow.
4. The device can only be used during open operations with direct visualization of the organ and vessels of interest. When using the device, there may be a potential inability to visualize the entire vascular supply.

### Conclusions:

The testing demonstrates that the SPY Fluorescent Imaging System is equivalent to predicate devices and intra-operatively enables surgeons to visually assess blood flow and related tissue perfusion during organ transplant procedures.

---

2 Kubota, K et al. *Intra-operative assessment of reconstructed vessels in living-donor liver transplantation, using a novel fluorescence imaging technique.* Journal of Hepatobiliary Pancreatic Surgery. 2006; 13: 100-104.

3 Sekijima, M et al. *An intra-operative fluorescent imaging system in organ transplantation.* Transplantation Proceedings. 2006; 36: 2188-2900.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JAN 10 2008

Novadaq Technologies, Inc.  
% Ms. Allison Manners  
VP, Regulatory and Clinical  
Affairs  
2585 Skymark Avenue, Suite 306  
Mississauga, Ontario, Canada  
L4W 4L5

Re: K073130

Trade/Device Name: SPY<sup>®</sup> Fluorescent Imaging System  
Regulation Number: 21 CFR 892.1600  
Regulation Name: Angiographic x-ray system  
Regulatory Class: II  
Product Code: IZI  
Dated: November 5, 2007  
Received: November 13, 2007

Dear Ms. Manners:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set



forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K 073130

**Device Name:** SPY® Fluorescent Imaging System

**Indications for Use:**

The SPY® Fluorescent Imaging System is intended to intra-operatively enable surgeons to visually assess blood flow and related tissue perfusion during organ transplant procedures.

Prescription Use     X      
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

(Division Sign-Off)

## Division of General, Restorative, and Neurological Devices

Page 1 of 1

Slack Number 14073130