

K080950

Otto Bock®

QUALITY FOR LIFE

APR - 2 2009

Summary of Safety and Effectiveness STIWELL med4**A Manufacturer:**

Otto Bock Healthcare Product GmbH
Kaiserstrasse 39
1070 Vienna
Austria

Telephone: +43 1 523 37 86 692
Fax Number: +43 1 523 22 64
Contact: Juergen Weiss
Regulatory Affairs Manager

B Contact:

Otto Bock Healthcare Product GmbH
Kaiserstrasse 39
1070 Vienna
Austria

Telephone: +43 1 523 37 86 692
Fax Number: +43 1 523 22 64
Contact: Juergen Weiss
Regulatory Affairs Manager

C Product Name/ Classification Name:

Common or Usual Name: Powered Muscle Stimulator
Proprietary or Trade Name: STIWELL med4
Model Number: 900101S

Class	Regulation No.	Device Classification Name	Product Code	Speciality
II	890.5850	Stimulator, Muscle, Powered	IPF	Physical Medicine
II	882.5890	Stimulator, Nerve, Transcutaneous, For Pain Relief	GZI	Neurology
II	882.5050	Device, Biofeedback	HCC	Neurology
II	882.5810	Stimulator, Neuromuscular, External Functional	GZI	Neurology
II	876.5320	Stimulator, Electrical, Non- Implantable, For Incontinence	KPI	Gastroenterology/ Urology


Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 64 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400
Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

D Predicate Devices:

510(k) No.	Predicate Device	Manufacturer	Regulation Number	Product Code(s)	Clearance Date
K032954	Elpha Models II 500, II 1000, II 2000, II 3000 A	Danmeter A/S	890.5850	IPF, GZJ	08/09/2004
K053434	Myotrac Infinity Electrical Stimulator	Thought Technology Ltd.	876.5320	KPI, HCC, IPF	03/15/2006
K040849	Mentamove	Cole & Associates	890.5850	IPF, HCC	10/01/2004
K940301	Compex 2	Biodex Medical Systems, Inc.	890.5850	IPF, GZJ	12/28/1994
K032988	Elpha 4 Conti	Danmeter A/S	876.5320	KPI	12/03/2003
K031900	Handmaster	Ness- Neuromuscular Electrical Stimulation Systems	882.5810	GZI	08/08/2003

E Description:

The STIWELL med4 is powered by rechargeable batteries. It has four stimulation channels and two EMG measurement channels.

The STIWELL med4 is intended for stationary use in a hospital as well as home use by the patient. The physician/therapist has the flexibility to adjust the programs and monitors the progress of the therapy. Statistics regarding the completed treatments can be retrieved from the STIWELL med4 and from the PC.

In order to gain a proper understanding of STIWELL med4, it is important to read the manual before beginning to use the STIWELL med4.

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 64 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400
Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

F Substantial Equivalence:**Powered Muscle Stimulator:**

Characteristics/ Specific.	New Device	Predicate Device	Predicate Device	Predicate Device	Predicate Device
Basic Unit Characteristics					
510(k) Number	K080950	K032954	K053434	K040849	K940301
Manufacturer	Otto Bock	Danmeter A/S	Thought Technology Ltd.	Mentamove North America, LLC	Comper SA
Device Name, Model	STWELL med4	Elpha II3000	Myotrac Infinity	Mentamove	Comper 2
Power Source(s)	Battery Pack Li-Ion 11.1V	Battery NMH or Alkaline 9V	Battery Pack NiMH rechargeable	Battery Pack NiMH rechargeable	Battery Pack NiMH 7.2 V
Method of Line Current Isolation	Medical Class II Power Adapter – Moscott (12.6VDC-15.1V)	N/A	Medical Class II Power Adapter (6VDC-15W)	Battery Charger (6VDC-2.1W)	Battery Charger (11VDC-7.2W)
Patient Leakage Current (normal condition)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)
Patient Leakage Current (single fault condition)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)
Number of Output Modes	1	1	1	1	1
Number of Output Channels	4	2	2	1	4
Number of EMG (input) Channels	2	N/A	2	1	N/A
EMG sensitivity	1µV	N/A	<0.1µV	1µV	N/A
EMG Sampling Rate	3kHz	N/A	2,048kHz	Unknown	N/A
EMG detection (bipolar/ monopolar)	Bipolar	N/A	Bipolar	Bipolar	N/A
EMG range (µV)	1-2000µV	N/A	0-2000µV	2-2000µV	N/A
EMG bandwidth	70-480Hz	N/A	10Hz-1kHz	Unknown	N/A
EMG signal processing (eg. RMS)	AVR (Average Rectified Value)	N/A	RMS (Root Mean Square)	Unknown	N/A
Synchronous or Alternating?	Alternating	Unknown	Alternating	N/A	Synchronous
Method of Channel Isolation	Transformer, Inductive couplers	Unknown	Unknown	Unknown	Transformer
Regulated Current or Regulated Voltage?	Regulated Current	Regulated Current	Regulated Current	Regulated Current	Regulated Current
Software/Firmware/ Microprocessor Control?	Yes	Yes	Yes	Yes	Yes
PC Software?	Yes	No	Yes	No	No
Automatic Overload Trip?	Yes	Yes	Yes	Yes	Yes
Automatic No-Load Trip?	Yes	Yes	Yes	Yes	Yes
Automatic Shut Off?	Yes (10min)	Yes	Unknown	Unknown	Unknown
Patient Override Control?	Yes (Stop Button)	Yes	Yes	Yes	Yes
Indicator Display: On/Off Status?	Yes	Yes	Yes	Yes	Yes
Indicator Display: Low Batt.?	Yes	Yes	Yes	Yes	Yes
Indicator Display: Voltage/Current Level?	Yes	Yes	Yes	Yes	Yes
Timer Range (minutes)	2 – 120min	5 – 95min and cont.	1 – 120min	Unknown	3 – 64min
Compliance with Voluntary Standards? (if yes, specify)	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown	Unknown	Unknown
Compliance with 21 CFR 888?	Yes	Unknown	Unknown	Unknown	Unknown
Weight	440g	158g	330g	Unknown	420g
Dimensions (WxHxD) in (mm)	175x95x30	114x80x31	102x152x51	Unknown	148x80x30
Housing Material and Constr.	Plastics	Plastics	Plastics	Plastics	Aluminum
Output Specifications					
Waveform	Biphasic symmetrical	Biphasic asymmetrical	Biphasic asymmetrical	Biphasic asymmetrical	Biphasic asymmetrical
Shape	Rectangular	Rectangular with discharge	Rectangular with discharge	Rectangular	Rectangular, Triangular
Maximum Output Voltage (500Ω)	50V	50V	50V	Unknown	50V
Maximum Output Voltage (2kΩ)	115V	164V	Unknown	Unknown	200V
Maximum Output Voltage (10kΩ)	N/A	Unknown	Unknown	Unknown	N/A ^{a1}
Maximum Output Current (500Ω)	100mA	100mA	100mA	Unknown	100mA
Maximum Output Current (2kΩ)	58mA	82mA	Unknown	Unknown	100mA
Maximum Output Current (10kΩ)	N/A	Unknown	Unknown	Unknown	N/A
Pulse Width (specify units)	50µs – 400µs	100µs – 400µs	50µs – 400µs	320µs	15µs – 999µs
Frequency (Hz)	1 – 140Hz	2 – 100Hz	2 – 100Hz	3, 1kHz GF, 33Hz AMF	1 – 200Hz
For interlaminar modes only: Beat Frequency (Hz)	N/A	N/A	N/A	N/A	N/A
For multiphasic waveforms only: Symmetrical phases?	N/A	N/A	N/A	N/A	N/A
For multiphasic waveforms only: Phase Duration (including units)	N/A	N/A	N/A	N/A	N/A
Net Charge (µC per pulse): (500Ω)	0µC Some positive and negative impulse	0µC Output capacitor	0µC Output capacitor	0µC Some positive and negative impulse	0µC Some positive and negative impulse
Maximum Phase Charge (µC) (500Ω)	40µC	40µC	60µC	Unknown	80µC
Maximum Current Density (mA/cm ²) (500Ω)	12.5mA/cm ²	5.9mA/cm ²	20.4mA/cm ²	Unknown	5.0mA/cm ²
Maximum Power Density (W/cm ²) (500Ω)	7.9mW/cm ²	1.7mW/cm ²	6.5mW/cm ²	Unknown	39.9mW/cm ²
Burst Mode (i.e., pulse trans): Pulses per burst	N/A	N/A	N/A	N/A	N/A
Burst Mode (i.e., pulse trans): Bursts per second	N/A	N/A	N/A	N/A	N/A
Burst Mode (i.e., pulse trans): Burst duration (seconds)	N/A	N/A	N/A	N/A	N/A
Burst Mode (i.e., pulse trans): Duty Cycle (Line (b) • Line (c))	N/A	N/A	N/A	N/A	N/A
ON Time (seconds)	1 – 20s	2 – 20s	2 – 20s	Default 8s	Unknown
OFF Time (seconds)	1 – 50s	2 – 20s	2 – 50s	Default 12s	Unknown
Additional features	N/A	N/A	N/A	N/A	N/A


Otto Bock Healthcare Products GmbHKaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 64 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400

Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

Functional Electrical Stimulation Programs:

Characteristics/ Specifications	New Device	Predicate Device	New Device	Predicate Device	New Device	Predicate Device	New Device	Predicate Device
Basic Unit Characteristics								
Program	FES 1 Grasp/Release	Exercise Open	FES 2 Grasp/Release (EMG)	Exercise Open	FES 3 Open/Close	Exercise	FES 4 Open/Close (EMG)	Palmar Grasp
510(a) Number	K080950	K031900	K080950	K031900	K080950	K031900	K080950	K031900
Manufacturer	Otto Bock	Ness	Otto Bock	Ness	Otto Bock	Ness	Otto Bock	Ness
Device Name	STWELL med4	Handmaster	STWELL med4	Handmaster	STWELL med4	Handmaster	STWELL med4	Handmaster
Power Source(s)	Battery Pack Li-Ion 11,1V	Battery Pack Ni-Cd 9,6V	Battery Pack Li-Ion 11,1V	Battery Pack Ni-Cd 9,6V	Battery Pack Li-Ion 11,1V	Battery Pack Ni-Cd 9,6V	Battery Pack Li-Ion 11,1V	Battery Pack Ni-Cd 9,6V
Method of Line Current Isolation	Medical Class II Power Adapter	N/A	Medical Class II Power Adapter	N/A	Medical Class II Power Adapter	N/A	Medical Class II Power Adapter	N/A
Leakage Current (normal condition)	N/A (Battery)	N/A	N/A (Battery)	N/A	N/A (Battery)	N/A	N/A (Battery)	N/A
Leakage Current (single fault cond.)	N/A (Battery)	N/A	N/A (Battery)	N/A	N/A (Battery)	N/A	N/A (Battery)	N/A
No. Output Mod.	1	1	1	1	1	1	1	1
No. Output Chan.	3	2	3	2	3	3	3	5
Simulated Muscl.	Wrist extensors Finger flexors Thumb flexor	Finger extensors Thumb extensor	Wrist extensors Finger flexors Thumb flexor	Finger extensors Thumb extensor	Finger/thumb extensors Finger flexors Thumb flexor	Finger extensors Thumb extensor Finger flexors Thumb flexor Thenar muscle	Finger/thumb extensors Finger flexors Thumb flexor	Finger extensors Thumb extensor Finger flexors Thumb flexor Thenar muscle
No. of EMG Chan.	0	0	1	0	0	0	1	0
EMG sensitivity	N/A	N/A	1µV	N/A	N/A	N/A	1µV	N/A
EMG Sampl. Rate	N/A	N/A	3kHz	N/A	N/A	N/A	3kHz	N/A
EMG detection	N/A	N/A	Bipolar	N/A	N/A	N/A	Bipolar	N/A
EMG range (µV)	N/A	N/A	1-2000µV	N/A	N/A	N/A	1-2000µV	N/A
EMG bandwidth	N/A	N/A	70-480Hz	N/A	N/A	N/A	70-480Hz	N/A
EMG signal proc.	N/A	N/A	AVR	N/A	N/A	N/A	AVR	N/A
Synch. or Altern.?	Alternating	Alternating	Alternating	Alternating	Alternating	Alternating	Alternating	Alternating
Meth. Chn. Isol.	Transformer, Inductive couplers	Transformer, Opt. isolator	Transformer, Inductive couplers	Transformer, Opt. isolator	Transformer, Inductive couplers	Transformer, Opt. isolator	Transformer, Inductive couplers	Transformer, Opt. isolator
RC or RV?	Regulated Current	Regulated Voltage	Regulated Current	Regulated Voltage	Regulated Current	Regulated Voltage	Regulated Current	Regulated Voltage
Firmware?	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
PC Software?	No	No	No	No	No	No	No	No
Aut. Overload Trip	Yes	Unknown	Yes	Unknown	Yes	Unknown	Yes	Unknown
Aut. No-Load Trip	Yes	Unknown	Yes	Unknown	Yes	Unknown	Yes	Unknown
Aut. Shut Off	Yes (10min)	No	Yes (10min)	No	Yes (10min)	No	Yes (10min)	No
Override Control?	Yes (Stop Button)	Yes	Yes (Stop Button)	Yes	Yes (Stop Button)	Yes	Yes (Stop Button)	Yes
Display: On/Off?	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Display: Low Bat?	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Display: V/C Level	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Timer Range	15-60min	Max. 90min	15-60min	Max. 90min	15-60min	Max. 90min	15-60min	Max. 45min
Compliance with Voluntary Standards?	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown
Compliance with 21 CFR 898?	Yes	Unknown	Yes	Unknown	Yes	Unknown	Yes	Unknown
Weight	440g	685g	440g	685g	440g	685g	440g	685g
Dimensions (WxHxD) in [mm]	175x95x30	172x91x40 (Control unit)	175x95x30	172x91x40 (Control unit)	175x95x30	172x91x40 (Control unit)	175x95x30	172x91x40 (Control unit)
Housing Material and Construction	Plastics	Plastics	Plastics	Plastics	Plastics	Plastics	Plastics	Plastics
Output Specifications								
Waveform	Biphasic symmetrical	Biphasic symmetrical	Biphasic symmetrical	Biphasic symmetrical	Biphasic symmetrical	Biphasic symmetrical	Biphasic symmetrical	Biphasic symmetrical
Shape	Rectangular	Sinusoidal	Rectangular	Sinusoidal	Rectangular	Sinusoidal	Rectangular	Sinusoidal
Maximum Output Voltage (500Ω)	50V	106V	50V	106V	50V	106V	50V	106V
Maximum Output Voltage (2kΩ)	116V	163V	115V	163V	115V	163V	115V	163V
Maximum Output Voltage (10kΩ)	N/A	Unknown	N/A	Unknown	N/A	Unknown	N/A	Unknown
Maximum Output Current (500Ω)	100mA	212mA	100mA	212mA	100mA	212mA	100mA	212mA
Maximum Output Current (2kΩ)	58mA	92mA	58mA	92mA	58mA	92mA	58mA	92mA
Maximum Output Current (10kΩ)	N/A	Unknown	N/A	Unknown	N/A	Unknown	N/A	Unknown
Pulse Width	50-400µs	100-350µs	50-400µs	100-350µs	50-400µs	100-350µs	50-400µs	100-500µs
Frequency [Hz]	1-140Hz Default: 35Hz	38Hz	1-140Hz Default: 35Hz	38Hz	1-140Hz Default: 35Hz	38Hz	1-140Hz Default: 35Hz	18Hz
Net Charge (µC per pulse) (500Ω)	0µC	0µC	0µC	0µC	0µC	0µC	0µC	0µC
Maximum Phase Charge (µC) (500Ω)	40µC	75µC	40µC	75µC	40µC	75µC	40µC	106µC
Maximum Current Density (mA/cm²) (500Ω)	12,5mA/cm²	15,9mA/cm²	12,5mA/cm²	15,9mA/cm²	12,5mA/cm²	15,9mA/cm²	12,5mA/cm²	15,9mA/cm²
Maximum Power Density (W/cm²) (500Ω)	7,9mW/cm²	0,39 mW/cm²	7,9mW/cm²	0,39 mW/cm²	7,9mW/cm²	0,39 mW/cm²	7,9mW/cm²	0,39 mW/cm²
Burst Mode: Pulses per burst	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Burst Mode: Bursts per second	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Burst Mode: Burst duration (seconds)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Burst Mode: Duty Cycle	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
ON Time [sec.]	1 - 20s	Unknown	1 - 20s	Unknown	1 - 20s	Unknown	1 - 20s	Trigger controlled
OFF Time [sec.]	1 - 30s	Unknown	Trigger controlled (min. 1 - 30s)	Unknown	1 - 30s	Unknown	Trigger controlled (min. 1 - 30s)	Trigger controlled
Additional features	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 64 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400

Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAUTWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

Incontinence Programs:

Characteristics/ Specifications	New Device	Predicate Device	New Device	Predicate Device	New Device	Predicate Device	New Device	Predicate Device
Basic Unit Characteristics								
Program	U1 Urge	P2, P3 Urge	U2 Stress	P1 Stress	U3 Mixed	P4 Mixed	U4 Feedback	EMG Script
510(k) Number	K080950	K032988	K080950	K032988	K080950	K032988	K080950	K053434
Manufacturer	Otto Bock	Donmeter A/S	Otto Bock	Donmeter A/S	Otto Bock	Donmeter A/S	Otto Bock	Thought Technology Ltd.
Device Name	STWELL med4	Elpha 4 Cont.	STWELL med4	Elpha 4 Cont.	STWELL med4	Elpha 4 Cont.	STWELL med4	Myotrac Infinity
Power Source(s)	Battery Pack Li-Ion 11.1V	Battery NMH or Alkaline 9V	Battery Pack Li-Ion 11.1V	Battery NMH or Alkaline 9V	Battery Pack Li-Ion 11.1V	Battery NMH or Alkaline 9V	Battery Pack Li-Ion 11.1V	Battery Pack NMH
Method of Line Current Isolation	Medical Class II Power Adapter	N/A	Medical Class II Power Adapter	N/A	Medical Class II Power Adapter	N/A	Medical Class II Power Adapter	Medical Class II Power Adapter
Leakage Current (normal condition)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)
Leakage Current (single fault cond.)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)
No. Output Mod.	1	1	1	1	1	1	0	0
No. Output Chom.	1	1	1	1	1	1	0	0
No. of EMG Chom.	0	0	0	0	0	0	1	1
EMG sensitivity	N/A	N/A	N/A	N/A	N/A	N/A	1µV	<0.1µV
EMG Samp. Rate	N/A	N/A	N/A	N/A	N/A	N/A	5KHz	2,048Hz
EMG detection	N/A	N/A	N/A	N/A	N/A	N/A	Bipolar	Bipolar
EMG range (µV)	N/A	N/A	N/A	N/A	N/A	N/A	1-2000µV	0-2000µV
EMG bandwidth	N/A	N/A	N/A	N/A	N/A	N/A	70-480Hz	10Hz-1kHz
EMG signal proc.	N/A	N/A	N/A	N/A	N/A	N/A	AVR	RMS
Synch. or Assem.?	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Math. Chom. Isol.	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
RC or RV?	Regulated Current	Regulated Current	Regulated Current	Regulated Current	Regulated Current	Regulated Current	N/A	N/A
Firmware?	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Aut. Overload Trip	Yes	Yes	Yes	Yes	Yes	Yes	N/A	N/A
Aut. No-Load Trip	Yes	Yes	Yes	Yes	Yes	Yes	N/A	N/A
Aut. Shut Off	Yes (10min)	Yes	Yes (10min)	Yes	Yes (10min)	Yes	Yes (10min)	Unknown
Override Control?	Yes (Stop Button)	Yes	Yes (Stop Button)	Yes	Yes (Stop Button)	Yes	Yes (Stop Button)	Yes
Display: On/Off?	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Display: Low Batt?	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Timer Range	5 - 25min	5 - 95min	5 - 25min	5 - 95min	5 - 25min	5 - 95min	2-4min	1-120min
Compliance with Voluntary Standards?	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown
Compliance with 21 CFR 890?	Yes	Unknown	Yes	Unknown	Yes	Unknown	Yes	Unknown
Weight	440g	158g	440g	158g	440g	158g	440g	330g
Dimensions (WxHxD) in (mm)	175x95x30	114x50x31	175x95x30	114x50x31	175x95x30	114x50x31	175x95x30	102x152x51
Housing Material and Construction	Plastics	Plastics	Plastics	Plastics	Plastics	Plastics	Plastics	Plastics
Output Specifications								
Waveform	Biphasic symmetrical	Pulsed biphasic asymmetrical	Biphasic symmetrical	Pulsed biphasic asymmetrical	Biphasic symmetrical	Pulsed biphasic asymmetrical	N/A	N/A
Shape	Rectangular	Rectangular with discharge	Rectangular	Rectangular with discharge	Rectangular	Rectangular with discharge	N/A	N/A
Maximum Output Voltage (500Ω)	50V	50V	50V	50V	50V	50V	N/A	N/A
Maximum Output Voltage (2kΩ)	115V	150V	115V	150V	115V	150V	N/A	N/A
Maximum Output Voltage (10kΩ)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Maximum Output Current (500Ω)	100mA	100mA	100mA	100mA	100mA	100mA	N/A	N/A
Maximum Output Current (2kΩ)	58mA	80mA	58mA	80mA	58mA	80mA	N/A	N/A
Maximum Output Current (10kΩ)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Pulse Width	500µs	50-400µs	300µs	50-400µs	350-500µs	50-400µs	N/A	N/A
Frequency (Hz)	10Hz	1-10Hz	35Hz	20-60Hz	5-60Hz	1-60Hz	N/A	N/A
Net Charge (µC per pulse) (500Ω)	0µC	0µC	0µC	0µC	0µC	0µC	N/A	N/A
Maximum Phase Charge (µC) (500Ω)	50µC	40µC	30µC	40µC	60µC	40µC	N/A	N/A
Surface area of the electrodes (cm²)	Axial Guard (Surface) PR100 21.2cm²	Vaginal: 4.9cm² Anal: 3 cm²	Axial Guard (Surface) PR100 21.2cm²	Vaginal: 4.9cm² Anal: 3 cm²	Axial Guard (Surface) PR100 21.2cm²	Vaginal: 4.9cm² Anal: 3 cm²	N/A	N/A
Maximum Current Density (mA/cm²) (500Ω)	4.7 mA/cm²	Vaginal: 20.4mA/cm² Anal: 33.3mA/cm²	4.7 mA/cm²	Vaginal: 20.4mA/cm² Anal: 33.3mA/cm²	4.7 mA/cm²	Vaginal: 20.4mA/cm² Anal: 33.3mA/cm²	N/A	N/A
Maximum Power Density (W/cm²) (500Ω)	23.5 µW/cm²	Vaginal: 2.3mW/cm² Anal: 3.8mW/cm²	105 µW/cm²	Vaginal: 2.3mW/cm² Anal: 3.8mW/cm²	0.42 mW/cm²	Vaginal: 2.3mW/cm² Anal: 3.8mW/cm²	N/A	N/A
Burst Mode: Pulses per burst	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Burst Mode: Bursts per second	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Burst Mode: Burst duration (seconds)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Burst Mode: Duty Cycle	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
ON Time (sec.)	13s	5-10s	9-12s	2-20s	9-11s	2-20s	N/A	N/A
OFF Time (sec.)	3s	1-5s	7-9s	2-20s	3-8s	2-20s	N/A	N/A
Additional features	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 84 · e-mail: vertrieb.oustria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8884 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400

Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

TENS:

Characteristics/ Specifications	New Device	Predicate Device	Predicate Device	Predicate Device
Basic Unit Characteristics				
510(k) Number	K060950	K032954	K032954	K945301
Manufacturer	Otto Bock	Donimeter A/S	Donimeter A/S	Compez SA
Device Name, Model	STIMWELL med4	Elpho II1000	Elpho II3000	STIMWELL med4
Power Source(s)	Battery Pack Li-Ion 11.1V	Battery NMH or Alkaline 9V	Battery Pack NMH or Alkaline 9V	Battery Pack Li-Ion 11.1V
Method of Line Current Isolation	Medical Class II Power Adapter – Mascoi (12.6VDC-15.1W)	N/A	N/A	Medical Class II Power Adapter – Mascoi (12.6VDC-15.1W)
Patient Leakage Current (normal condition)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)
Patient Leakage Current (single fault condition)	N/A (Battery)	N/A (Battery)	N/A (Battery)	N/A (Battery)
Number of Output Modes	1	1	1	1
Number of Output Channels	4	2	2	4
Number of EMG (input) Channels	N/A	N/A	N/A	N/A
EMG sensitivity	N/A	N/A	N/A	N/A
EMG Sampling Rate	N/A	N/A	N/A	N/A
EMG detection (bipolar/ monopolar)	N/A	N/A	N/A	N/A
EMG range (µV)	N/A	N/A	N/A	N/A
EMG bandwidth	N/A	N/A	N/A	N/A
EMG signal processing (eg. RMS)	N/A	N/A	N/A	N/A
Synchronous or Alternating?	Alternating	Unknown	Unknown	Synchronous
Method of Channel Isolation	Transformer, Inductive couplers	Unknown	Unknown	Transformer
Regulated Current or Regulated Voltage?	Regulated Current	Regulated Current	Regulated Current	Regulated Current
Software/Firmware/ Microprocessor Control?	Yes	Yes	Yes	Yes
Automatic Overload Trip?	Yes	Yes	Yes	Yes
Automatic No-Load Trip?	Yes	Yes	Yes	Yes
Automatic Shut Off?	Yes (10min)	Yes	Yes	Unknown
Patient Override Control?	Yes (Stop Button)	Yes	Yes	Yes
Indicator Display: On/Off Status?	Yes	Yes	Yes	Yes
Indicator Display: Low Battery?	Yes	Yes	Yes	Yes
Indicator Display: Voltage/Current Level?	Yes	Yes	Yes	Yes
Timer Range (minutes)	10 – 120min	45min and cont.	5 – 95min and cont.	3 – 54min
Compliance with Voluntary Standards? (if yes, specify)	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown
Compliance with 21 CFR 908?	Yes	Unknown	Unknown	Unknown
Weight	440g	185g	158g	420g
Dimensions (Width/D) in [mm]	175x95x30	117x60x31	114x60x31	148x80x30
Housing Material and Construction	Plastics	Plastics	Plastics	Aluminum
Output Specifications				
Waveform (e.g., pulsed monophasic, biphasic)	Biphasic symmetrical	Biphasic asymmetrical	Biphasic asymmetrical	Biphasic symmetrical
Shape (e.g., rectangular, spike, rectified sinusoidal)	Rectangular	Rectangular with discharge	Rectangular with discharge	Rectangular
Maximum Output Voltage (500Ω)	50V	20V	50V	50V
Maximum Output Voltage (2kΩ)	115V	76V	164V	200V
Maximum Output Voltage (10kΩ)	N/A	Unknown	Unknown	N/A
Maximum Output Current (500Ω)	100mA	40mA	60mA	100mA
Maximum Output Current (2kΩ)	58mA	38mA	52mA	100mA
Maximum Output Current (10kΩ)	N/A	Unknown	Unknown	N/A
Pulse Width (specify units)	150µs, 200µs	150µs, 180µs, 200µs	50 – 400µs	Unknown
Frequency [Hz]	2 – 100Hz	2 – 120Hz	0.5 – 150Hz	Unknown
For interlaminar modes only: Beat Frequency [Hz]	N/A	N/A	N/A	N/A
For multiphasic waveforms only: Symmetrical phases?	N/A	N/A	N/A	N/A
For multiphasic waveforms only: Phase Duration (including units); state range, if applicable; both phases if asymmetrical	N/A	N/A	N/A	N/A
Net Charge (µC per pulse); (500Ω) // zero, state method of achieving zero net charge	0µC Same positive and negative impulse	0µC Output capacitor	0µC Output capacitor	0µC Same positive and negative impulse
Maximum Phase Charge (µC) (500Ω)	20µC	8µC	24µC	Unknown
Maximum Current Density (mA/cm²) (500Ω)	12.5mA/cm²	2.1mA/cm²	3.2mA/cm²	5.0mA/cm²
Maximum Power Density (W/cm²) (500Ω) using smallest electrode conductive surface area	1.0mW/cm²	0.09mW/cm²	1.4mW/cm²	Unknown
Burst Mode (i.e., pulse trains): Pulses per burst	8	7	7	Unknown
Burst Mode (i.e., pulse trains): Bursts per second	2	2	0.5 – 5	Unknown
Burst Mode (i.e., pulse trains): Burst duration (seconds)	100ms	87.5ms	87.5ms	Unknown
Burst Mode (i.e., pulse trains): Duty Cycle [Line (b) ± Line (a)]	20%	17.5%	4.4% – 44%	Unknown
ON Time (seconds)	Continuous or Burst	Continuous or Burst	Continuous or Burst	Continuous
OFF Time (seconds)	N/A	N/A	N/A	N/A
Additional features (if applicable)	N/A	N/A	N/A	N/A

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 • 1070 Wien • Austria • Telefon (+43-1) 523 37 86 • Telefax (+43-1) 523 22 64 • e-mail: vertrieb.austria@ottobock.com • www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 • Telefax (+43-1) 526 79 85

HG WIEN FN 242378p • UID-Nr.: ATU 57529204 • ARA-Nr.: 8864 • ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400

Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

Biofeedback:

Characteristics/ Specifications	New Device	Predicate Device	Predicate Device
Basic Unit Characteristics			
510(k) Number	K080950	K053434	K040849
Manufacturer	Otto Bock	Thought Technology Ltd.	Mentamove North America, LCC
Device Name, Model	STIWELL med4	Myotrac Infinity	Mentamove
Power Source(s)	Battery Pack Li-Ion 11.1V	Battery Pack NiMH rechargeable	Battery Pack NiMH rechargeable
Method of Line Current Isolation	Medical Class II Power Adapter – Mascot (12.6VDC-15, 1W)	Medical Class II Power Adapter (6VDC-15W)	Battery Charger (6VDC-2, 1W)
Patient Leakage Current (normal condition)	N/A (Battery)	N/A (Battery)	N/A (Battery)
Patient Leakage Current (single fault condition)	N/A (Battery)	N/A (Battery)	N/A (Battery)
Number of Output Modes	1	0	0
Number of Output Channels	1	0	0
Number of EMG (input) Channels	2	2	1
EMG sensitivity	1µV	<0, 1µV	1µV
EMG Sampling Rate	3kHz	2,048kHz	Unknown
EMG detection (bipolar/ monopolar)	Bipolar	Bipolar	Bipolar
EMG range (µV)	1-2000µV	0-2000µV	2-2000µV
EMG bandwidth	70-480Hz	10Hz-1kHz	Unknown
EMG signal processing (eg. RMS)	AVR (Average Rectified Value)	RMS (Root Mean Square)	Unknown
Synchronous or Alternating?	N/A	N/A	N/A
Method of Channel Isolation	N/A	N/A	N/A
Regulated Current or Regulated Voltage?	N/A	N/A	N/A
Software/Firmware/ Microprocessor Control?	Yes	Yes	Yes
Automatic Overload Trip?	N/A	N/A	N/A
Automatic No-Load Trip?	N/A	N/A	N/A
Automatic Shut Off?	Yes (10min)	Unknown	Unknown
Patient Override Control?	Yes (Stop Button)	Yes	Yes
Indicator Display: On/Off Status?	Yes	Yes	Yes
Indicator Display: Low Battery?	Yes	Yes	Yes
Indicator Display: Voltage/Current Level?	N/A	N/A	N/A
Timer Range (minutes)	5-30min	1-120min	Unknown
Compliance with Voluntary Standards? (if yes, specify)	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-10	Unknown	Unknown
Compliance with 21 CFR 898?	Yes	Unknown	Unknown
Weight	440g	330g	Unknown
Dimensions (WxHxD) in [mm]	175x95x30	102x152x51	Unknown
Housing Material and Construction	Plastics	Plastics	Plastics
Output Specifications			
Waveform (e.g., pulsed monophasic, biphasic)	N/A	N/A	N/A
Shape (e.g., rectangular, spike, rectified sinusoidal)	N/A	N/A	N/A
Maximum Output Voltage (500Ω)	N/A	N/A	N/A
Maximum Output Voltage (2kΩ)	N/A	N/A	N/A
Maximum Output Voltage (10kΩ)	N/A	N/A	N/A
Maximum Output Current (500Ω)	N/A	N/A	N/A
Maximum Output Current (2kΩ)	N/A	N/A	N/A
Maximum Output Current (10kΩ)	N/A	N/A	N/A
Pulse Width (specify units)	N/A	N/A	N/A
Frequency [Hz]	N/A	N/A	N/A
For interventional modes only: Beat Frequency [Hz]	N/A	N/A	N/A
For multiphasic waveforms only: Symmetrical phases?	N/A	N/A	N/A
For multiphasic waveforms only: Phase Duration (including units)	N/A	N/A	N/A
Net Charge (µC per pulse): (500Ω)	N/A	N/A	N/A
Maximum Phase Charge (µC): (500Ω)	N/A	N/A	N/A
Maximum Current Density [mA/cm²] (500Ω)	N/A	N/A	N/A
Maximum Power Density [W/cm²] (500Ω)	N/A	N/A	N/A
Burst Mode (i.e., pulse trains): Pulses per burst	N/A	N/A	N/A
Burst Mode (i.e., pulse trains): Bursts per second	N/A	N/A	N/A
Burst Mode (i.e., pulse trains): Burst duration (seconds)	N/A	N/A	N/A
Burst Mode (i.e., pulse trains): Duty Cycle [Line (b) x Line (c)]	N/A	N/A	N/A
ON Time [seconds]	N/A	N/A	N/A
OFF Time [seconds]	N/A	N/A	N/A
Additional features (if applicable)	N/A	N/A	N/A

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 84 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400

Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

G Functional and Safety Testing:

Tests have been carried out according to the following standards:

No.	Title	Version	Comments
IEC 60601-1	Medical electrical equipment Part 1: General requirements for safety	1996-03	Medical electrical devices (Testing)
IEC 60601-1-1	Medical electrical equipment Part 1-1: General requirements for safety – Safety requirements for medical electrical systems	2002-09-01	Medical electrical devices (Testing)
IEC 60601-1-2	Medical electrical equipment Part 1-2: General requirements for safety – Electromagnetic compatibility – requirements and tests	2001	EMC (Testing)
IEC 60601-1-4	Medical electrical equipment Part 1-4: General requirements for safety – Programmable electrical medical systems	2001-06-01	Programmable Systems
IEC 60601-2-10	Medical electrical equipment Part 2-10: Particular requirements for the safety of nerve and muscle stimulators	2003-04	Medical electrical devices; EMC (Testing)
IEC 721-3-2	Classification of Environmental Conditions Part 3-2: Classes of Influencing Factors - Transport	1981	Transport (Testing)
ISO 2248	Packaging - Complete, filled transport packages - Vertical impact test by dropping	1985	Transport (Testing)
ISO 2876	Packaging - Complete, filled transport packages – Rolling test	1985	Transport (Testing)
IEC 60512-8	Electromechanical components for electronic equipment; basic testing procedures and measuring methods Part 8: Connector test (mechanical) and mechanical test on contacts and terminations	3.0/1994-05	Specific sections were used for verification testing
IEC 60512-9	Electromechanical components for electronic equipment; basic testing procedures and measuring methods Part 9: Miscellaneous tests	2.0/1992-05	Specific sections were used for verification testing
IEC 68-1 +Corr.+A1	Environmental testing Part 1: General and Guidance	1998	Specific sections were used for verification testing
IEC 68-2-1 +A1+A2	Environmental testing Part 2-1: Tests, Tests A: Cold	1994	Specific sections were used for verification testing

Otto Bock Healthcare Products GmbHKaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 88 · Telefax (+43-1) 523 22 64 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400
Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

No.	Title	Version	Comments
IEC 68-2-2 +IEC68-2-2A+A2	Environmental testing Part 2-2: Tests, Tests B: Dry heat	1994	Specific sections were used for verification testing
IEC 68-2-30	Environmental testing - Part 2-30: Tests - Test Db: Damp heat, cyclic (12 h + 12 h cycle)	2000	Specific sections were used for verification testing
IEC 68-2-33	Environmental testing - Part 2-33: Tests. Guidance on change of temperature tests	2000	Specific sections were used for verification testing
IEC 68-2-36	Environmental testing Part 2-36, Tests Fdb: Vibrations, noise (broadband), repeatability: medium	1973	Specific sections were used for verification testing
IEC 68-2-78	Environmental testing Part 2-78: Tests, Test Cb: Damp heat, steady state, primarily intended for devices	2001	Specific sections were used for verification testing
21 CFR 898	Performance Standard for electrode lead wires and patient cables	1997	Electrical safety (Testing)

In addition to tests according to the above mentioned standards several functional and safety tests defined by the manufacturer have been conducted.

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 64 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400
Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAATWW BLZ 20 111, Konto-Nr. 30001045517



QUALITY FOR LIFE

H Indications for use:

The STIWELL med4 is a neuromuscular electronic stimulator indicated for use under medical supervision for adjunctive therapy in the treatment of medical diseases and conditions.

As a powered muscle stimulator the STIWELL med4 is indicated for the following conditions:

- Relaxation of muscle spasm
- Prevention or retardation of disuse atrophy
- Increasing local blood circulation
- Muscle re-education
- Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis
- Maintaining or increasing range of motion

As a transcutaneous electrical nerve stimulator for pain relief the STIWELL med4 is indicated for the following conditions:

- Symptomatic relief and management of chronic (long-term), intractable pain
- Adjunctive treatment in the management of post-surgical pain and post traumatic acute pain

As a biofeedback device the STIWELL med4 is indicated for the following conditions:

- Biofeedback, relaxation and muscle re-education purposes

As an external functional neuromuscular stimulator the STIWELL med4 is indicated for the following conditions:

- Helps to relearn voluntary motor functions of the extremities

As a nonimplanted electrical continence device the STIWELL med4 is indicated for the following conditions:

- Acute and ongoing treatment of stress, urge or mixed urinary incontinence and where the following results may improve urinary control: Inhibition of the detrusor muscles through reflexive mechanisms and strengthening of pelvic floor muscles
- Incontinence treatment for assessing EMG activity of the pelvic floor and accessory muscles such as the abdominal and the gluteus muscles

I Conclusion:

The STIWELL med4 stimulation device is safe and effective for its intended use. The STIWELL med4 is substantially equivalent to the predicate devices.

J Date summary submitted:

21st January 2008

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien · Austria · Telefon (+43-1) 523 37 86 · Telefax (+43-1) 523 22 64 · e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Kundenservice: Telefon (+43-1) 526 95 48 · Telefax (+43-1) 526 79 85

HG WIEN FN 242378p · UID-Nr.: ATU 57529204 · ARA-Nr.: 8864 · ERA-Nr.: 50411

Bank Austria Creditanstalt Aktiengesellschaft, Filiale Westbahnstraße Swift Code: BKAUATWW BLZ 12 000, Konto-Nr. 09743146400
Erste Bank, Kommerzcenter Wien West 2 Swift Code: GIBAAATWW BLZ 20 111, Konto-Nr. 30001045517



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Otto Bock, Austria GmbH
% W.F. Jackson Associates, Ltd
Mr. William Jackson
2247 Jennifer Lane
North St. Paul, Minnesota 55109-2851

Re: K080950

APR - 2 2009

Trade/Device Name: STIWELL med4
Regulation Number: 21 CFR 890.5850
Regulation Name: Powered muscle stimulator
Regulatory Class: II
Product Code: IPF, GZI, KPI, GZJ, HCC
Dated: March 2, 2009
Received: March 25, 2009

Dear Mr. Jackson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometrics' (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a long horizontal flourish extending to the right.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K080950

Device Name: STIWELL med4

Indications for Use:

The STIWELL med4 is a neuromuscular electronic stimulator indicated for use under medical supervision for adjunctive therapy in the treatment of medical diseases and conditions.

As a powered muscle stimulator the STIWELL med4 is indicated for the following conditions:

- Relaxation of muscle spasm
- Prevention or retardation of disuse atrophy
- Increasing local blood circulation
- Muscle re-education
- Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis
- Maintaining or increasing range of motion

As a transcutaneous electrical nerve stimulator for pain relief the STIWELL med4 is indicated for the following conditions:

- Symptomatic relief and management of chronic (long-term), intractable pain
- Adjunctive treatment in the management of post-surgical pain and post traumatic acute pain

As a biofeedback device the STIWELL med4 is indicated for the following conditions:

- Biofeedback, relaxation and muscle re-education purposes

As an external functional neuromuscular stimulator the STIWELL med4 is indicated for the following conditions:

- Helps to relearn voluntary motor functions of the extremities

As a nonimplanted electrical continence device the STIWELL med4 is indicated for the following conditions:

- Acute and ongoing treatment of stress, urge or mixed urinary incontinence and where the following results may improve urinary control: Inhibition of the detruser muscles through reflexive mechanisms and strengthening of pelvic floor muscles
- Incontinence treatment for assessing EMG activity of the pelvic floor and accessory muscles such as the abdominal and the gluteus muscles

Prescription Use X
(Per 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(Per 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDHR, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of General, Restorative,
and Neurological Devices

Page 1 of 1

510(k) Number

K080950