



JAN 21 2009

**510(k) Summary
COULTER 6C Cell Control**

1 Submitted By:

Lourdes Coba
Staff Regulatory Affairs Specialist
Beckman Coulter, Inc
11800 SW 147 Avenue, M/C 31-B06
Miami, Florida 33196-2500
Telephone (305) 380-4191
FAX (305) 380-4344

K081822

2 Date Submitted:

June 26, 2008

3 Device Name(s):

3 1 Proprietary Names

COULTER® 6C Cell Control

3 2 Classification Name

Hematology quality control mixture
(21 CFR § 864.8625)

4 Predicate Device:

Candidate(s)	Predicate	Manufacturer	Docket Number
COULTER® 6C Cell Control	COULTER® 5C® Cell Control (Cleared as COULTER® PX Cell Control)	Beckman Coulter, Inc	K912133 & K060464

5 **Description:**

6C Cell Control is a reference product prepared from treated, stabilized human erythrocytes in an isotonic medium. 6C Cell Control also contains a stabilized, platelet-sized component, and fixed erythrocytes to simulate leukocytes and nucleated red blood cells. By design, 6C Cell Control confirms and monitors instrument accuracy and precision performance by providing measurements for counting, sizing, hemoglobin determination, NRBC enumeration and White Blood Cell differentiation using VCSn technology.

6 **Intended Use:**

6C Cell Control is a hematology quality control material used to monitor the performance of COULTER hematology analyzers listed in the TABLE OF EXPECTED RESULTS in conjunction with specific COULTER reagents.

The assigned values and expected ranges on the TABLE OF EXPECTED RESULTS can be used to monitor instrument performance. This product can also be used to establish your own laboratory mean.

7 **Comparison to Predicate(s):**

COULTER 6C Cell Control is essentially identical to the current COULTER 5C Cell Control with the addition of 2 parameters: NRBC% and NRBC#. This control is for use with the UniCel® DxH 800 analyzer (pending 510(k) submission). An NRBC analog was added to the 5C control product formulation.

8 **Summary of Performance Data:**

The data in the Premarket Notification on safety and effectiveness supports a finding of substantial equivalence to products already in commercial distribution. Stability studies of 6C Cell Control support the Beckman Coulter stability claims of 18 events within 16 days (open vial) and 95 days (closed vial).

This summary of safety and effectiveness is being submitted in accordance with the requirements of the Safe Medical Device Act of 1990 and the implementing regulation 21 CFR 807.92.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Beckman Coulter, Inc
c/o Ms Nancy Nadler
Staff Regulatory Affairs Specialist
11800 SW 147th Avenue
Miami, FL 33196

JAN 21 2009

Re k081822

Trade/Device Name Coulter® 6C Cell Control
Regulation Number 21 CFR 864.8625
Regulation Name Hematology Quality Control Mixture
Regulatory Class Class II
Product Code JPK
Dated January 05, 2009
Received January 06, 2009

Dear Ms Nadler

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to, registration and listing (21 CFR Part 807), labeling (21 CFR Parts 801 and 809), and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometrics (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,



Maria M. Chan, Ph.D.
Acting Division Director
Division of Immunology and Hematology Devices
Office of In Vitro Diagnostic Device Evaluation
and Safety
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K081822

Device Name COULTER® 6C Cell Control

Indications For Use

6C Cell Control is a hematology quality control material used to monitor the performance of COULTER hematology analyzers listed in the TABLE OF EXPECTED RESULTS in conjunction with specific COULTER reagents

The assigned values and expected ranges on the TABLE OF EXPECTED RESULTS can be used to monitor instrument performance. This product can also be used to establish your own laboratory mean.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH Office of In Vitro Diagnostic Devices (OIVD)

Division Sign-Off

Office of In Vitro Diagnostic Device
Evaluation and Safety

Page 1 of 1

510(k)

K081822