



3150 NW 107<sup>th</sup> Avenue, Miami FL 33172  
Tel: 305.599.7174  
Fax: 305.592.4621

1092246

**510(k) Summary**  
**NIPRO Needleless Transfer Device**

**NOV 12 2009**

**807.92(a) (1)**

Contact Person: Jessica Oswald, Regulatory Affairs Specialist  
Date of summary preparation: July 22, 2009

**807.92(a) (2)**

Trade Name: NIPRO Needleless Transfer Device  
Common Name: TRANSFER NEEDLE  
Classification Name: Set, IV, Fluid transfer  
Regulation Number: 21 CFR part 880.5440  
Regulatory Class: II  
Panel: 80  
Product Code: LHI

**807.92(a) (3)**

Legally marketed substantial equivalent device:  
K972117: MSI Transflow

**807.92(a) (4)**

Description of device:

The Nipro Needleless Transfer Device with or without 20mm vial closure adapters is a sterile, single-use needleless liquid transfer device with the intent of transferring sterile water from its vial into an evacuated vial containing lyophilized product requiring reconstitution.

This device consists of a single piece injection molded main body with 2 capped spikes on both ends. The device is provided with and without 2 adapters used for setting at the mouth of the solution bottle in order to stabilize the container at the time of mixture.

The device is sterile, single-use only, latex-free, non-toxic and non-pyrogenic.

**807.92(a) (5)**

Indications for Use:

The NIPRO Needleless Transfer Device is a sterile device for the aseptic transfer of solutions from container to container or reconstituting medications.

**807.92(a) (6)**

Comparison of technological characteristics:

The NIPRO Needleless Transfer Device is substantially equivalent to the predicate device in the following technological characteristics –

- Design
- Physical characteristics
- Basic Scientific Technology
- Intended Use

**807.92(b) (1)**

Non-clinical tests performed and included in this submission include:

- Dimensional
- Mechanical
- Operational
- Chemical
- Biocompatibility

**807.92(b) (3)**

Conclusions drawn from non-clinical and clinical tests:

The results of the non-clinical tests and the comparison of technological characteristics with the predicate device demonstrate that the NIPRO Needleless Transfer Device performs equivalent to the predicate device and is safe and effective when used as intended.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

Ms. Jessica Oswald  
Regulatory Affairs Specialist  
Nipro Medical Corporation  
3150 NW 107<sup>th</sup> Avenue  
Miami, Florida 33172

NOV 12 2009

Re: K092246  
Trade/Device Name: NIPRO Needleless Transfer Device  
Regulation Number: 21CFR 880.5440  
Regulation Name: Intravascular Administration Set  
Regulatory Class: II  
Product Code: LHI  
Dated: October 6, 2009  
Received: October 9, 2009

Dear Ms Oswald:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page – Ms. Oswald

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Susan Runner".

Susan Runner, D.D.S., M.A.  
Acting Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known):

Device Name: NIPRO Needleless Transfer Device

Indications For Use:

The NIPRO Needleless Transfer Device is a sterile device for the aseptic transfer of solutions from container to container or reconstituting medications.

Prescription Use √

(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

NIPRO Needleless Transfer Device  
4 Indications for Use

510(k) Number: K092246

-4.1-