

1092320

## 510(k) Summary—Elecsys PreciControl ThyroAB

NOV - 3 2009

---

|                     |   |
|---------------------|---|
| <b>Introduction</b> | In accordance with 21 CFR 807.92, Roche Diagnostics hereby submits official notification as required by Section 510(k) of the Federal Food, Drug and Cosmetics Act of our intention to market the device described in this Premarket Notification [510(k)]. |
|---------------------|---|

---

|   |  |
|---|--|
| <b>Submitter<br/>Name, Address,<br/>Contact</b> | <p>Roche Diagnostics<br/>9115 Hague Road<br/>Indianapolis, IN 46250</p> <p>Contact Person: Sarah Baumann<br/>Phone: 317-521-3952<br/>Fax: 317-521-2324<br/>Email: <a href="mailto:sarah.baumann@roche.com">sarah.baumann@roche.com</a></p> <p>Secondary Contact: Stephanie Greeman<br/>Phone: 317-521-2458<br/>Fax: 317-521-2324<br/>Email: <a href="mailto:stephanie.greeman@roche.com">stephanie.greeman@roche.com</a></p> <p>Date Prepared: July 30, 2009</p> |
|---|--|

---

|                               |   |
|-------------------------------|---|
| <b>Submission<br/>Purpose</b> | <p>PreciControl ThyroAB is used for quality control of specified Elecsys immunoassays on the Elecsys and cobas e immunoassay analyzers. This product contains control material for numerous Elecsys assays in one convenient solution.</p> <p>Changes to PreciControl ThyroAB consist of the addition of anti-thyroperoxidase (Anti-TPO) and anti-thyroglobulin (Anti-TG) antibodies to extend the current functionality.</p> |
|-------------------------------|---|

---

|                    |  |
|--------------------|--|
| <b>Device Name</b> | <p>Proprietary name: Elecsys PreciControl ThyroAB<br/>Common name: PreciControl ThyroAB<br/>Classification name: Multi-Analyte Controls, All Kinds (assayed and Unassayed)</p> |
|--------------------|--|

---

*Continued on next page*

## 510(k) Summary—Elecsys PreciControl ThyroAB, Continued

---

|                           |  |
|---------------------------|--|
| <b>Device Description</b> | The Elecsys PreciControl ThyroAB is a lyophilized product consisting of antibodies in a human serum matrix. During manufacture, the antibodies are spiked into the matrix at the desired concentration levels. |
|---------------------------|--|

---

|                     |  |
|---------------------|--|
| <b>Intended Use</b> | Elecsys PreciControl ThyroAB is used for quality control of the Elecsys Anti-TSHR, Anti-TPO and Anti-Tg immunoassays on the Elecsys and cobas e immunoassay analyzers. |
|---------------------|--|

---

|                         |   |
|-------------------------|---|
| <b>Predicate Device</b> | The modified Elecsys PreciControl ThyroAB is substantially equivalent to other products in commercial distribution intended for similar use. We claim equivalency to the currently marketed Elecsys PreciControl ThyroAB (K080092). |
|-------------------------|---|

---

|                                       |   |
|---------------------------------------|---|
| <b>Device Comparison—Similarities</b> | The table below presents the similarities between the modified Elecsys PreciControl ThyroAB and the predicate device, Elecsys PreciControl ThyroAB (K080092). |
|---------------------------------------|---|

| Characteristic  | Predicate Device<br>Elecsys PreciControl<br>ThyroAB<br>(K080092)   | Elecsys PreciControl<br>ThyroAB |
|-----------------|--|---------------------------------|
| Analyzer System | Elecsys and cobas e immunoassay analyzers:<br>-Elecsys 2010<br>-MODULAR ANALYTICS E170<br>-cobas e 411<br>-cobas e 601 | Same                            |
| Format          | Lyophilized  | Same                            |
| Matrix          | Human Serum  | Same                            |

*Continued on next page*

## 510(k) Summary—Elecsys PreciControl ThyroAB, Continued

**Device  
Comparison—  
Differences**

The table below presents the differences between the modified Elecsys PreciControl ThyroAB and the predicate device, Elecsys PreciControl ThyroAB (K080092).

| Characteristic           | Predicate Device<br>Elecsys PreciControl<br>ThyroAB<br>(K080092)   | Elecsys PreciControl<br>ThyroAB   |
|--------------------------|--|---|
| Intended use             | Elecsys PreciControl ThyroAB is used for quality control of the Elecsys Anti-TSHR immunoassay on the Elecsys and cobas e immunoassay analyzers.  | Elecsys PreciControl ThyroAB is used for quality control of the Elecsys Anti-TSHR, Anti-TPO and Anti-Tg immunoassays on the Elecsys and <b>cobas e</b> immunoassay analyzers.                                   |
| Analyte concentration    | Anti-TSHR (IU/L):<br>Level 1 = 4<br>Level 2 = 16   | Anti-TSHR (IU/L):<br>Level 1 = 4<br>Level 2 = 16<br><br>Anti-TPO (IU/mL):<br>Level 1 = 35<br>Level 2 = 100<br><br>Anti-TG (IU/mL):<br>Level 1 = 100<br>Level 2 = 200  |
| Antibody source and type | Anti-TSHR: Human monoclonal  | Anti-TSHR: Human monoclonal<br>Anti-TPO: Sheep polyclonal<br>Anti-TG: Sheep polyclonal  |
| Handling                 | Dissolve carefully the contents of one bottle by adding exactly 2.0 mL of distilled water and allow to stand closed for 15 minutes to reconstitute. Mix carefully, avoiding the formation of foam. | Dissolve carefully the contents of one bottle by adding exactly 2.0 mL of distilled or deionized water and allow to stand closed for 30 minutes to reconstitute. Mix carefully, avoiding the formation of foam. |

*Continued on next page*

## 510(k) Summary—Elecsys PreciControl ThyroAB, Continued

**Device Comparison—Differences** (continued) The table below presents the differences between the modified Elecsys PreciControl ThyroAB and the predicate device, Elecsys PreciControl ThyroAB (K080092).

| Characteristic | Predicate Device<br>Elecsys PreciControl<br>ThyroAB<br>(K080092)   | Elecsys PreciControl<br>ThyroAB  |
|----------------|--|--|
| Stability      | <u>Unopened</u> :<br>Store at 2-8°C until expiration date<br><br><u>Reconstituted</u> :<br>on the analyzer at 20-25°C:<br>up to 3 hrs<br>at -20°C: 3 months (freeze only once)<br><br><u>After Thawing</u> :<br>use only once. | <u>Unopened</u> :<br>Store at 2-8°C until expiration date<br><br><u>Reconstituted</u> :<br>on the analyzer at 20-25 °C:<br>up to 5 hrs<br>at -20°C: 1 month (freeze only once)<br>or at 2-8°C for 3 day<br>(for Anti-TG & Anti-TPO only) |

**Performance Characteristics** The Elecsys PreciControl ThyroAB was evaluated for value assignment, stability, and duration of reconstitution.



## DEPARTMENT OF HEALTH & HUMAN SERVICES

---

Roche Diagnostics  
c/o Ms. Sarah Baumann  
Regulatory Affairs Consultant  
9115 Hague Road, PO Box 50410  
Indianapolis, Indiana 46250-4016

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Mail Center – WO66-0609  
Silver Spring, MD 20993-0002

NOV - 3 2009

Re: k092320  
Trade Name: Elecsys PreciControl-ThyroAB  
Regulation Number: 21 CFR §862.1660  
Regulation Name: Quality Control Material (assayed and unassayed).  
Regulatory Class: Class I, reserved  
Product Codes: JJY  
Dated: September 30, 2009  
Received: October 01, 2009

Dear Ms. Baumann:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (301) 796-5760. For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-5680 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to be 'CCH', with a long horizontal flourish extending to the right.

Courtney C. Harper, Ph.D.  
Director  
Division of Chemistry and Toxicology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Indication for Use

510(k) Number (if known): K092320

Device Name: Elecsys PreciControl ThyroAB

Indication For Use:

Elecsys PreciControl ThyroAB is used for quality control of the Elecsys Anti-TSHR, Anti-TPO and Anti-Tg immunoassays on the Elecsys and **cobas e** immunoassay analyzers.


Prescription Use   X    
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use         
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)

  
\_\_\_\_\_  
Division Sign-Off  
Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k)\_k092320 K092320