

DEC 18 2009

510(K) SUMMARY
(as required by 807.92(c))

Submitter of 510(k): Dynarex Corporation
10 Glenshaw Street
Orangeburg, NY 10972
USA
Phone: (845) 365-8200
Fax: (845) 365-8201

Contact Person: Daniel Consaga

Date of Summary: December 17, 2009

Trade/Proprietary Name: Dynarex Sterile Lubricating Jelly

Classification Name: Lubricant, Patient

Product Code: KMJ

Intended Use:

Dynarex Sterile Lubricating Jelly is a medical device intended for medical purposes, to lubricate body orifices to facilitate entry of diagnostic or therapeutic devices when a sterile field is required.

Device Description:

Dynarex Sterile Lubricating Jelly is a clear, greaseless, water soluble jelly. Typical packaging configurations for the Dynarex Sterile Lubricating Jelly is 3 g or 5 g foil packs, and 2 oz or 4 oz tubes. Other sizes may become available.

Predicate Device:

E-Z Lubricating Jelly, 510k Number K041060 is manufactured by Chester Labs, Inc. 1900 Section Rd., Cincinnati, OH 45237.

Substantial Equivalence:

The Dynarex Sterile Lubricating Jelly provides effective lubrication during the insertion of diagnostic and therapeutic devices into body orifices. Its function and performance are similar to the predicate device as presented in this 510(k)

Safety and Effectiveness of the device:

This device is as safe and effective as the predicate device cited above based on .

510(K) SUMMARY (as required by 807.92(c))

Safety and Effectiveness of the device:

This device is as safe and effective as the predicate device cited above based on the following;

Summary comparing technological characteristics with other predicate device:

Dynarex Sterile Lubricating Jelly is similar in terms of intended use and technological characteristics to predicate devices reviewed as a lubricating jelly to facilitate the entry of a diagnostic or therapeutic device. The device is substantially equivalent with respect to indications for use and other physical characteristics to predicate devices in terms of section 510(k) substantial equivalency.

Please find below a tabulated comparison supporting that this device is substantially equivalent to the predicate device in commercial distribution.

TECHNOLOGICAL CHARACTERISTICS	Dynarex Lubricating Jelly	E-Z Lubricating Jelly
Purified water	Yes	Yes
Carbomer thickeners	Yes	Yes
Methylparaben	Yes	Yes
Labeled water soluble	Yes	Yes
Labeled non-staining	Yes	Yes
Labeled alcohol and fragrance free	Yes	Yes
Container material	Plastic/Film Laminate	Plastic
Sterile	Yes	Yes
Physical Tests	Viscosity & pH	Viscosity & pH
	Preservative effectiveness	Preservative effectiveness
	Sterility Test Study	Sterility Test Study
Biocompatibility Testing	ISO 10993	ISO 10993
	In-Vitro Cytotoxicity - Pass*	In-Vitro Cytotoxicity - Pass
	Implantation - Pass	Implantation - Pass
	Irritation & Hypersensitivity - Pass	Irritation & Hypersensitivity - Pass
	Systemic Toxicity - Pass	Systemic Toxicity - Pass
In-Vitro Cytotoxicity - Pass* Grade 2 Result Response	The two preservatives (Methylparaben and Propylparaben) are well known to cause this type result. Sodium Hydroxide also can cause this phenomena, but the concentration in the gel of sodium hydroxide is very low (less than 1%) so it is believed that its' contribution to this grade 2 result is very low if any.	

Sterilization:

The Dynarex Sterile Lubricating Jelly is sterilized by gamma irradiation under parameters that have been validated according to ISO/AAMI 11137 requirements (Sterilization of health care products - Requirements for validation and routine control - radiation sterilization).

Method: Vdmax²⁵ with an SAL of 10⁻⁶.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room W-066-0609
Silver Spring, MD 20993-0002

Mr. Daniel Consaga
Regulatory Affairs Engineer
Dynarex Corporation
10 Glenshaw Street
Orangeburg, New York 10962

DEC 18 2009

Re: K092488
Trade/Device Name: Dynarex Sterile Lubricating Jelly
Regulation Number: 21 CFR 880.6375
Regulation Name: Patient Lubricant
Regulatory Class: I
Product Code: KMJ
Dated: October 26, 2009
Received: December 8, 2009

Dear Mr. Consaga:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

A handwritten signature in black ink, appearing to read "Anthony D. Watson". The signature is stylized with a large, looped "A" and a cursive "Watson".

Anthony D. Watson, B.S., M.S., M.B.A.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

Device Name: Dynarex Sterile Lubricating Jelly

Indications for Use:

Dynarex Sterile Lubricating Jelly is a medical device intended for medical purposes, to lubricate body orifices to facilitate entry of diagnostic or therapeutic devices when a sterile field is required.


Prescription Use _____
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use X _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K 092488