

510(k) Summary

OCT - 5 2010

Introduction	According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.
Submitter name, address, contact	<p>Roche Diagnostics 9115 Hague Road, P.O. Box 50416 Indianapolis, IN 46250-0416 317-521-3577</p> <p>Contact Person: Kelly Colleen O'Maine Adams Phone: 317-521-3577 Fax: 317-521-2324 Email: colleen.adams@roche.com</p> <p>Secondary Contact: Stephanie Greeman Phone: 317-521-2458 Fax: 317-521-2324 Email: stephanie.greeman@roche.com</p> <p>Date Prepared: September 27, 2010</p>
Device Name	<p>Proprietary name: Elecsys CA 19-9 CalCheck 5 Common name: CA 19-9 CalCheck 5 Classification name: Single (specified) analyte controls (assayed and unassayed)</p>
Predicate device	The Elecsys CA 19-9 CalCheck 5 is substantially equivalent to other products in commercial distribution intended for similar use. We claim equivalency to the currently marketed Elecsys C-Peptide CalCheck 5 (K100810) and CA 19-9 CalCheck (K051185).
Device Description	The Elecsys CA 19-9 CalCheck 5 is a lyophilized product consisting of human CA 19-9 in human serum matrix. During manufacture, the analyte is spiked into the matrix at the desired concentration levels.
Intended use	The Elecsys CA 19-9 CalCheck 5 is an assayed control for use in calibration verification and for use in the verification of the assay range established by the Elecsys CA 19-9 reagent on the indicated Elecsys and cobas e immunoassay analyzers. For in vitro diagnostic use only.

Continued on next page

510(k) Summary, Continued

Comparison Table

The table below compares Elecsys CA 19-9 CalCheck 5 with the predicate devices, Elecsys C-Peptide CalCheck 5 (K100810) and CA 19-9 CalCheck (K051185). The first predicate shows that CA 19-9 CalCheck 5 is substantially equivalent to another CalCheck 5 product, Elecsys C-Peptide CalCheck 5. The Elecsys CA 19-9 CalCheck 5 is also substantially equivalent to the second predicate, CA 19-9 CalCheck, with several key similarities, especially the analyte. The shaded fields indicate similar characteristics between the candidate device and one or both predicate devices.

Characteristic	Elecsys C-Peptide CalCheck 5 (K100810)	Elecsys CA 19-9 CalCheck 5 (Candidate Device)	Elecsys CA 19-9 CalCheck (K051185)
Intended Use	The Elecsys C-Peptide CalCheck 5 is an assayed control for use in calibration verification and for use in the verification of the serum and plasma assay range established by the Elecsys C-Peptide reagent on the indicated Elecsys and cobas e immunoassay analyzers.	The Elecsys CA 19-9 Calcheck 5 is an assayed control for use in calibration verification and for use in the verification of the assay range established by the Elecsys CA 19-9 reagent on the indicated Elecsys and cobas e immunoassay analyzers. For in vitro diagnostic use only.	For use in the verification of the calibration established by the Elecsys CA 19-9 reagent on the indicated Elecsys and cobas e immunoassay analyzers.
Analyte	C-Peptide	CA 19-9	CA 19-9
Levels	Five	Five	Three
Assay Measuring Range	0.1 – 40 ng/mL	0.600 – 1000 U/mL	0.600 – 1000 U/mL
Check Target Values	Check 1: <0.2 ng/mL Check 2: 5.0 ng/mL Check 3: 20 ng/mL Check 4: 30 ng/mL Check 5: 40 ng/mL	Check 1: ≤ 3U/mL Check 2: 35 U/mL Check 3: 500 U/mL Check 4: 800 U/mL Check 5: 1000 U/mL	Check 1: < 5 U/mL Check 2: 179 U/mL Check 3: 676 U/mL
Format	Lyophilized	Lyophilized	Lyophilized
Handling	Reconstitute Check 1, Check 2, Check 3, Check 4, and Check 5 with exactly 1.0 mL distilled or deionized water. Allow to stand closed for 15 minutes, then mix gently by inversion.	Reconstitute Check 1, Check 2, Check 3, Check 4, and Check 5 with exactly 1.0 mL distilled or deionized water. Allow to stand closed for 15 minutes, then mix gently by inversion.	Reconstitute Check 1, Check 2, and Check 3 with exactly 1.0 mL distilled or deionized water. Allow to stand closed for 15 minutes, then mix gently by inversion.
Stability	<u>Unopened:</u> • Store at 2-8°C until expiration date <u>Reconstituted:</u> • 20-25°C: 4 hours	<u>Unopened:</u> • Store at 2-8°C until expiration date <u>Reconstituted:</u> • 20-25°C: 4 hours	<u>Unopened:</u> • Store at 2-8°C until expiration date <u>Reconstituted:</u> • 20-25°C: 4 hours
Matrix	Equine serum matrix	Human serum matrix	Human serum matrix

Continued on next page

**Performance
Characteristics**

The Elecsys CA 19-9 CalCheck 5 was evaluated for value assignment and stability.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food & Drug Administration
10903 New Hampshire Avenue
Building 66
Silver Spring, MD 20993

Roche Diagnostics
c/o Ms. Kelly Colleen O'Maine Adams, MTSC
Regulatory Affairs Consultant
9115 Hague Road
P.O.Box 50416
Indianapolis, IN, 46250-0416

OCT 05 2010

Re: k101365

Trade/Device Name: Elecsys CA 19-9 CalCheck 5
Regulation Number: 21 CFR§862.1660
Regulation Name: Quality Control Material, Assayed and Unassayed
Regulatory Class: Class I (Reserved)
Product Code: JJX
Dated: August 27, 2010
Received: August 31, 2010

Dear Ms. Adams:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



for

Maria M. Chan, Ph.D.

Director

Division of Immunology and Hematology Devices

Office of *In Vitro* Diagnostic Device Evaluation and Safety

Center for Devices and Radiological Health

Enclosure

Indication for Use

OCT - 5 2010

510(k) Number (if known): K101365

Device Name: Elecsys CA 19-9 CalCheck 5

Indication For Use:

The Elecsys CA 19-9 CalCheck 5 is an assayed control for use in calibration verification and for use in the verification of the assay range established by the Elecsys CA 19-9 reagent on the indicated Elecsys and **cobas e** immunoassay analyzers. For in vitro diagnostic use only.

Prescription Use X
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use ____
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)

Deena Philip

Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K101365