K102201

510(k) SUMMARY

DEC 1 5 2010

Summary Date:

December 4, 2010

Submitter Data:

Orthofix Inc.

1720 Bray Central Drive McKinney, TX 75609

214-937-2000 214-937-2736 (fax)

Primary Contact:

Daria Chew

darlachew@orthofix.com

Device Trade Name: Contours Lapidus Plating System

Common Name:

bone plate

Classification Name: Single/multiple component metallic bone fixation appliances and

accessories. (21 CFR Parts 888,3030

Product Code:

HRS

Legally Marketed

Predicate Device:

Darco ® Locking Plate System: K061808 / August 17, 2006

Device Description:

The Contours Lapidus Plating System consists of bone plates, screws and manual, orthopedic, surgical instruments. The plates (left and right versions) are low-profile and anatomically contoured for the 1st metatarsal and cuneiform. Each plate has a trapezoidal screw hole geometry with three angled, threaded holes at both the

distal and proximal ends as well as a compression hole for

intraoperative compression across the joint. The Contours Lapidus plate and bone screws are made from titanium alloy conforming to

ASTM F136.

Indications for Use:

The Contours Lapidus Plating System is intended for revision

procedures and joint fusion in the small bones of the foot.

Biomechanical Testing:

In order to demonstrate that the Contours Lapidus Plating System has the mechanical properties necessary to perform its intended use and to perform as well as the predicate device, Orthofix has conducted mechanical and functional testing of the Contours Lapidus Plating System. This testing includes bending strength, bending stiffness, fatigue life, screw torque and cadaveric functional testing. The results of the testing demonstrated the Contours Lapidus Plating System to meet or exceed all testing requirements and to perform as well as the predicate device.

1<102201

p-2/2

Technological Characteristics:

The Contours Lapidus Plating System is considered to be substantially equivalent in design, intended use and material to the predicate device. However, there are certain design differences, but these do not raise new questions regarding safety and effectiveness. Refer to the Table of Technological

Characteristics for a summary.

Features	Contours Lapidus Plating System	
Indications for Use	The Contours Lapidus Plating System is indicated for revision procedures and joint fusion in the small bones of the foot.	
Plate/Screw Material	Implant grade Titanium (Ti6A14V ELI)	
Plate Size	50mm x 16mm x 2mm (left and right versions)	
Plate Geometry	Low-profile, anatomically shaped for 1 st metatarsal and cuneiform bones.	
Fixation Method, Screw Holes	Trapezoidal screw hole placement. Three in proximal end, three in distal end plus compression hole.	
Screw Type	Low profile, locking, non-locking and compression	
Screw Length	Locking and non-locking screws: 12mm to 32mm lengths (in 2mm increments); compression screws are offered in sizes ranging from 16mm to 20mm lengths (in 2mm increments), 25mm to 40mm lengths in 5mm increments, and 42mm to 50mm lengths (in 2 mm increments)	

Sterilization:

The Contours Lapidus Plating System components are supplied NON-STERILE and require sterilization prior to use.

Substantial Equivalence:

Substantial equivalence is based upon design, dimension, material characterization, and biomechanical testing of the device in comparison to the predicate. The Contours Lapidus Plating System is substantially equivalent in design and function to the Darco Locking Plate System which was cleared under K061808 /

August 17, 2006.







Food and Drug Administration 10903 New Hampshire Avenue Document Control Room --WO66-G609 Silver Spring, MD 20993-0002

DEC 1 5 2010

Orthofix, Inc. % Ms. Mary E. Biggers, RAC Regulatory Affairs Consultant 1720 Bray Central Drive McKinney, Texas 75069

Re: K102201

Trade/Device Name: Contours Lapidus Plating System (bone plate)

Regulation Number: 21 CFR 888.3030

Regulation Name: Single/multiple component metallic bone fixation appliances and

accessories

Regulatory Class: Class II

Product Code: HRS

Dated: December 4, 2010 Received: December 8, 2010

Dear Ms. Biggers:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm.

Sincerely yours,

Mark N. Melkerson

Director

Division of Surgical, Orthopedic and Restorative Devices Office of Device Evaluation Center for Devices and Radiological Health

AB 2. Rho

Enclosure

K102201

p-1/1

INDICATION FOR USE STATEMENT

	Page of
510(k) Number (if known):	
Device Name:	Contours Lapidus Plating System (bone plate)
Indications for Use:	•
The Contours Lapidus Plating S joint fusion in the small bones o	System is intended for revision procedures and for the foot.
Prescription Use: X (Per 21 CFR 801 109)	Or Over-The-Counter (Optional Format 1-2-96)
(PLEASE DO NOT WRITE BELO NEEDED)	W THIS LINE-CONTINUE ON ANOTHER PAGE IF
Concurrence of CD	RH, Office of Device Evaluation (ODE)
Division o	Sign-Off) Surgical, Orthopedic, rative Devices

510(k) Number K102201

PAGE 13