

## 510(k) SUMMARY

DEC 16 2010

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR § 807.92.

The Assigned 510(k) Number is: K102419

Date: December 9, 2010

**Submitted by:** Wallac Oy, subsidiary of PerkinElmer  
940 Winter Street  
Waltham, MA 02451 USA

**Contact Person:** Susan K. Hamann  
**Primary:** Tel: 781-663-5872  
Fax: 781-663-5983

**Secondary:** Kay A. Taylor  
Tel: 317-418-1735  
Fax: 317-536-3064

**Trade Name:** GSP® Neonatal IRT kit (3306-001U)

**Common Name:** GSP Neonatal IRT kit  
**Regulation:** 21 CFR 862.1725

**Classification Name:** Trypsin Test System

**Product Code:** JNO

**Predicate Device:** AutoDELFIA® Neonatal IRT kit,  
510(k) Number (K0003668)

**Device Description:** The GSP Neonatal IRT assay is a solid phase, two-site fluoroimmuno-metric assay based on the direct sandwich technique in which two monoclonal antibodies (derived from mice) are directed against two separate antigenic determinants on the IRT molecule. Calibrators, controls or test specimens containing IRT are reacted simultaneously with immobilized monoclonal antibodies directed against a specific antigenic site on the IRT molecule and europium-labeled monoclonal antibodies (directed against a different antigenic site) in assay buffer. The assay buffer elutes IRT from dried blood on filter paper disks. The complete assay requires only one incubation step.

DELFLA Inducer dissociates europium ions from the labeled antibody into solution where they form highly fluorescent chelates with components of the DELFLA Inducer. The fluorescence in each well is then measured. The fluorescence of each sample is proportional to the concentration of IRT in the sample.

**Intended Use:**

The GSP Neonatal IRT kit is intended for the quantitative determination of IRT in blood specimens dried on filter paper as an aid in screening newborns for cystic fibrosis using the GSP® instrument.

**Substantial Equivalence:**

The GSP Neonatal IRT kit is substantially equivalent to our currently marketed AutoDELFLA IRT kit (K0003668). There are the following similarities and differences between the two kits:

Table 1. Characteristics of the two kits.

Characteristic	GSP Neonatal IRT kit  (New Device)	AutoDELFIA Neonatal IRT kit  (Predicate Device)
<b>Similarities</b>		
Intended User	Adequately trained laboratory personnel in laboratories performing newborn screening	Same
Intended Use / Indications for Use	The GSP Neonatal IRT kit is intended for the quantitative determination of IRT in blood specimens dried on filter paper as an aid in screening newborns for cystic fibrosis using the GSP® instrument.	The AutoDELFIA Neonatal IRT kit is intended for the quantitative determination of human IRT in blood specimens dried on filter paper as an aid in screening newborns for cystic fibrosis using the 1235 AutoDELFIA automatic immunoassay system.
Chemical Principle	<p>The GSP Neonatal IRT assay is a solid phase, two-site fluorimetric assay based on the direct sandwich technique in which two monoclonal antibodies (derived from mice) are directed against two separate antigenic determinants on the IRT molecule. Calibrators, controls, or test specimens containing IRT are reacted simultaneously with immobilized monoclonal antibodies directed against a specific antigenic site on the IRT molecule and europium-labeled monoclonal antibodies (directed against a different antigenic site) in assay buffer. The assay buffer elutes IRT from dried blood on filter paper disks. The complete assay requires only one incubation step.</p> <p>DELFIA Inducer dissociates europium ions from the labeled antibody into solution where they form highly</p>	Same

Characteristic	GSP Neonatal IRT kit  (New Device)	AutoDELFIA Neonatal IRT kit  (Predicate Device)
	fluorescent chelates with components of the DELFIA Inducer. The fluorescence in each well is then measured. The fluorescence of each sample is proportional to the concentration of IRT in the sample.	
Detection principle	Time-resolved fluorescence	Same
Specimen	Dried blood on filter paper disks with a diameter of approximately 3.2 mm (1/8 inch)	Same
Antibodies	Two different mouse monoclonal antibodies	Same
Calibrator and Control Matrix	Human blood derivative with a hematocrit of 50-55% and spotted onto filter paper cassettes (Whatman, no. 903)  (Washed RBCs in buffer containing BSA and protease inhibitors)	Same  (Washed RBCs in buffer containing saccharose)
Calibration	Calibrated using gravimetric methods	Same
Controls	3 levels (approx. values 30, 70 and 110 ng/mL blood)	Same
Assay buffer	IRT Assay Buffer, ready for use  Containing blockers and BSA item 1	Same  Containing BSA item 2
Calibrators	6 levels  (approx. values 0, 25, 50, 100, 250, 500 ng/mL blood)	Same
Expected Values	The measurement of IRT from dried blood spots is used as a means of identifying a population of newborns who are at increased risk of having CF and should be selected for 2nd tier testing. The identification is based on the use of a fixed cut-off value or population percentile. The IRT cut-off levels must be	Same

Characteristic	GSP Neonatal IRT kit (New Device)	AutoDELFIA Neonatal IRT kit (Predicate Device)
	determined by each newborn screening laboratory to meet the desired sensitivity and specificity of the screen and should be evaluated periodically.	
Coated Plates	Anti-IRT Microtitration Strips, 8 X 12 wells coated with antibodies directed against a specific site on the IRT molecule (mouse monoclonal)  Microtitration plate raw material item 1.	Same  Microtitration plate raw material item 2

Characteristic	GSP Neonatal IRT kit (New Device)	AutoDELFIA Neonatal IRT kit (Predicate Device)	
Differences			
Instrument	GSP Instrument	1235 AutoDELFIA Instrument	
Dissociation solution	DELFIA Inducer	Enhancement Solution	
Antibody Cross-Reactions in the Assay	$\alpha$ 2-macroglobulin	0.000%	$\alpha$ 2-macroglobulin < 4 ng/ml blood
	$\alpha$ 1-antitrypsin	0.000%	$\alpha$ 1-antitrypsin < 4 ng/ml blood
	Phospholipase A2	0.014%	Phospholipase A2 < 4 ng/ml blood
	Chymotrypsin	0.959%	Chymotrypsin < 4 ng/ml blood
	Human IgG	0.000%	Human IgG < 4 ng/ml blood
	Pepsinogen	-0.056%	(Uro)Pepsinogen < 4 ng/ml blood
	Complement Factor I	0.000%	Complement Factor I NA
Measuring Range	9 to 500 ng/mL blood	4(as defined by LoB) to 500 (as defined by upper calibrator) ng/mL blood	
	Linearity Range: 9 to500 ng/mL blood	Linearity Range: No claims for linearity in labeling.	
Tracer	Anti-IRT-Eu tracer stock solution, approximate concentration of ~40 $\mu$ g/mL mouse monoclonal,	Anti-IRT-Eu tracer stock solution, approximate concentration of ~50 $\mu$ g/mL mouse monoclonal,	

	ready for use.  Tracer antibody labeling with europium-chelate 1	ready for use.  Tracer antibody labeling with europium-chelate 2  Contains mouse IgG as blocker.
Analytical Sensitivity / Limit of Blank,  Limit of Detection  Limit of Quantitation	Limit of Blank 0.76 ng/mL blood  Limit of Detection 2.2 ng/mL blood  Limit of Quantitation 2.2 ng/mL blood	Limit of Blank < 4 ng/mL blood
Precision (Total Variation using a full calibration curve on each plate)	10.9 ng/mL blood CV% 7.3 22.2 ng/mL blood CV% 7.2 28.5 ng/mL blood CV% 7.0 40.0 ng/mL blood CV% 8.2 50.2 ng/mL blood CV% 8.0 61.6 ng/mL blood CV% 7.8 93.5 ng/mL blood CV% 7.2 302.3 ng/mL blood CV% 7.4 449 ng/mL blood CV% 7.5	42.6 ng/mL blood CV% 9.3 98.8 ng/mL blood CV% 10.0 266 ng/mL blood CV% 9.6



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Wallac Oy, Subsidiary of Perkin Elmer Inc.  
c/o Ms. Susan K. Hamann  
Regulatory Affairs Manager  
940 Winter Street  
Waltham, MA 02451

Food & Drug Administration  
10903 New Hampshire Avenue  
Building 66  
Silver Spring, MD 20993

DEC 16 2010

Re: k102419  
Trade Name: GSP Neonatal IRT kit (3306-001U)  
Regulation Number: 21 CFR §862.1725  
Regulation Name: Trypsin Test System  
Regulatory Class: Class I exempt, exceeds the limitation to exemption in 862.9(c)(2)  
Product Codes: JNO  
Dated: November 24, 2010  
Received: November 26, 2010

Dear Ms. Hamann:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (301) 796-5760. For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or ( 301 ) 796-5680 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in dark ink, appearing to read 'CCH', with a long horizontal flourish extending to the right.

Courtney C. Harper, Ph.D.  
Director  
Division of Chemistry and Toxicology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure



## Indications for Use Form

510(k) Number (if known):    K102419   

Device Name: GSP® Neonatal IRT Kit (3306-001U)

Indications for Use:

The GSP Neonatal IRT kit is intended for the quantitative determination of human immunoreactive trypsin(ogen) in blood specimens dried on filter paper as an aid in screening newborns for Cystic Fibrosis using the GSP® instrument.

Prescription Use    X     
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use             
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)



Division Sign-Off  
Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k)    K102419