

## 510(k) Summary

**Manufacturer:** Königsee Implantate GmbH  
OT Aschau, Am Sand 4,  
07426 Allendorf - GERMANY

**Dated Prepared:** January 10, 2011

**Device Trade Name:** Spinal Fusion Carrier (SFC™) VBR

**Contact:** Justin Eggleton  
Musculoskeletal Clinical Regulatory Advisers, LLC  
1331 H Street NW, 12<sup>th</sup> Floor  
Washington, DC 20005  
Office: 202.552.5800  
Fax: 202.552.5798

**Classification:** 21 CFR §888.3060

**Class:** II

**Product Code:** MQP

### Indications For Use:

SFC™ has been designed to be used in the thoracolumbar spine (T1 to L5) to replace a collapsed, damaged, or unstable vertebral body due to tumor or trauma (i.e. fracture). SFC™ is intended to be used with supplemental spinal fixation systems that are cleared by the FDA for use in thoracic and/or lumbar spine such as posterior pedicle screw and rod system, anterior plate systems, and anterior screw and rod systems. The SFC™ may be used with autograft or allograft.

### Device Description:

The SFC™ acts as an expandable spacer to maintain proper vertebral body spacing and angulation following vertebrectomy. The SFC™ VBR is manufactured from Ti6Al4V.

### Predicate Device(s):

The SFC™ was shown to be substantially equivalent to previously cleared devices (Innovotec SEC VBR, K091743; Aesculap Hydrolift VBR, K083186; Synthes Synex II, K061891) and has the same indications for use, design, function, and materials used.

### Performance:

Testing performed on this device indicates that the SFC™ is substantially equivalent to predicate devices. ASTM F2077 and ASTM F2267 performance standards (static axial compression, static axial torsion, static compression shear, dynamic axial compression,

dynamic axial torsion, subsidence, and expulsion) were adhered to and all applicable requirements were met.

**Conclusions:**

The subject and predicate devices share the same indications for use, design, function, and materials of manufacture. The non-clinical test results demonstrate that any minor differences do not impact device performance as compared to the predicates. The AccuLIF Cage was shown to be substantially equivalent to the cited predicate devices.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

OCT 12 2011

Konigsee Implantate GmbH  
% Musculoskeletal Clinical  
Regulatory Advisers, LLC  
Mr. Justin Eggleton  
1331 H Street NW, 12<sup>th</sup> Floor  
Washington, District of Columbia 20005

Re: K110153

Trade/Device Name: Konigsee Implantate Spinal Fusion Carrier (SFC<sup>TM</sup>) VBR  
Regulation Number: 21 CFR 888.3060  
Regulation Name: Spinal intervertebral fixation orthosis  
Regulatory Class: Class II  
Product Code: MQP  
Dated: September 26, 2011  
Received: September 27, 2011

Dear Mr. Eggleton:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K110153

Device Name: Königsee Implantate Spinal Fusion Carrier (SFC™) VBR

SFC™ has been designed to be used in the thoracolumbar spine (T1 to L5) to replace a collapsed, damaged, or unstable vertebral body due to tumor or trauma (i.e. fracture). SFC™ is intended to be used with supplemental spinal fixation systems that are cleared by the FDA for use in thoracic and/or lumbar spine such as posterior pedicle screw and rod system, anterior plate systems, and anterior screw and rod systems. The SFC™ may be used with autograft or allograft.


Prescription Use ✓  
(Part 29 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(29 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K110153