

JUL 29 2011

**510(k) SUMMARY**

Care back pain relief stimulator, K ( )

**Date of Submission:** 07/27/2011**510(k) Submitter's Name:** Shenzhen Dongdixin Technology Co., Ltd.

No.3bldg. xiliyangguang industrial Estae

**Address:** xilixiaobaimang Nanshan district Shenzhen, China  
518108**Telephone:** +86(755) 27652471**Fax:** +86(755) 27652674**510(k) Correspondent  
contact:** Kang Jian Ping**E-mail:** microcong@gmail.com

### **1. Proposed Device:**

Trade Name: Care back pain relief stimulator  
Model: TENS 7000,LT3074  
Classification Name: Stimulator, Nerve, Transcutaneous.  
Over-The-Counter  
Regulation Number: 21CFR 882.5890  
Product Code: NUH  
Device Class: II

### **2. Predicate Device:**

Predicate Device: Low Back Pain Relief System  
510(k) Number: K060222  
Manufacturer: Gemore Technology Co.,Ltd

### **3. Description of Proposed Device:**

Care back pain relief stimulator, which includes models TENS 7000 and LT3074 are non-invasive devices which are tented for over the counter use in temporary relief of pain associated with sore and aching muscles in the lower back due to strain form exercise or normal household and working activities.

The devices contain the following main parts: TENS stimulation unit, Support Belt, electrode, and snap cable.

The electrode is permanently sewn up in the support belt. Each belt is flexible, and is available in range of sizes to ensure good patient contact. A male snap connector is placed within the electrode and is connected via the female snap connector to a short lead wire. The lead wire has a female pin connection at the distal end which accepts the lead wire connection from the stimulator.

The device can be worn on the low back part of user so as to place the electrode on the treatment location of low back.

### **4. Proposed Device Intended Use Statement:**

Care back pain relief stimulator is intended for temporary relief of pain associated with sore and aching muscles in the lower back due to strain form exercise or

normal household and working activities.

#### 5. Technological Characteristics:

Both the stimulator and the Predicate device stimulator have the same intended use and fundamental technology. A side-by-side comparison of the technological characteristics of the care back pain relief and the cited predicate devices is included in the 510(k) submission.

Parameter	Care back pain relief stimulator Model: TENS7000,LT3074	Low Back Pain Relief System, model GM320PP
Indications for use	Care back pain relief stimulator is intended for temporary relief of pain associated with sore and aching muscles in the lower back due to strain from exercise or normal household and working activities.	The model GM320PP is intended for temporary relief of pain associated with sore and aching muscles in the lower back due to strain from exercise or normal household and work activities.
Power Source	9V Battery for TENS7000 4.5V Battery for LT3074	9V Battery
Max Output Voltage (V) ±20%	40.0V@500Ω for TENS 7000	40.0V@500Ω
	45.0V@500Ω for LT3074	
	100.0V@2kΩ for TENS 7000	120.0V@2kΩ
	50.0V@2kΩ for LT3074	
	110.0V@10kΩ for TENS 7000	150.0V@10kΩ
	60.0V@10kΩ for LT3074	
Max Output Current (mA) ±20%	80.0mA @500Ω for TENS 7000	80.0 mA @500Ω
	90.0mA @500Ω for LT3074	
	50.0mA @2kΩ for TENS 7000	60.0 mA @2kΩ
	25.0mA @2KΩ for LT3074	
	11.0mA @10kΩ for TENS 7000	15.0 mA @10kΩ
	6.0mA @2kΩ for LT3074	
Pulse Width Range	50-300us for TENS 7000	150-260 us
	50-250us for LT3074	
Frequency	2-150Hz for TENS 7000	2-120Hz
	2-110Hz for LT3074	
Timer Range (minutes)	5 - 60 minutes or Continuous for TENS7000	15min, 30min, 60min and continuous
	30 minutes for LT3074	
Maximum Average Current (average absolute value), mA <i>(refer to the remark 1)</i>	3.6mA,500Ω, for TENS 7000	2.5mA,500Ω
	2.5mA,500Ω, for LT3074	
RMS current <i>(refer to the remark 1)</i>	17.0 mA.r.m.s for TENS 7000	20.0 mA.r.m.s
	14.2 mA.r.m.s for LT3074	

Maximum Current Density <i>(refer to the remark 1)</i>	0.56mAmps/cm <sup>2</sup> ,500Ω for TENS 7000	0.56mAmps/cm <sup>2</sup> ,500Ω
	0.14mAmps/cm <sup>2</sup> ,500Ω for LT3074	
Maximum Power Density <i>(refer to the remark 1)</i>	4.8mWrms/cm <sup>2</sup> ,500Ω for TENS 7000	3.3mWrms/cm <sup>2</sup> ,500Ω
	1.0mWrms/cm <sup>2</sup> ,500Ω for LT3074	
Compliance with Voluntary Standards?	IEC60601-1, IEC60601-1-2 IEC60601-2-10	IEC60601-1, IEC60601-1-2 IEC60601-2-10
Compliance with 21 CFR 898?	Yes	Yes
Weight (grams.)	150 grams with battery for TENS7000	140 grams with battery
	80 grams with battery for LT3074	
Dimensions (mm.) H x W x T	101*61*24.5cm for TENS7000	108*61.5*25
	85*72*30cm for LT3074	
Housing Materials & Construction	Enclosure: ABS,94V-1,80°C,UL Approved	Enclosure: ABS,94V-1,80°C,UL Approved
Materials of electrodes:	Rubber or fabric	Gel

## 6. Non-Clinical Tests Performed:

Care back pain relief stimulator did not conduct, nor rely upon, clinical tests to determine substantial equivalence. Care back pain relief stimulator has received test reports certifying that the Care back pain relief stimulator was tested and found to be in conformity with voluntary standards includes IEC60601, IEC60601-1-2, IEC60601-2-10, ISO 10993 and ISO14971. The software verification has been carried out according to the FDA Guidance of the content of Premarket submissions for the software contained in Medical Devices. This test reports have provided as part of this premarket notification.

## 7. Conclusions:

The Care back pain relief Stimulator, which includes models TENS 7000 and LT3074, has the same intended use and technological characteristics as the cited predicate device. Moreover, bench testing, safety report and Risk Analysis Report documentation supplied in this submission demonstrates that the difference in the submitted models could maintain the same safety and effectiveness as that of the cited predicate device. In the other words, those engineering difference do not affect the intended use or alter the fundamental scientific technology of the device. Thus, the care back pain relief stimulator is substantially equivalent to the predicate device.

Remark 1

**TENS 7000 sample calculations:**

Maximum Average Current= $I_{max} \times \text{pulse rate} \times \text{pulse width} = 80\text{mA} \times 150\text{Hz} \times 300\text{us} \times 10^{-6} = 3.6\text{mA}$

$$I_{r.m.s} = \sqrt{\frac{5}{5} \times 150\text{Hz} \times \frac{300\text{us}}{1000000}} \times 80\text{mA} = 17.0\text{mA}_{r.m.s}$$

$$\text{Maximum Current Density} = \frac{I_{r.m.s}}{S} = \frac{17.0\text{mA}_{r.m.s}}{30\text{cm}^2} = 0.56\text{mA}_{r.m.s} / \text{cm}^2$$

$$V_{r.m.s} = \sqrt{\frac{5}{5} \times 150\text{Hz} \times \frac{300\text{us}}{1000000}} \times 40\text{V} = 8.5\text{V}_{r.m.s}$$

$$\text{Maximum Voltage Density} = \frac{V_{r.m.s}}{S} = \frac{8.5\text{V}_{r.m.s}}{30\text{cm}^2} = 0.28\text{V}_{r.m.s} / \text{cm}^2$$

$$\text{Maximum Power Density} = \frac{V_{r.m.s} \times I_{r.m.s}}{S} = \frac{8.5\text{V}_{r.m.s} \times 17\text{mA}_{r.m.s}}{30\text{cm}^2} = 4.8\text{mW}_{r.m.s} / \text{cm}^2$$

**LT3074 sample calculations:**

Maximum Average Current= $I_{max} \times \text{pulse rate} \times \text{pulse width} = 90\text{mA} \times 100\text{Hz} \times 250\text{us} \times 10^{-6} = 2.5\text{mA}$

$$I_{r.m.s} = \sqrt{\frac{5}{5} \times 100\text{Hz} \times \frac{250\text{us}}{1000000}} \times 90\text{mA} = 14.2\text{mA}_{r.m.s}$$

$$\text{Maximum Current Density} = \frac{I_{r.m.s}}{S} = \frac{14.2\text{mA}_{r.m.s}}{100\text{cm}^2} = 0.14\text{mA}_{r.m.s} / \text{cm}^2$$

$$V_{r.m.s} = \sqrt{\frac{5}{5} \times 100\text{Hz} \times \frac{250\text{us}}{1000000}} \times 45\text{V} = 7.1\text{V}_{r.m.s}$$

$$\text{Maximum Voltage Density} = \frac{V_{r.m.s}}{S} = \frac{7.1\text{V}_{r.m.s}}{100\text{cm}^2} = 0.07\text{V}_{r.m.s} / \text{cm}^2$$

$$\text{Maximum Power Density} = \frac{V_{r.m.s} \times I_{r.m.s}}{S} = \frac{7.1\text{V}_{r.m.s} \times 14.2\text{mA}_{r.m.s}}{100\text{cm}^2} = 1.0\text{mW}_{r.m.s} / \text{cm}^2$$



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Handelhaus Dittman GMBH  
Shenzen Dongdixin Technology Co. LTD.  
c/o Mr. Kang Jiang Ping  
R.A. Director  
No. 3 Building Xiliyangguang  
Industrial Estate Xillixaobaim  
Shenzen  
China 518108

JUL 29 2011

Re: K110390

Trade/Device Name: ETG 255 Combo Stimulator, TEN 260 Stimulator  
Regulation Number: 21 CFR 882.5890  
Regulation Name: Transcutaneous electrical nerve stimulator for pain relief  
Regulatory Class: Class II  
Product Codes: GZJ, IPF  
Dated: November 30, 2011  
Received: February 11, 2011

Dear Mr. Jiang Ping:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in cursive script, appearing to read "Malvina B. Eydelman".

Malvina B. Eydelman, M.D.  
Director  
Division of Ophthalmic, Neurological,  
Ear, Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number ( ):

Device Name:

**ETG 255 Combo Stimulator, TEN 260 Combo Stimulator**

### Indications for Use:

ETG 255 Combo Stimulator,

#### For TENS mode

Symptomatic relief of chronic intractable pain, relief of acute post-surgical and posttraumatic pain.

#### For EMS mode

1. Relaxation of muscle spasm.
2. Increase of local blood flow circulation
3. Prevention or retardation of disuse atrophy
4. Muscle re-education
5. Maintaining or increasing range of motion.
6. Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis

#### TEN 260 Stimulator

Symptomatic relief of chronic intractable pain, relief of acute post-surgical and posttraumatic pain.

Prescription Use ✓  
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Ophthalmic, Neurological and Ear,  
Nose and Throat Devices

510(k) Number K110390