

### 510(k) SUMMARY

#### Houston Medical Robotics, Inc. Euclid™ Tier 1 Mini Access System

This 510(k) Summary is submitted in accordance with 21 CFR Part 807, section 807.92.

**Sponsor's Name and Address:** Houston Medical Robotics, Inc.  
17225 El Camino Real, Suite 350  
Houston, Texas 77058

**Contact Person:** Darla J. Elkin  
Elkin RC, LLC  
42 North Chantson Circle  
The Woodlands, TX 77382  
Telephone: (281) 450-8163  
Facsimile: (281) 298-7277

**Date Summary Prepared:** May 20, 2011

**Device Trade Name:** Euclid™ Tier 1 Mini Access System

**Product Code:** ITX

**Regulation Number:** 21CFR 892.1570

**Classification:** Class II

**Common Name:** Diagnostic ultrasonic transducer

**Predicate Devices:** The WAND™ (K081697)  
Ultra-Pro II Needle Guidance System (K093713)

#### Device Description:

The Euclid™ Tier 1 Mini Access System is designed to provide the medical practitioner with the capability to accurately and reliably insert a guidewire into a vessel. The device consists of (1) the Euclid™ Tier 1 Mini Base Kit containing the sterile disposable device assembly consisting of a needle (18/21 Gauge) and guidewire (0.18" diameter x 40 cm length), a Euclid™ Ultrasound Imaging Screen Overlay, and IFU; and (2) the Euclid™ Tier 1 Mini Transducer Adapter Kit which contains a Euclid™ Transducer Adapter and IFU.

#### Intended Use:

The Euclid™ Tier 1 Mini Access System is used to facilitate the placing of a needle and guidewire into a targeted anatomical location.

## Comparison of the Technological Characteristics of the New Device and Predicate Devices:

| Features                               | Euclid™ Tier 1 Mini Access System   | The WAND Microaccess Safety Introducer  | Ultra-Pro II Needle Guidance System  |
|--|---|---|--|
| <b>K#</b>                              | K111426   | K093022   | K093713  |
| <b>Disposable Components</b>           | Sterile Disposable: <ul style="list-style-type: none"> <li>• Needle</li> <li>• Guidewire</li> </ul> | Sterile Disposable: <ul style="list-style-type: none"> <li>• Needle</li> <li>• Guidewire</li> <li>• Sheath</li> </ul> | Sterile Disposable: <ul style="list-style-type: none"> <li>• Disposable Needle Guide</li> <li>• Cover</li> </ul> |
| <b>Reusable Component</b>              | None  | None  | Bracket to hold Ultrasound Transducer  |
| <b>Guidewire</b>                       | Nitinol 0.018"  | Nitinol 0.018"  | N/A  |
| <b>Angle Needle Path</b>               | Yes (Adjustable & Fixed During Use)   | No  | Yes (Fixed)  |
| <b>Placement</b>                       | Manual  | Manual  | Manual   |
| <b>Ultrasound Guided Procedure</b>     | Yes<br>Couples with commercially available ultrasound   | Yes<br>Labeling suggests use of commercially available image guidance.  | Yes<br>Couples with commercially available ultrasound.   |
| <b>Vessel Depth Setting</b>            | Yes   | No  | Yes  |
| <b>Needle</b>                          | 18 /21 GA thin-walled stainless steel   | 21 GA thin-walled stainless steel   | N/A  |
| <b>Ultrasound Visibility</b>           | Yes   | Yes   | No   |
| <b>Visual Display of Vessel Access</b> | Yes, external Ultrasound Screen   | Yes, Flashback window   | Yes, external Ultrasound Screen  |
| <b>EtO Sterilized</b>                  | Yes   | Yes   | Yes  |

## Performance Testing

Results of bench studies conducted on the Euclid™ Tier 1 Mini Access System demonstrate the System to be as safe and as effective as the predicate device based on the biocompatibility of the materials used, sterilization validation, bench testing and verification and validation in animal models. The following studies were conducted on the Euclid™ Tier 1 Mini Access System and acceptance criteria were met:

- Functional Verification**  
 Benchtop testing (Guidewire insertion) to the functional extremes (5mm and 60mm);
- Shipping Integrity**  
 Simulated shipping conditions followed by device and packaging inspection, guidewire insertion in benchtop model to functional extremes (5mm and 60mm), pouch seal strength testing, and dye penetration testing;
- Accelerated Aging**  
 Aging to 6 months followed by functional testing (guidewire insertions to 5mm and 60mm), pouch seal strength testing, and dye penetration testing;

- **Animal Validation**  
Guidewire insertion in animal model by end user; and
- **Biocompatibility Testing**  
Materials tested per ISO 10993, including: Systemic Toxicity, Intracutaneous Toxicity, Implantation, Cytotoxicity, Hemolysis, Pyrogenicity, and Sensitization.

### **Conclusion**

The Euclid™ Tier 1 Mini Access System is substantially equivalent to the predicate devices. The indication for the devices is substantially equivalent. The technological design and functional characteristics of placement of a guidewire using ultrasound guidance with a sterile disposable device are all substantially equivalent to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

MAR - 1 2012

Houston Medical Robotics, Inc.  
c/o Ms. Darla J. Elkin  
Elkin RC, LLC  
42 North Chantsong, Circle  
The Woodlands, TX. 77382

Re: K111426  
Trade Name: Euclid™ Tier 1 Mini Access System  
Regulation Number: 21 CFR 892.1570  
Regulation Name: Diagnostic Ultrasonic Transducer  
Regulatory Class: II (two)  
Product Code: ITX  
Dated: February 24, 2012  
Received: February 27, 2012

Dear Ms. Elkin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Ms. Darla J. Elkin

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



*for* Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K111426

Device Name: Euclid™ Tier 1 Mini Access System


Indications for Use:

The Euclid™ Tier 1 Mini Access System is used to facilitate the placing of a needle and guidewire into a targeted anatomical location.

Prescription Use X AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

  
(Division Sign-Off) Concurrency of CDRH, Office of Device Evaluation (ODE)  
Division of Cardiovascular Devices

510(k) Number K111426