

## 510(K) SUMMARY OF SAFETY AND EFFECTIVENESS

**Trade Name:** MicroPlex Coil System – Vector SEP 20 2011

**Generic Name:** Neurovascular Embolization Device

**Classification:** Class II, 21 CFR 882.5950

**Submitted By:** MicroVention, Inc  
1311 Valencia Avenue  
Tustin, California 92780 U.S.A.

**Contact:** Laraine Pangelina

**Predicate Device:** MicroPlex Coil System – Cosmos 10 (K082461, K093919, K103758)  
MicroPlex Coil System – Cosmos 18 (K090891, K093358)

### Device Description:

The MCS Vector consists of an implantable coil made of bare platinum alloy. The Vector implantable coils has a 3D shape in various loop sizes and lengths. The coil is attached to a V-Trak delivery pusher. The proximal end is inserted into a hand held battery powered V-Grip Detachment Controller (sold separately). The implant segment detaches upon activation of the Detachment Controller.

The table below provides information about the physical properties of the MCS Vector with a comparison to the predicate devices.

Feature	MCS – Cosmos 10	MCS – Cosmos 18	MCS – Vector
Coil shape	3D	Same	Same
Coil implant diameter	2-12mm	6-24mm	3-15mm
Coil restrained length	2-45cm	17-68cm	6-60cm
Deliver pusher length	185cm	Same	Same
Main coil wire material	Platinum/Tungsten alloy	Same	Same
Coupler material	Platinum/Iridium	Same	Same
Adhesive material	Ultraviolet cure	Same	Same
Implant to pusher material	Polyolefin elastomer	Same	Same
Stretch resistant filar material	Polyolefin elastomer	Same	Same
MRI compatibility	Yes	Yes	Yes
Method of supply	Sterile, single use	Same	Same
Packaging configuration	Dispenser coil, pouch, shipping carton	Same	Same

**Indications for Use:**

Intended for the endovascular embolization of intracranial aneurysms and other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistula. The MCS is also intended for vascular occlusion of blood vessels within the neurovascular system to permanently obstruct blood flow to an aneurysm or other vascular malformation and for arterial and venous embolization in the peripheral vasculature.

**Bench Test Summary:**

Test	Result
Dimensional Measurement	Met same criteria as predicate
Simulated Use: Introduction, Tracking, Deployment, Frame movement, Microcatheter movement, Microcatheter manipulation, Compartmentalization, Detachment, Overall performance	Met same criteria as predicate
Spring Constant	Met same criteria as predicate
Weld Tensile Strength	Met same criteria as predicate
Detachment Zone Tensile Strength	Met same criteria as predicate

**Summary of Substantial  
Equivalence:**

The MCS Vector is substantially equivalent to the predicate devices with regard to intended use, patient population, device design, materials, processes, and operating principal.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - W066-0607  
Silver Spring, MD 20903-0002

Microvention, Inc.  
c/o Ms. Laraine Pangelina  
1311 Valencia Ave  
Tustin, CA 92780

SEP 20 2011

Re: K111451

Trade/Device Name: Microplex Coil System-Vector  
Regulation Number: 21 CFR 882.5950  
Regulation Name: Device, Neurovascular embolization  
Regulatory Class: Class II  
Product Code: HCG  
Dated: August 9, 2011  
Received: August 10, 2011

Dear Ms. Pangelina:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in dark ink, appearing to read 'Malvina B. Eydelman'.

Malvina B. Eydelman, M.D.

Director

Division of Ophthalmic, Neurological,  
and Ear, Nose and Throat Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known): K111451

Device Name: MicroPlex Coil System – Vector

Indications for Use: *Intended for the endovascular embolization of intracranial aneurysms and other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae. The MCS is also intended for vascular occlusion of blood vessels within the neurovascular system to permanently obstruct blood flow to an aneurysm or other vascular malformation and for arterial and venous embolizations in the peripheral vasculature.*

Prescription Use X  
(Per 21 CFR 801.109)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(Optional Format 1-2-96)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Jeffrey Toy  
(Division Sign-Off)  
Division of Ophthalmic, Neurological and Ear,  
Nose and Throat Devices

510(k) Number K111451