

**Section 5 – 510(k) Summary or 510(k) Statement**

K121150

Page 1 of 8

SEP 21 2012

**I. General Information**

Submitter:

Alma Lasers, Inc.

485 Half Day Road Suite 100

Buffalo Grove IL 60089

Contact Person:

Kathy Maynor - Regulatory

352-586-3113 (cell)

Summary Preparation Date: June 28, 2012

**II. Names**

Device Name(s): Accent Pixel RF Tips

Primary Classification Name(s): Electrosurgical cutting and coagulation device and accessories

**III. Predicate Devices**

- Alma Lasers Family of Accent RF Systems [Accent, Accent XL, Accent Elite] (K072699, K070004, K101147)
- Syneron Medical Matrix RF Applicator (K073572, K090025)
- EndyMed RF Applicator (K083461, K101510, K120513)

#### **IV. Product Description**

The Accent Pixel RF Tips include cleanable, sterilizable, multiple use radiofrequency (RF) energy delivery devices (accessory) intended for use with the Alma Lasers Family of Accent™ RF Systems [Accent, Accent XL, Accent Elite]. These Pixel RF Tips are also available as single use, disposable products.

The Accent Pixel RF Tips are comprised of the following main components:

- Pixel<sup>RF</sup> Stationary Tip
- Pixel<sup>RF</sup> Roller Tip

#### **V. Intended Use & Indications for Use**

##### **Intended Use**

The Alma Lasers family of Accent RF Systems (Accent, Accent XL, Accent Elite) with Unipolar and Bipolar Handpieces are intended for use in dermatologic and general surgical procedures.

##### **Indications for Use**

The Unipolar and Bipolar handpieces, when used as a combined treatment, are indicated for the non-invasive treatment of wrinkles and rhytids.

The Pixel RF Tips (stationary and rolling), when used with the unipolar handpiece, are indicated for dermatological procedures requiring ablation and resurfacing of the skin.

#### **VI. Rationale for Substantial Equivalence**

The Accent Pixel RF Tips shares the same or similar indications for use, device operation, overall technical and functional capabilities, and therefore is substantially equivalent to the predicate devices. Additionally, animal and human clinical data were submitted to substantiate product performance.

K 121150

Page 3 of 8

Characteristic	K0XXXX Alma Lasers, Ltd. Family of Accent RF Pixel Tips with Pixel <sup>RF</sup> Module/ Handpiece	K072699, K070004, K101147 Alma Lasers Family of Accent RF Systems [Accent, Accent XL, Accent Elite] (w/ UniPolar & BiPolar Handpieces)	K073572, K090025 Syneron Medical Matrix RF Applicator	K083461, K101510, K120513 EndyMed RF applicator
Product Code Regulation	General & Plastic Surgery • GEL, OUH - 21 CFR 878.4400	General & Plastic Surgery • NUV, GEL, ISA 21 CFR 878.4810	General & Plastic Surgery • GEL, OUH - 21 CFR 878.4400	General & Plastic Surgery GEL, OUH - 21 CFR 878.4400
Intended Use	Intended for use in dermatologic and general surgical procedures	Intended for use in dermatologic and general surgical procedures	Intended for use in dermatologic and general surgical procedures	Intended for use in dermatologic and general surgical procedures

Characteristic	K0XXXX Alma Lasers, Ltd. Family of Accent RF Pixel Tips with Pixel <sup>RF</sup> Module/ Handpiece	K072699, K070004, K101147 Alma Lasers Family of Accent RF Systems [Accent, Accent XL, Accent Elite] (w/ UniPolar & BiPolar Handpieces)	K073572, K090025 Syneron Medical Matrix RF Applicator	K083461, K101510, K120513 EndyMed RF applicator
Indications for Use	<p>The Alma Lasers family of Accent RF Systems (Accent, Accent XL, Accent Elite) with Unipolar and Bipolar Handpieces are intended for use in dermatologic and general surgical procedures.</p> <p>The Unipolar and Bipolar handpieces, when used as a combined treatment, are indicated for the non-invasive treatment of wrinkles and rhytids.</p> <p>The Pixel RF Tips (stationary and rolling), when used with the unipolar handpiece, are indicated for dermatological procedures requiring ablation and resurfacing of the skin.</p> <p><b>Supported by animal histology data – section 19 – and human clinical data – section 20.</b></p>	<p>Indicated for: The non-invasive treatment of <b>wrinkles and rhytids</b> using a combined treatment with Unipolar and Bipolar handpieces</p>	<p>The Matrix RF Applicator is intended for dermatological procedures requiring ablation and resurfacing of the skin, and for the treatment of facial wrinkles.</p>	<p><i>The Glow by EndyMed</i> is a noninvasive device intended for use in Dermatologic and General Surgical procedures.</p> <p>The TC applicator is indicated for mild to moderate facial wrinkles and rhytides.</p> <p>The FSR applicator is indicated for Dermatological procedures requiring ablation and resurfacing of the skin.</p>

Characteristic	K0XXXX Alma Lasers, Ltd. Family of Accent RF Pixel Tips with Pixel <sup>RF</sup> Module/ Handpiece	K072699, K070004, K101147 Alma Lasers Family of Accent RF Systems [Accent, Accent XL, Accent Elite] (w/ UniPolar & BiPolar Handpieces)	K073572, K090025 Syneron Medical Matrix RF Applicator	K083461, K101510, K120513 EndyMed RF applicator
RF Power (or Energy)	<ul style="list-style-type: none"> <li>UniPolar               <ul style="list-style-type: none"> <li>Recommended Range - Stationary Tips 45-60 W</li> <li>Recommended Range - Roller Tips 35-55 W</li> <li>Maximum for Accent System Up to 200 W (2.08 W/mm<sup>2</sup>) (Nominal 35-80 W - Pixel<sup>RF</sup> Tips - range of 0.15 to 2.08 W/mm<sup>2</sup>)</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>UniPolar &amp; UniLarge Up to 200 W (1.51 W/mm<sup>2</sup>) (Nominal 120-140 W)</li> <li>BiPolar &amp; BiPolarS Up to 420 W (1.51 W/mm<sup>2</sup>) (Nominal 250-290 W)</li> </ul>	<ul style="list-style-type: none"> <li>Bipolar - Up to 25 Joules per cm<sup>2</sup> (Tunable SelectPulse™ fractional technology - 3 programs which emulate the tissue effects of common skin resurfacing fractional laser technologies)</li> </ul>	<ul style="list-style-type: none"> <li>Bipolar up to 65 watts Up to 62mj/pin</li> </ul>

K121150

Page 6 of 8

Characteristic	K0XXXX Alma Lasers, Ltd. Family of Accent RF Pixel Tips with Pixel <sup>RF</sup> Module/ Handpiece	K072699, K070004, K101147 Alma Lasers Family of Accent RF Systems [Accent, Accent XL, Accent Elite] (w/ UniPolar & BiPolar Handpieces)	K073572, K090025 Syneron Medical Matrix RF Applicator	K083461, K101510, K120513 EndyMed RF applicator
	Tip Matrix	Tip Matrix	Tip Matrix	Tip Matrix
Treatment Area	Stationary Pixel <sup>RF</sup> Tips	Not applicable	64 Matrix Points – per pulse	112 matrix spots
	• ~110 Pixel points (1.0 pixels/mm <sup>2</sup> )	(K072699 & K070004)	0.4 points/mm <sup>2</sup>	(1.3 matrix spots/mm <sup>2</sup> )
Treatment Area	Roller Pixel <sup>RF</sup> Tips			
	• ~132 Pixel points per cm (1.3 pixels/mm <sup>2</sup> )			
Treatment Area	Stationary Pixel <sup>RF</sup> Tips (Unipolar)	• UniPolar & UniLarge	12 x 12 mm – disposable tip	10 x 15 mm
	12 mm dia.	13 to 19 mm dia.	‘spot’ size (unipolar)	
Treatment Area	Roller Pixel <sup>RF</sup> Tips (Unipolar)	• BiPolarS & BiPolar		
	10 x 22.7 mm circumference ‘roller’	8.95 to 13 mm dia.		
Handpiece Dimensions	Pixel <sup>RF</sup> Handpiece - 160 x 155 mm	• UniPolar - 169 x 205 mm	5.9 x 6.3 inches (15 x 16 cm)	Unknown
		• BiPolar - 167 x 203 mm		

K121150

Page 7 of 8

Characteristic	K0XXXX Alma Lasers, Ltd. Family of Accent RF Pixel Tips with Pixel <sup>RF</sup> Module/ Handpiece	K072699, K070004, K101147 Alma Lasers Family of Accent RF Systems [Accent, Accent XL, Accent Elite] (w/ UniPolar & BiPolar Handpieces)	K073572, K090025 Syneron Medical Matrix RF Applicator	K083461, K101510, K120513 EndyMed RF applicator
Patient Contact Materials	<ul style="list-style-type: none"> <li>Stationary Pixel<sup>RF</sup> Tips – Aluminum</li> <li>Roller Pixel<sup>RF</sup> Tips – Aluminum, Sapphire</li> </ul>	Handpieces: <ul style="list-style-type: none"> <li>Anodized aluminum and Delrin plastic</li> </ul>	Unknown	Unknown
How provided	<b>Disposable and Reusable Stationary &amp; Roller Pixel<sup>RF</sup> Tips:</b> <ul style="list-style-type: none"> <li>Reusable, cleanable, sterilizable Tips and disposable tips also available</li> </ul>	Not applicable	Disposable tip – up to 125 pulses	Disposable tip – up to 300 pulses
Electrical Reqs	110–120 V; 50–60 Hz; 5A	110–120 V; 50–60 Hz; 5A	90–230 VAC, 50–60 Hz, 7.5A	100–230 VAC; 1–3 A; 50–60 Hz
System Dimensions (“)	26 x 17 x 16 AccentXL; 21 x 17 x 38 Accent; 26 x 18 x 16 Accent Beauty	26 x 17 x 16 AccentXL 21 x 17 x 38 Accent	20.87 x 15.35 x 11.42 (Tabletop system)	115 x 30 x 19 cm; 45 x 12 x 7 in.
System Weight	110(XL);	55 lbs	15.4 lbs (7 kg)	33 kg; 72.5 lb.

**VII. Safety and Effectiveness Information**

The review of the indications for use and technical characteristics provided demonstrates that the Accent Pixel RF Tips are substantially equivalent to the predicate devices. Additionally, animal and human clinical data were submitted to substantiate performance.

**VIII. Conclusion**

The Accent Pixel RF Tips were found to be substantially equivalent to the predicate devices.

The Accent Pixel RF Tips share the identical indications for use, similar design features, and functional features with, and thus are substantially equivalent to, the predicate devices.





DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Alma Lasers Incorporated  
% Ms. Kathy Maynor  
Regulatory Consultant  
485 Half Day Road, Suite W No. 100  
Buffalo Grove, Illinois 60089

SEP 21 2012

Re: K121150

Trade/Device Name: Pixel RF Tips  
Regulation Number: 21 CFR 878.4400  
Regulation Name: Electrosurgical cutting and coagulation device and accessories  
Regulatory Class: Class II  
Product Code: GEI  
Dated: August 31, 2012  
Received: September 11, 2012

Dear Ms. Maynor:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "For Paul DeLeon", is written over the typed name and title of Mark N. Melkerson.

Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use Statement**

510(k) Number (if known): ~~K12~~

K121150

Device Name: Pixel RF Tips

Indications for Use:

**Intended Use**

The Alma Lasers Family of Accent RF Systems [Accent, Accent XL, Accent Elite] is intended for use in dermatologic and general surgical procedures.

**Indications for Use**

The Unipolar and Bipolar handpieces, when used as a combined treatment, are indicated for the non-invasive treatment of wrinkles and rhytids.

The Pixel RF Tips (stationary and rolling), when used with the unipolar handpiece, are indicated for dermatological procedures requiring ablation and resurfacing of the skin.

Neil R. Gade for rxm  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K121150

Prescription Use ☒

Over-The-Counter Use

AND/OR

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page

1 of 1