Hutchinson Santé G-DermTM Surgical Glove

August 28, 2012 Traditional 510(k)

5. 510(k) SUMMARY

[per 21 CFR 807.92]

SEP 1 3 2012

5.1 **Submitter Information**

Name:

Hutchinson Santé S.N.C.

Address:

Rue Marret et Paturel

F-60140 Liancourt

FRANCE

Phone:

+33 (0) 3 44 73 87 00

Fax:

+33 (0) 3 44 73 87 09

Contact Name: Dr Raffi Krikorian, Ph.D.

Date of summary: August 28, 2012

5.2 **Device Identification**

Trade name:

G-DermTM Powder Free Surgical Glove made of styrene-

ethylene/butylene-styrene (SEBS) and styrene-

ethylene/propylene-styrene-ethylene/propylene (SEPSEP)

synthetic copolymers

Common name:

Surgical glove

Classification name: Surgeon's glove

Product Code:

KGO

Regulation number:

21 CFR 878.4460

Device Class:

Class I (general controls)

5.3 **Identification of Predicate Devices**

- Elastyfree by ECI Medical Technologies Inc (K020918)
- Safeskin Tactylon PF powder-free surgical gloves (K994081)
- Tactylon by Tactyl Technologies (K955419)

5.4 **Device Description**

The G-DermTM Powder Free Surgical Glove is made of a blend of styreneethylene/butylene-styrene (SEBS) and styrene-ethylene/propylene-styreneethylene/propylene (SEPSEP) synthetic copolymers. It is coated with a polyurethane based coating on the inner side to facilitate donning.

5.5 Indications for Use

The G-DermTM Surgical Glove is a disposable sterile powder-free medical device made of styrene-ethylene/butylene-styrene (SEBS) and styrene-ethylene/propylene-styrene-ethylene/propylene (SEPSEP) synthetic copolymers that is intended to be worn by operating room personnel to protect a surgical wound from contamination.

5.6 Comparison to Predicate Devices

	G-Derm ^{IM}	Elastyfree	Safeskin	Tactylon
- w A	40.75	9 50 220 120	Tactylon	•
Regulation No.	21 CFR	21 CFR	21 CFR	21 CFR
	878.4460	878.4460	878.4460	878.4460
Device Class	Class I	Class I	Class I	Class !
Technology	Solvent borne material	Solvent borne material	Solvent borne material	Solvent borne material
Patient	Styrenic	Styrenic	Styrenic	Styrenic
contacting	thermoplastic	thermoplastic	thermoplastic	thermoplastic
physical barrier material	elastomers	elastomers	elastomer	elastomer
User contacting	Styrenic	Styrenic	Styrenic	Styrenic
physical barrier	thermoplastic	thermoplastic	thermoplastic	thermoplastic
material	elastomers with	elastomers with	elastomers with	elastomer
	a polyurethane	a polyurethane	a polymer	
	coating	coating	coating	
Standards met	ASTM D3577	ASTM D3577	ASTM D3577	ASTM D3577
	ASTM D5151	ASTM D5151	ASTM D5151	ASTM D5151
	ASTM D6124	ASTM D6124	ASTM D6124	ASTM D6124
Sterilization	Sterile (SAL 10	Sterile (SAL	Sterile. Gamma	Sterile. (method
	⁶). Gamma	10 ⁻⁶). Gamma	irradiation	unknown)
	irradiation	irradiation		

5.7 Performance Data

The G-DermTM Surgical Glove possesses the following technological characteristics:

Characteristics	Standard	
Thickness: 0.21 ±0.02 mm	Meets ASTM D3577	
Length: 280 mm min.	Meets ASTM D3577	
Physical Properties	Meets ASTM D3577,	

Type 2	
Meets ASTM D6124	
Meets ASTM D5151	
ISO 11737	
ISO 10993-1	
ISO 10993-5	
ISO 10993-11	
ISO 10993-10	
ISO 10993-10	

5.8 Clinical Data

Clinical data is not needed for medical glove 510(k) submissions.

It can be concluded that the G-DermTM Surgical Glove will perform according to the performance standards referenced above, FDA requirements and the labeling claims for this product. Consequently, this device is substantially equivalent to currently marketed devices.



Food and Drug Administration 10903 New Hampshire Avenue Document Control Room –WO66-G609 Silver Spring, MD 20993-0002

Hutchinson Santé S.N.C. Dr. Raffi Krikorian Vice President Science & Regulatory Rue Marret et Paturel Liancourt, France F-60140 SEP 1 3 2012

Re: K121335

Trade/Device Name: G-DermTM Powder Free Surgical Glove made of Styrene-

Ethylene/Butylene-Styrene (SEBS) and Styrene-Ethylene/ Propylene-Styrene-Ethylene/Propylene (SEPSEP) Synthetic

Copolymers

Regulation Number: 21 CFR 878.4460 Regulation Name: Surgeon's Glove

Regulatory Class: I Product Code: KGO Dated: August 28, 2012 Received: August 29, 2012

Dear Dr. Krikorian:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898.

In addition, FDA may publish further announcements concerning your device in the <u>Federal</u> <u>Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm.

Sincerely yours,

Anthony D. Watson, B.S., M.S., M.B.A.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

4. INDICATIONS FOR USE STATEMENT

510(k) Number (if known):K121335
Device Name: G-Derm TM Powder Free Surgical Glove made of styrene-ethylene/butylene-styrene (SEBS) and styrene-ethylene/propylene-styrene-ethylene/propylene (SEPSEP) synthetic copolymers
Indications for Use:
The G-Derm TM Surgical Glove is a disposable sterile powder-free medical device made of styrene-ethylene/butylene-styrene (SEBS) and styrene-ethylene/propylene-styrene-ethylene/propylene (SEPSEP) synthetic copolymers that is intended to be worn by operating room personnel to protect a surgical wound from contamination.
Prescription Use AND/OR Over-The-Counter Use X (21 CFR 801 Subpart C)
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)
Majur Jangulur
(Division Sign-Off) Division of Anesthesiology, General Hospital Infection Control, Dental Devices
510(k) Number: <u>K 12 1335</u> Page 1 of 1
1 age 1 of 1