



K123680

FEB 21 2013

## 510(k) Summary

owner's name:	White Peaks Dental Systems GmbH & Co. KG
address:	Langeheide 9 45239 Essen Germany
Establishment registration no.:	3009409473
Contact person:	Oliver Puckert
Phone / Fax:	+49-281-206458-0 / +49-281-206458-13
Name of the device:	CopranColor
Trade or proprietary name:	CopranColor
Classification name:	Powder, Porcelain
Device classification:	2
Regulation / Product code:	21 CFR 872.6660 / EIH

### Predicate device equivalence:

CopranColor equals the predicate device with respect to the indications for use and the fundamental technology, relating to the use of Zirkonium dioxide blanks, the application process (CAD/CAM technology), the biocompatibility, coloring system and other technical aspects.

### Device Description

CopranColor is a device that may be used as an accessory to zirkonium dioxide dental restorations such as Copran ZR /Origin YZ.

Together with the zirkonium dioxide blanks the CopranColor can be used by professional dental technicians for fabrication of all ceramic restorations.



### **Indications for Use**

*CopranColor is intended to be used as the coloring agent for the shading of white zirconia ceramic materials that are intended to be used by professional dental technicians to fabricate all-ceramic dental restorations for the sole use of particular patients.*

### **Comparison of Technology Characteristics**

The CopranColor is a liquid including colorants that is used prior to the final sintering of the restoration enabling the user to adjust the restoration to match the natural color of the patient's teeth.

The colouring liquid has no negative effect to the physical material properties.

This includes the variants:

- Copran Color (used for white Zirconium)
- Copran Ultra-T Color ( used for ceramic blending)
- Copran Monolith Color (used for milled teeth)

### **Determination of Substantial Equivalence**

CopranColor in combination with the zirconia blanks has the same technological characteristics as the predicate device ZENOTEC Zr / ZENOSTAR Color Zr.

The new device has the following similarities to the previously cleared predicate device:

- Same intended use
- Same operating principles
- Same technology
- Same manufacturing process

The materials used in the new devices are all biocompatible as the predicate device.

### **Conclusion drawn from nonclinical / clinical test:**

Based on the above and the fact that main component CopranZR is cleared for the US market since 2009 we consider the device CopranColor being safe and effective and we conclude the device is substantially equivalent to the predicate device: ZENOTEC Zr / ZENOSTAR Color Zr.



**Summary of Device Testing:**

Bench testing was performed to ensure that the CopranColor met its specifications. All tests were verified to meet acceptance criteria. Biocompatibility testing was performed to verify the equivalent safety of the materials that are used.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

February 21, 2013

Mr. Oliver Puckert  
CEO  
White Peakes Dental Systems GmbH & Co. KG  
Langeheide 9  
45239 Essen  
GERMANY

Re: K123680  
Trade/Device Name: CopranColor  
Regulation Number: 21 CFR 872.6660  
Regulation Name: Porcelain Powder for Clinical Use  
Regulatory Class: II  
Product Code: EIH  
Dated: November 13, 2012  
Received: November 30, 2012

Dear Mr. Puckert:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mary S. Runner  
*Mary S. Runner, DDS, MA* 2013.02.21  
12:03:50 -05'00'

FOR Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

***Indications for Use***

510(k) Number (if known): K123680

Device Name: CopranColor \_\_\_\_\_

CopranColor is intended to be used as the coloring agent for the shading of white zirconia ceramic materials that are intended to be used by professional dental technicians to fabricate all-ceramic dental restorations for the sole use of particular patients.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER  
PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page   1   of   1  

Mary S. Runner  
*Susan Runner, DDS, MA* 2013:02.21 09:56:39  
05'00"

(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K123680