

**MAY 22 2014****SECTION 5 - 510(K) SUMMARY**

<b>VIGILEO™ MONITOR 510(k) SUMMARY</b>	
<b>510(k) Submitter</b>	Edwards Lifesciences
<b>Contact Person</b>	Sally L. Maher One Edwards Way Irvine, CA 92614-5686
<b>Date Prepared</b>	May 6, 2014
<b>Trade Name</b>	Vigileo™ Arterial Pressure Cardiac Output (APCO)/Oximetry Monitor
<b>Common Name</b>	Cardiac Output/Oximetry Computer
<b>Classification Name</b>	Single-Function, Preprogrammed Diagnostic Computer (21 CFR 870.1435, product codes DXG, DQE )
<b>Regulation Class/ Product Code</b>	Class II/ DXG, DQE
<b>Predicate Device</b>	K103094 (cleared 17 March 2011), Edwards Lifesciences, Vigileo™ Arterial Pressure Cardiac Output/Oximetry Monitor
<b>Device Description</b>	<p>The Vigileo™ Arterial Pressure Cardiac Output (APCO)/Oximetry Monitor (Vigileo™ Monitor) is a microprocessor-based instrument. When used with an Edwards FloTrac sensor, the Vigileo™ Monitor is capable of continuously monitoring the following parameters:</p> <ul style="list-style-type: none"> <li>• Cardiac Output (CO);</li> <li>• Cardiac Index (CI);</li> <li>• Stroke Volume (SV);</li> <li>• Stroke Volume Variation (SVV);</li> <li>• Systemic Vascular Resistance (SVR);</li> <li>• Systemic Vascular Resistance Index (SVRI);</li> <li>• Oxygen Delivery (DO<sub>2</sub>); and</li> <li>• Oxygen Delivery Index (DO<sub>2</sub>I).</li> </ul> <p>When used with Edwards Oximetry catheters, the Vigileo™ Monitor is capable of continuously monitoring the following parameters:</p> <ul style="list-style-type: none"> <li>• Central venous oxygen saturation (ScvO<sub>2</sub>); and,</li> <li>• Mixed venous oxygen saturation (SvO<sub>2</sub>).</li> </ul>

<b>VIGILEO™ MONITOR 510(k) SUMMARY, con't.</b>	
<b>Indications for Use/ Intended Use</b>	The Vigileo™ APCO/Oximetry Monitor is indicated for continuously measuring hemodynamic parameters such as cardiac output and oximetry to assess oxygen delivery and consumption. When connected to an Edwards oximetry catheter, the monitor measures oximetry in adults and pediatrics. The monitor also displays parameters, such as stroke volume and stroke volume variation, used to assess fluid status and vascular resistance. The Vigileo™ APCO/Oximetry Monitor may be used in all setting in which critical care is provided.
<b>Comparative Analysis</b>	Verification and validation testing was conducted to compare the performance and functionality of the subject and predicate devices. This testing regimen included side-by-side bench, pre-clinical studies, and comparative analysis of clinical data. The Vigileo™ Monitor has been shown to be safe, effective, and substantially equivalent to the cited predicate device for its intended use critical care environments.
<b>Functional/ Safety Testing</b>	The Vigileo™ Monitor has successfully passed functional and performance testing, including software verification and validation, mechanical and electrical testing, bench studies, pre-clinical animal studies, comparison testing of clinical cases, and clinical utility.
<b>Conclusion</b>	The Vigileo™ Monitor has been shown to be safe, effective, and is substantially equivalent to the cited predicate devices for their intended use in the OR and ICU environments.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

May 22, 2014

Edwards Lifesciences LLC  
Mr. Stephen M. Enos  
Director, Regulatory Affairs  
One Edwards Way  
Irvine, CA 92614

Re: K131588  
Edwards Lifesciences™ Vigileo™ Arterial Pressure Cardiac Output/Oximetry  
Monitor  
Regulation Number: 21 CFR 870.1435  
Regulation Name: Single-function, Preprogrammed Diagnostic Computer  
Regulatory Class: Class II (two)  
Product Code: DXG, DQE  
Dated: May 6, 2014  
Received: May 7, 2014

Dear Mr. Enos:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

Page 2 – Mr. Stephen M. Enos

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>

Sincerely yours,

A stylized, blocky signature of Bram D. Zuckerman, M.D., written in a dark ink. The signature is positioned over the FDA logo.

for Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K131588

Device Name

Edwards Lifesciences™ Vigileo™ Arterial Pressure Cardiac Output/Oximetry Monitor

Indications for Use (Describe)

The Vigileo™ APCO/Oximetry Monitor is indicated for continuously measuring hemodynamic parameters such as cardiac output and oximetry to assess oxygen delivery and consumption. When connected to an Edwards oximetry catheter, the monitor measures oximetry in adults and pediatrics.

The monitor also displays parameters, such as stroke volume and stroke volume variation, used to assess fluid status and vascular resistance. The Vigileo™ APCO/Oximetry Monitor may be used in all settings in which critical care is provided.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)



**FD** Date: 2014.05.22 11:18:12  
**0400'**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*