



**510(k) Summary**  
**ArthroCare® Corporation**  
**ArthroCare® Coblator IQ™ Perc-D® SpineWands®**

**General Information**

**Submitter Name/Address:** ArthroCare Corporation  
7000 West William Cannon Dr.  
Austin, TX 78735

**Establishment Registration Number:** 3006524618

**Contact Person:** Shirley Hyink  
Sr. Manager, Regulatory Affairs

**Date Prepared:** July 2, 2013

**Device Description**

**Trade Name:** ArthroCare® Coblator IQ™ Perc-D®  
SpineWand® (Coblator IQ DLR and DLG  
SpineWands)

**Generic/Common Name:** Electrosurgical Device and Accessories

**Classification Name:** Electrosurgical Cutting and Coagulation  
Device and Accessories (21 CFR 878.4400)

**Predicate Devices**

K100353 - ArthroCare Coblator IQ Perc-D SpineWand (May 13, 2010)

**Product Description**

The Coblator IQ Perc-D SpineWands consist of two Wands, the Coblator IQ DLR SpineWand and the Coblator IQ DLG SpineWand, intended for use only with the Coblator IQ Controller.

These wands are single use, disposable, high frequency electrosurgery devices intended for use in percutaneous spinal disc decompression surgery. The Wands consist of a stainless steel shaft with a straight electrode configuration design. The proximal end of the devices is housed in a handle composed of polycarbonate. A cable integrated into the handle connects the Wand to the ArthroCare Coblator IQ Controller.

The Wand electrode is activated using the CIQ Controller System energy which starts the Coblation process to ablate the disc tissue. The Coblation process involves use of a conductive media to create a plasma layer to ablate tissue.

**Intended Use/Indications For Use**

The ArthroCare Coblator IQ Perc-D Spine Wands (Coblator IQ DLR and Coblator IQ DLG Spine Wands) are indicated for ablation, coagulation, and decompression of disc material to treat symptomatic patients with contained herniated lumbar and lumbosacral discs. The Wands are designed to be used exclusively with the ArthroCare Coblator IQ Controller.

**Non-Clinical Data**

This submission seeks only to clarify the indication for use statement of the predicate device and there have been no other device, design or material changes since the previous clearance of these predicate devices (K100353), new testing was deemed to be unnecessary. The performance testing previously conducted as part of the design control activities include device acceptance criteria testing, biological testing, and packaging testing.

**Clinical Data**

No clinical data are included in this submission.

**Summary**

The purpose of this submission is to clarify the current indication for use statement to specifically denote use of the Wands within the lumbar and lumbosacral spine. Other than the change to the indication of use statement the Wands are identical to the predicate devices in:

- Performance specification,
- Materials,
- Technological characteristics,
- Sterilization, and
- Principle of operation

The proposed Indications for Use do not significantly affect the safety or efficacy of the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

August 22, 2013

Shirley Hyink  
Senior Manager, Regulatory Affairs  
ArthroCare Corporation  
7000 West William Cannon Drive  
Austin, Texas 78735

Re: K132099

Trade/Device Name: ArthroCare<sup>®</sup> Coblator IQ<sup>™</sup> Perc-D<sup>®</sup> SpineWand<sup>®</sup>  
(Coblator IQ DLR and DLG SpineWands)  
Regulation Number: 21 CFR 878.4400  
Regulation Name: Electrosurgical cutting and coagulation device and accessories  
Regulatory Class: Class II  
Product Code: GEI  
Dated: July 3, 2013  
Received: July 10, 2013

Dear Ms. Hyink:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Mark N. Melkerson -S**

Mark N. Melkerson  
Acting Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

### Indications for Use

510(k) Number (if known): K132099

**Device Name:** ArthroCare® Coblator IQ™ Perc-D® SpineWand® (Coblator IQ DLR and DLG SpineWands)

### Indications for Use:

The ArthroCare Coblator IQ Perc-D Spine Wands (Coblator IQ DLR and Coblator IQ DLG SpineWands) are indicated for ablation, coagulation, and decompression of disc material to treat symptomatic patients with contained herniated lumbar and lumbosacral discs. The Wands are designed to be used exclusively with the ArthroCare Coblator IQ Controller.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Joshua C. Nipper -S**

(Division Sign-off)  
Division of Surgical Devices  
510(k) Number K132099