



NOV - 5 2013

510(k) SUMMARY
GRYPHON™ ANCHOR w/PROKNOT™ TECHNOLOGY

| | | |
|-------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|
| Date Summary Prepared | July 17, 2013 | |
| Submitter's Name and Address | DePuy Mitek a Johnson & Johnson company 325 Paramount Drive Raynham, MA 02767 | |
| Contact Person | Yayoi Fujimaki Regulatory Affairs Senior Associate DePuy Mitek a Johnson & Johnson company 325 Paramount Drive Raynham, MA 02767, USA | Telephone: 508-828-3541 Facsimile: 508-977-6911 e-mail: yfujimal@its.jnj.com |
| Name of Medical Device | Trade Name: GRYPHON™ Anchor w/PROKNOT™ Technology Common Name: fastener, fixation, biodegradable, soft tissue | |
| Device Classification | <ul style="list-style-type: none"> ▪ MAI - Single/multiple component metallic bone fixation appliances and accessories, classified as Class II, regulated per 21 CFR 888.3030. ▪ HWC - Smooth or threaded metallic bone fixation fastener, Classified as Class II, regulated per CFR 888.3040. ▪ Orthopedic panel | |
| Predicate Device | <ul style="list-style-type: none"> ▪ BIORAPTOR® Knotless Suture Anchor (K121018, K093428) ▪ GRYPHON™ P BR Anchor (K090124, K100012) ▪ GRYPHON™ PEEK Anchor (K103712) | |
| Indications for Use | Shoulder: Bankart Repair, SLAP Lesion Repair, Capsular Shift or Capsulolabral Reconstruction Hip: Capsular Repair, Acetabular Labral Repair | |
| Device Description | The proposed device is a suture-anchor to be used for soft tissue fixation to bone. The Gryphon Anchor is a cannulated, ribbed anchor, made of either non-absorbable PEEK (Polyetheretherketone) or absorbable Biocryl® Rapide™ (composite of β -TCP and PLGA copolymer). Orthocord suture (absorbable polydioxanone (PDS®) and non-absorbable polyethylene braided composite suture) is preloaded on the anchor. The suture incorporates a pre-tied sliding knot (ProKnot knot). The device is provided as sterile; the device is for single patient use only. | |

| | |
|----------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Comparison of Technological Characteristics | Substantial equivalence to the predicate devices has been justified by similarity of indications, design, material, operation principle and device performance data. Performance testing ensured that there is no new issue of safety and efficacy. |
| Safety and Performance | <p>Non-clinical Testing</p> <p>Performance requirement of the proposed device is to secure soft tissue to bone to heal. Fixation force testing was conducted under <i>in vitro</i> condition throughout two times of healing period, and substantial equivalence of fixation performance to the predicate device has been confirmed. Material biocompatibility has been also confirmed. The proposed device raises no new issue of safety and efficacy.</p> |
| Substantial Equivalence | Based on technological characteristics comparison and performance evaluation, the proposed device is concluded to be substantially equivalent to the predicate devices. |



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

November 5, 2013

DePuy Mitek, a Johnson and Johnson company
Ms. Yayoi Fujimaki
Regulatory Affairs Senior Associate
325 Paramount Drive
Raynham, Massachusetts 02767

Re: K132241

Trade/Device Name: GRYPHON™ ANCHOR w/PROKNOT™ TECHNOLOGY
Regulation Number: 21 CFR 888.3030
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories
Regulatory Class: Class II
Product Code: MAI, HWC
Dated: September 25, 2013
Received: September 26, 2013

Dear Ms. Fujimaki:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

ErinFDKeith

for

Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K132241

Device Name: GRYPHON™ ANCHOR w/PROKNOT™ TECHNOLOGY

Indications for Use:

GRYPHON ANCHOR w/PROKNOT TECHNOLOGY is indicated for followings.

Shoulder: Bankart Repair, SLAP Lesion Repair, Capsular Shift or Capsulolabral Reconstruction

Hip: Capsular Repair, Acetabular Labral Repair

Prescription Use x

AND/OR

Over-The-Counter Use

(Part 21 CFR 801 Subpart D)

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Casey L. Hanley, Ph.D.

Division of Orthopedic Devices

Page 1 of 1