

MAR 06 2014

510(k) SUMMARY
Amendia's Apollo Suture Anchor and Titan Screws

Date: September 11, 2013
Contact: Daniel Lanois, General Manager
Trade Name: Apollo Suture Anchor System and Titan Screws
Common Name: Screw, Fixation, Bone
Product Class: Class II
Classification: 21 CFR §888.3040 – Smooth or threaded metallic bone fixation fastener
Product Code: MBI
Panel Code: 87

Name/Address of Sponsor

Amendia, LLC
1755 West Oak Pkwy
Marietta, GA 30062
770-575-5200

Purpose:

The purpose of this submission is clearance of the Apollo Suture Anchor System and Titan Screws as a new medical device that is substantially equivalent to the predicate device.

Device Description

Apollo Family

The Apollo Medial Suture Anchor, Lateral Anchor, and Suture Tac Delivery Systems are delivery systems for anchors for use in fixation of ligament, tendon, bone, or soft tissue to bone in knee, shoulder, foot/ankle, elbow, and hand/wrist procedures. These anchors consist of cannulated anchors with integrated suture attachment or separate suture punch eyelet. The Anchors are provided loaded on individual inserters with and without integrated sutures, sterile, for single use only.

Titan Family

The Titan Interference and Tenodesis Screws are interference screws for use in fixation of ligament, tendon, bone, or soft tissue to bone in knee, shoulder, foot/ankle, elbow, and hand/wrist procedures. The screws are provided sterile, for single use only.

Screw and anchor implants are made from either a titanium alloy (6Al4V) per ASTM F136, or PEEK (Zeniva ZA-500) per ASTM F2026 from Solvay Advanced Polymers.

Predicate Device

The Apollo Suture Anchor System and Titan Screws are substantially equivalent to Smith and Nephew Biosure PK Interference Screw (K083635), the DePuy Mitek Healix PEEK Anchor (K120449), the Biomet PEEK Allthread Knotless Suture Anchor and Hitch PEEK Anchor (K071569 and K060693), and the Arthrex Tenodesis Screw (K041356).

Intended Use / Indications for Use

The Apollo Suture Anchors and Titan Screws are indicated for use in fixation of ligament, tendon, bone, or soft tissue to bone in knee, shoulder, foot/ankle, elbow, and hand/wrist procedures. The screws are intended for use in the following procedures:

<p>Indications – Apollo Medial Suture Anchor The Apollo Medial Suture Anchor is intended for:</p> <p>Shoulder</p> <ul style="list-style-type: none"> <input type="checkbox"/> Rotator Cuff Repair <input type="checkbox"/> Bankart Repair <input type="checkbox"/> SLAP Lesion Repair <input type="checkbox"/> Biceps Tenodesis <input type="checkbox"/> Acromioclavicular Separation Repair <input type="checkbox"/> Deltoid Repair <input type="checkbox"/> Capsular Shift or Capsulolabral Reconstruction <p>Foot/Ankle</p> <ul style="list-style-type: none"> <input type="checkbox"/> Lateral Stabilization <input type="checkbox"/> Medial Stabilization <input type="checkbox"/> Achilles Tendon Repair <p>Knee</p> <ul style="list-style-type: none"> <input type="checkbox"/> Medial Collateral Ligament Repair <input type="checkbox"/> Lateral Collateral Ligament Repair <input type="checkbox"/> Posterior Oblique Ligament Repair <input type="checkbox"/> Iliotibial Band Tenodesis <p>Elbow</p> <ul style="list-style-type: none"> <input type="checkbox"/> Biceps Tendon Reattachment • <input type="checkbox"/> Ulnar or Radial Collateral Ligament Reconstruction <p>Hip</p> <ul style="list-style-type: none"> <input type="checkbox"/> Capsular Repair <input type="checkbox"/> Acetabular Labral Repair 	<p>Indications – Apollo Lateral Anchor The Apollo Lateral Anchor is indicated for:</p> <p>Shoulder</p> <ul style="list-style-type: none"> <input type="checkbox"/> Rotator Cuff Repair <input type="checkbox"/> Bankart Repair <input type="checkbox"/> SLAP Lesion Repair <input type="checkbox"/> Biceps Tenodesis <input type="checkbox"/> Acromioclavicular Separation Repair <input type="checkbox"/> Deltoid Repair <input type="checkbox"/> Capsular Shift or Capsulolabral Reconstruction <p>Wrist/Hand</p> <ul style="list-style-type: none"> <input type="checkbox"/> Scapholunate Ligament Reconstruction <input type="checkbox"/> Ulnar/Radial Collateral Ligament Reconstruction <p>Foot/Ankle</p> <ul style="list-style-type: none"> <input type="checkbox"/> Lateral Stabilization <input type="checkbox"/> Medial Stabilization <input type="checkbox"/> Achilles Tendon Repair/Reconstruction <input type="checkbox"/> Hallux Valgus Reconstruction <input type="checkbox"/> Mid and Forefoot Reconstruction <p>Elbow</p> <ul style="list-style-type: none"> <input type="checkbox"/> Biceps Tendon Reconstruction <input type="checkbox"/> Ulnar or Radial Collateral Ligament Reconstruction <input type="checkbox"/> Lateral Epicondylitis Repair (PEEK Anchor Only) <p>Knee</p> <ul style="list-style-type: none"> <input type="checkbox"/> Medial Collateral Ligament Repair <input type="checkbox"/> Lateral Collateral Ligament Repair <input type="checkbox"/> Posterior Oblique Ligament Repair <input type="checkbox"/> Joint Capsule Closure <input type="checkbox"/> Iliotibial Band Tenodesis <input type="checkbox"/> Patellar Ligament/Tendon Repair
<p>Indications – Apollo Labral Suture Tac Shoulder</p> <ul style="list-style-type: none"> <input type="checkbox"/> Rotator Cuff Repair <input type="checkbox"/> Bankart Repair <input type="checkbox"/> SLAP Lesion Repair 	

<input type="checkbox"/> Biceps Tenodesis <input type="checkbox"/> Acromioclavicular Separation Repair <input type="checkbox"/> Deltoid Repair <input type="checkbox"/> Capsular Shift or Capsulolabral Reconstruction Wrist <input type="checkbox"/> Scapholunate Ligament Reconstruction Elbow <input type="checkbox"/> Biceps Tendon Reattachment <input type="checkbox"/> Ulnar or Radial Collateral Ligament Reconstruction Hip <input type="checkbox"/> Capsular Repair <input type="checkbox"/> Acetabular Labral Repair Knee <input type="checkbox"/> Extracapsular Repair <input type="checkbox"/> Medial Collateral Ligament Repair <input type="checkbox"/> Lateral Collateral Ligament Repair <input type="checkbox"/> Posterior Oblique Ligament Repair <input type="checkbox"/> Joint Capsule Closure <input type="checkbox"/> Iliotibial Band Tenodesis Reconstruction <input type="checkbox"/> Patellar Ligament/Tendon Repair <input type="checkbox"/> Vastus Medialis Obliquus Muscle Advancement	
Indications –Interference Screws <p>The Titan Interference Screws are indicated for the reattachment of ligament, tendon, soft tissue, or bone to bone during cruciate ligament reconstruction surgeries of the knee. All screws with a diameter of 9 mm or less and a length of 23 mm or less are also intended for the use in the following procedures:</p> <p>Knee</p> <input type="checkbox"/> ACL repairs <input type="checkbox"/> PCL repairs <input type="checkbox"/> Extracapsular repairs o Medial collateral ligament o Lateral collateral ligament o Posterior oblique ligament <input type="checkbox"/> Patellar realignment and tendon repairs o Vastus medialis obliquus advancement <input type="checkbox"/> Iliotibial band tenodesis <p>Shoulder</p> <input type="checkbox"/> Capsular stabilization o Bankart repair o Anterior shoulder instability o SLAP lesion repairs o Capsular shift of capsulolabral reconstructions <input type="checkbox"/> Acromioclavicular separation repairs <input type="checkbox"/> Deltoid repairs <input type="checkbox"/> Rotator cuff tear repairs <input type="checkbox"/> Biceps tenodesis <p>Foot and Ankle</p> <input type="checkbox"/> Hallux valgus repairs <input type="checkbox"/> Medial or lateral instability repairs/reconstructions <input type="checkbox"/> Achilles tendon repairs/reconstructions <input type="checkbox"/> Midfoot reconstructions	Indications –Tenodesis Screws <p>The Titan Tenodesis Screws are intended to be used for fixation of tissue, including ligament or tendon to bone, or a bone/tendon to bone. See below for specific indications. The Tenodesis Screws are intended to provide soft tissue reattachment (i.e. fixation of ligament and tendon graft tissue). See below for specific indications.</p> <p>Shoulder</p> <input type="checkbox"/> Capsular stabilization o Bankart repair o Anterior shoulder instability o SLAP lesion repairs o Capsular shift of capsulolabral reconstructions <input type="checkbox"/> Acromioclavicular separation repairs <input type="checkbox"/> Deltoid repairs <input type="checkbox"/> Rotator cuff tear repairs <input type="checkbox"/> Biceps tenodesis <p>Foot and Ankle</p> <input type="checkbox"/> Hallux valgus reconstruction <input type="checkbox"/> Medial stabilization <input type="checkbox"/> Lateral stabilization <input type="checkbox"/> Achilles Tendon Repair <input type="checkbox"/> Midfoot reconstructions <input type="checkbox"/> Metatarsal ligament repair <input type="checkbox"/> Bunionectomy <input type="checkbox"/> Flexor Hallucis Longus for Achilles Tendon reconstruction <input type="checkbox"/> Tendon transfers in the foot and ankle <p>Knee</p> <input type="checkbox"/> Medial Collateral Ligament Repair <input type="checkbox"/> Lateral Collateral Ligament Repair

<input type="checkbox"/> Metatarsal ligament/tendon repairs/reconstructions <input type="checkbox"/> Bunionectomy <input type="checkbox"/> Flexor Hallucis Longus <input type="checkbox"/> Tendon transfers Elbow, Wrist, and Hand <input type="checkbox"/> Biceps tendon reattachment <input type="checkbox"/> Ulnar or radial collateral ligament reconstructions <input type="checkbox"/> Lateral epicondylitis repair <input type="checkbox"/> Scapholunate ligament reconstruction <input type="checkbox"/> Tendon transfers	<input type="checkbox"/> Patellar Tendon Repair <input type="checkbox"/> Posterior Oblique Ligament Repair <input type="checkbox"/> Iliotibial Band Tenodesis <input type="checkbox"/> Posterior Cruciate Ligament Repair Elbow <input type="checkbox"/> Biceps tendon reattachment <input type="checkbox"/> Ulnar or radial collateral ligament reconstruction Wrist and Hand <input type="checkbox"/> Scapholunate Ligament Reconstruction <input type="checkbox"/> Ulnar Collateral Ligament Reconstruction <input type="checkbox"/> Radial Collateral Ligament Reconstruction <input type="checkbox"/> Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty) <input type="checkbox"/> Carpal Ligament Reconstructions and repairs <input type="checkbox"/> Tendon transfer in the hand/wrist <input type="checkbox"/> Lateral Epicondylitis repair
---	--

Performance Data:

Performance testing included axial pullout, axial push out, insertion torque, and torque to failure tests per ASTM F543-07. Side-by-side testing to predicate devices was completed with test results demonstrating that the Apollo Suture Anchor System and Titan Screws are substantially equivalent to the predicate devices.

Summary:

The Apollo Suture Anchor System and Titan Screws are substantially equivalent to the predicate devices in regards to:

- Indications for Use
- Materials
- Mechanism of action
- Dimensions

There are no significant differences in technological characteristics compared to the predicate device that would raise any new types of safety or efficacy issues.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

March 6, 2014

Amendia, LLC
% Dr. Richard Jansen, Pharm. D.
President
Silver Pine Consulting
13540 Guild Ave.
Apple Valley, Minnesota 55124

Re: K133036

Trade/Device Name: Apollo Suture Anchor System and Titan Screws
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: MBI
Dated: January 17, 2014
Received: February 14, 2014

Dear Dr. Jansen,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Vincent ~~FD~~ Devlin -S

for
Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known): K133036

Device Name: Apollo Suture Anchor System and Titan Screws

Intended Use / Indications for Use

The Apollo Suture Anchors and Titan Screws are indicated for use in fixation of ligament, tendon, bone, or soft tissue to bone in knee, shoulder, foot/ankle, elbow, and hand/wrist procedures. The screws are intended for use in the following procedures:

<p>Indications – Apollo Medial Suture Anchor The Apollo Medial Suture Anchor is intended for:</p> <p>Shoulder</p> <p><input type="checkbox"/> Rotator Cuff Repair</p> <p><input type="checkbox"/> Bankart Repair</p> <p><input type="checkbox"/> SLAP Lesion Repair</p> <p><input type="checkbox"/> Biceps Tenodesis</p> <p><input type="checkbox"/> Acromioclavicular Separation Repair</p> <p><input type="checkbox"/> Deltoid Repair</p> <p><input type="checkbox"/> Capsular Shift or Capsulolabral Reconstruction</p> <p>Foot/Ankle</p> <p><input type="checkbox"/> Lateral Stabilization</p> <p><input type="checkbox"/> Medial Stabilization</p> <p><input type="checkbox"/> Achilles Tendon Repair</p> <p>Knee</p> <p><input type="checkbox"/> Medial Collateral Ligament Repair</p> <p><input type="checkbox"/> Lateral Collateral Ligament Repair</p> <p><input type="checkbox"/> Posterior Oblique Ligament Repair</p> <p><input type="checkbox"/> Iliotibial Band Tenodesis</p> <p>Elbow</p> <p><input type="checkbox"/> Biceps Tendon Reattachment</p> <p><input type="checkbox"/> Ulnar or Radial Collateral Ligament Reconstruction</p> <p>Hip</p> <p><input type="checkbox"/> Capsular Repair</p> <p><input type="checkbox"/> Acetabular Labral Repair</p>	<p>Indications – Apollo Lateral Anchor The Apollo Lateral Anchor is indicated for:</p> <p>Shoulder</p> <p><input type="checkbox"/> Rotator Cuff Repair</p> <p><input type="checkbox"/> Bankart Repair</p> <p><input type="checkbox"/> SLAP Lesion Repair</p> <p><input type="checkbox"/> Biceps Tenodesis</p> <p><input type="checkbox"/> Acromioclavicular Separation Repair</p> <p><input type="checkbox"/> Deltoid Repair</p> <p><input type="checkbox"/> Capsular Shift or Capsulolabral Reconstruction</p> <p>Wrist/Hand</p> <p><input type="checkbox"/> Scapholunate Ligament Reconstruction</p> <p><input type="checkbox"/> Ulnar/Radial Collateral Ligament Reconstruction</p> <p>Foot/Ankle</p> <p><input type="checkbox"/> Lateral Stabilization</p> <p><input type="checkbox"/> Medial Stabilization</p> <p><input type="checkbox"/> Achilles Tendon Repair/Reconstruction</p> <p><input type="checkbox"/> Hallux Valgus Reconstruction</p> <p><input type="checkbox"/> Mid and Forefoot Reconstruction</p> <p>Elbow</p> <p><input type="checkbox"/> Biceps Tendon Reconstruction</p> <p><input type="checkbox"/> Ulnar or Radial Collateral Ligament Reconstruction</p> <p><input type="checkbox"/> Lateral Epicondylitis Repair (PEEK Anchor Only)</p> <p>Knee</p> <p><input type="checkbox"/> Medial Collateral Ligament Repair</p> <p><input type="checkbox"/> Lateral Collateral Ligament Repair</p> <p><input type="checkbox"/> Posterior Oblique Ligament Repair</p> <p><input type="checkbox"/> Joint Capsule Closure</p> <p><input type="checkbox"/> Iliotibial Band Tenodesis</p> <p><input type="checkbox"/> Patellar Ligament/Tendon Repair</p>
<p>Indications – Apollo Labral Suture Tac Shoulder</p> <p><input type="checkbox"/> Rotator Cuff Repair</p> <p><input type="checkbox"/> Bankart Repair</p> <p><input type="checkbox"/> SLAP Lesion Repair</p> <p><input type="checkbox"/> Biceps Tenodesis</p> <p><input type="checkbox"/> Acromioclavicular Separation Repair</p> <p><input type="checkbox"/> Deltoid Repair</p> <p><input type="checkbox"/> Capsular Shift or Capsulolabral Reconstruction</p> <p>Wrist</p> <p><input type="checkbox"/> Scapholunate Ligament Reconstruction</p> <p>Elbow</p> <p><input type="checkbox"/> Biceps Tendon Reattachment</p> <p><input type="checkbox"/> Ulnar or Radial Collateral Ligament Reconstruction</p> <p>Hip</p> <p><input type="checkbox"/> Capsular Repair</p> <p><input type="checkbox"/> Acetabular Labral Repair</p>	

K133036

Knee <input type="checkbox"/> Extracapsular Repair <input type="checkbox"/> Medial Collateral Ligament Repair <input type="checkbox"/> Lateral Collateral Ligament Repair <input type="checkbox"/> Posterior Oblique Ligament Repair <input type="checkbox"/> Joint Capsule Closure <input type="checkbox"/> Iliotibial Band Tenodesis Reconstruction <input type="checkbox"/> Patellar Ligament/Tendon Repair <input type="checkbox"/> Vastus Medialis Obliquus Muscle Advancement	
Indications –Interference Screws The Titan Interference Screws are indicated for the reattachment of ligament, tendon, soft tissue, or bone to bone during cruciate ligament reconstruction surgeries of the knee. All screws with a diameter of 9 mm or less and a length of 23 mm or less are also intended for the use in the following procedures: Knee <input type="checkbox"/> ACL repairs <input type="checkbox"/> PCL repairs <input type="checkbox"/> Extracapsular repairs o Medial collateral ligament o Lateral collateral ligament o Posterior oblique ligament <input type="checkbox"/> Patellar realignment and tendon repairs o Vastus medialis obliquus advancement <input type="checkbox"/> Iliotibial band tenodesis Shoulder <input type="checkbox"/> Capsular stabilization o Bankart repair o Anterior shoulder instability o SLAP lesion repairs o Capsular shift of capsulolabral reconstructions <input type="checkbox"/> Acromioclavicular separation repairs <input type="checkbox"/> Deltoid repairs <input type="checkbox"/> Rotator cuff tear repairs <input type="checkbox"/> Biceps tenodesis Foot and Ankle <input type="checkbox"/> Hallux valgus repairs <input type="checkbox"/> Medial or lateral instability repairs/reconstructions <input type="checkbox"/> Achilles tendon repairs/reconstructions <input type="checkbox"/> Midfoot reconstructions <input type="checkbox"/> Metatarsal ligament/tendon repairs/reconstructions <input type="checkbox"/> Bunionectomy <input type="checkbox"/> Flexor Hallucis Longus <input type="checkbox"/> Tendon transfers Elbow, Wrist, and Hand <input type="checkbox"/> Biceps tendon reattachment <input type="checkbox"/> Ulnar or radial collateral ligament reconstructions <input type="checkbox"/> Lateral epicondylitis repair <input type="checkbox"/> Scapholunate ligament reconstruction <input type="checkbox"/> Tendon transfers	Indications –Tenodesis Screws The Titan Tenodesis Screws are intended to be used for fixation of tissue, including ligament or tendon to bone, or a bone/tendon to bone. See below for specific indications. The Tenodesis Screws are intended to provide soft tissue reattachment (i.e. fixation of ligament and tendon graft tissue). See below for specific indications. Shoulder <input type="checkbox"/> Capsular stabilization o Bankart repair o Anterior shoulder instability o SLAP lesion repairs o Capsular shift of capsulolabral reconstructions <input type="checkbox"/> Acromioclavicular separation repairs <input type="checkbox"/> Deltoid repairs <input type="checkbox"/> Rotator cuff tear repairs <input type="checkbox"/> Biceps tenodesis Foot and Ankle <input type="checkbox"/> Hallux valgus reconstruction <input type="checkbox"/> Medial stabilization <input type="checkbox"/> Lateral stabilization <input type="checkbox"/> Achilles Tendon Repair <input type="checkbox"/> Midfoot reconstructions <input type="checkbox"/> Metatarsal ligament repair <input type="checkbox"/> Bunionectomy <input type="checkbox"/> Flexor Hallucis Longus for Achilles Tendon reconstruction <input type="checkbox"/> Tendon transfers in the foot and ankle Knee <input type="checkbox"/> Medial Collateral Ligament Repair <input type="checkbox"/> Lateral Collateral Ligament Repair <input type="checkbox"/> Patellar Tendon Repair <input type="checkbox"/> Posterior Oblique Ligament Repair <input type="checkbox"/> Iliotibial Band Tenodesis <input type="checkbox"/> Posterior Cruciate Ligament Repair Elbow <input type="checkbox"/> Biceps tendon reattachment <input type="checkbox"/> Ulnar or radial collateral ligament reconstruction Wrist and Hand <input type="checkbox"/> Scapholunate Ligament Reconstruction <input type="checkbox"/> Ulnar Collateral Ligament Reconstruction <input type="checkbox"/> Radial Collateral Ligament Reconstruction <input type="checkbox"/> Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty) <input type="checkbox"/> Carpal Ligament Reconstructions and repairs <input type="checkbox"/> Tendon transfer in the hand/wrist <input type="checkbox"/> Lateral Epicondylitis repair

Prescription Use v
 (Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use
 (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Casey L. Hanley, Ph.D.
 Division of Orthopedic Devices