

FEB 21 2014

k133710

510(k) Summary of Safety and Effectiveness

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SDMA 1990 and 21 CFR 807.92.

Introduction

According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

Submitter name, Address, and Contact

Lin-Zhi International, Inc.
670 Almanor Avenue
Sunnyvale, CA 94085
Phone: (408) 732-3856
Fax: (408) 732-3849
e-mail: bclin@lin-zhi.com

Contact: Bernice Lin, Ph.D.
VP Operations

Submission Date

February 21, 2014

Multi-calibrator Set A, B, and C Device Name and Classification

Multi-calibrator Set A:

Classification Name: Drug Mixture Calibrator Materials
Class II, DKB (91 Toxicology),
21 CFR862.3200

Common Name: Benzoylcegonine, Methamphetamine, Methadone, and 6-acetylmorphine Calibrators

Proprietary Name: LZI Multiple Analyte Set A Drugs of Abuse Calibrators

Classification Name: Drug Mixture Control Materials
Class I, DIF (91 Toxicology),
21 CFR 862.3280

Common Name: Benzoylcegonine, Methamphetamine, Methadone, and 6-acetylmorphine Controls

Proprietary Name: LZI Multiple Analyte Set A Drugs of Abuse Controls

Multi-calibrator Set B:

Classification Name: Drug Mixture Calibrator Materials
Class II, DKB (91 Toxicology),
21 CFR862.3200

Common Name: MDMA, Morphine, Oxazepam, and Secobarbital Calibrators

Proprietary Name: LZI Multiple Analyte Set B Drugs of Abuse Calibrators

Classification Name: Drug Mixture Control Materials
Class I, DIF (91 Toxicology),
21 CFR 862.3280

Common Name: MDMA, Morphine, Oxazepam, and Secobarbital Controls

Proprietary Name: LZI Multiple Analyte Set B Drugs of Abuse Controls

Multi-calibrator Set C:

Classification Name: Drug Mixture Calibrator Materials
Class II, DKB (91 Toxicology),
21 CFR862.3200

Common Name: Benzoylecgonine, Methamphetamine, Methadone, Morphine,
Oxazepam, and Secobarbital Calibrators

Proprietary Name: LZI Multiple Analyte Set C Drugs of Abuse Calibrators

Classification Name: Drug Mixture Control Materials
Class I, DIF (91 Toxicology),
21 CFR 862.3280

Common Name: Benzoylecgonine, Methamphetamine, Methadone, Morphine,
Oxazepam, and Secobarbital Controls

Proprietary Name: LZI Multiple Analyte Set C Drugs of Abuse Controls

Previous Submission Information

There were no prior submissions for this subject device (k133710).

Legally Marketed Predicate Device(s)

The LZI Multiple Analyte Set A, B, and C Drugs of Abuse Calibrators and Controls (k133710) are substantially equivalent to the Multiple Analyte Drugs of Abuse Calibrators and Controls (k051088) manufactured by Lin-Zhi International, Inc. The LZI Multiple Analyte Set A, B, and C Drugs of Abuse Calibrators and Controls are identical or similar to its predicate in terms of intended use, method principle, device components, and clinical performance.

Device Description

All of the LZI Multiple Analyte Set A, B, and C Drugs of Abuse Calibrators and Controls are liquid and ready to use. These Calibrators and Controls do not have any especially unique technical characteristics. Each contains a known concentration of a specific drug analyte as a mixture.

The Negative DAU Calibrator is a processed, drug-free human urine matrix in human urine with sodium azide (0.09%) as preservative. The Low, Cutoff, Intermediate, and High Calibrators, as well as the 2 levels of Controls are prepared by spiking known concentrations of drug analyte into the Negative DAU Calibrator matrix. These five calibrators and two controls are sold as individual bottles. The various concentrations of each drug analyte in their corresponding calibrators and controls are summarized as follows:

Table 1: LZI Multiple Analyte Set A Calibrators and Controls:

	Low Calibrator	Cutoff Calibrator	Intermediate Calibrator	High Calibrator	Control Level 1	Control Level 2
Material	ng/mL	ng/mL	ng/mL	ng/mL	ng/mL	ng/mL
Benzoylcegonine	75	150	300	1000	112.5	187.5
d-Methamphetamine	250	500	750	1000	375	625
Methadone	150	300	600	1000	225	375
6-AM	5	10	20	40	7.5	12.5

* 5th Calibrator is the Negative Calibrator

Table 2: LZI Multiple Analyte Set B Calibrators and Controls:

	Low Calibrator	Cutoff Calibrator	Intermediate Calibrator	High Calibrator	Control Level 1	Control Level 2
Material	ng/mL	ng/mL	ng/mL	ng/mL	ng/mL	ng/mL
MDMA	250	500	750	1000	375	625
Morphine	1000	2000	4000	6000	1500	2500
Oxazepam	100	200	500	1000	150	250
Secobarbital	100	200	500	1000	150	250

* 5th Calibrator is the Negative Calibrator

Table 3: LZI Multiple Analyte Set C Calibrators and Controls:

	Low Calibrator	Cutoff Calibrator	Intermediate Calibrator	High Calibrator		Control Level 1	Control Level 2
Material	ng/mL	ng/mL	ng/mL	ng/mL		ng/mL	ng/mL
Benzoylcegonine	75	150	300	1000		112.5	187.5
d-Methamphetamine	250	500	750	1000		375	625
Methadone	150	300	600	1000		225	375
Morphine	1000	2000	4000	6000		1500	2500
Oxazepam	100	200	500	1000		150	250
Secobarbital	100	200	500	1000		150	250

* 5th Calibrator is the Negative Calibrator

Intended Use

Device Name: LZI Multiple Analyte Set A, B, and C Urine Drugs of Abuse Calibrators

The LZI Multiple Analyte Set A, B, and C Urine Drugs of Abuse Calibrators are intended for in vitro diagnostic use for the calibration of assays for the analytes currently listed in the package insert: Benzoylecgonine, Methamphetamine, Methadone, 6-acetylmorphine, MDMA, Morphine, Oxazepam, and Secobarbital. The calibrators are designed for prescription use with homogeneous enzyme immunoassays on automated clinical chemistry analyzers.

Device Name: LZI Multiple Analyte Set A, B, and C Urine Drugs of Abuse Controls

The LZI Multiple Analyte Set A, B, and C Urine Drugs of Abuse Controls are intended for in vitro diagnostic use to monitor the performance of assays for the analytes currently listed in the package insert: Benzoylecgonine, Methamphetamine, Methadone, 6-acetylmorphine, MDMA, Morphine, Oxazepam, and Secobarbital. The controls are designed for prescription use with homogeneous enzyme immunoassays on automated clinical chemistry analyzers.

Comparison to Predicate Device

The LZI Multiple Analyte Set A, B, and C Drugs of Abuse Calibrators and Controls (k133710) are substantially equivalent to the Multiple Analyte Drugs of Abuse Calibrators and Controls (k051088) manufactured by Lin-Zhi International, Inc. The LZI Multiple Analyte Set A Drugs of Abuse Calibrators and Controls are identical or similar to its predicate in terms of intended use, method principle, device components, and clinical performance

All spiked values of calibrators and controls were confirmed with GC/MS. Performance characteristics on precision, accuracy and stability are acceptable. The following table compares LZI Multiple Analyte Set A, B, and C Drugs of Abuse Calibrators and Controls with the predicate device.

Device Characteristics	Subject Device (k133710) LZI Multiple Analyte Set A DAU Calibrators and Controls	Subject Device (k133710) LZI Multiple Analyte Set B DAU Calibrators and Controls	Subject Device (k133710) LZI Multiple Analyte Set C DAU Calibrators and Controls	Predicate Device (k051088) LZI Multiple Analyte DAU Calibrators and Controls
Intended Use	Intended for in vitro diagnostic use for the calibration and validation of LZI DAU enzyme immunoassays to detect benzoyllecgonine, methamphetamine, methadone, and 6-acetylmorphine, in human urine. <i>The assays used with the multiple analyte calibrators and controls only provide a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas or liquid chromatography/mass spectrometry (GC/MS or LC/MS) is the preferred confirmatory method. Clinical consideration and professional judgement should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.</i>	Intended for in vitro diagnostic use for the calibration and validation of LZI DAU enzyme immunoassays to detect MDMA, morphine, oxazepam, and secobarbital in human urine. <i>The assays used with the multiple analyte calibrators and controls only provide a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas or liquid chromatography/mass spectrometry (GC/MS or LC/MS) is the preferred confirmatory method. Clinical consideration and professional judgement should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.</i>	Intended for in vitro diagnostic use for the calibration and validation of LZI DAU enzyme immunoassays to detect benzoyllecgonine, methamphetamine, methadone, morphine, oxazepam, and secobarbital in human urine. <i>The assays used with the multiple analyte calibrators and controls only provide a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas or liquid chromatography/mass spectrometry (GC/MS or LC/MS) is the preferred confirmatory method. Clinical consideration and professional judgement should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.</i>	Intended for in vitro diagnostic use for the calibration and validation of LZI DAU enzyme immunoassays to detect methamphetamine, opiate, phenyleclidine, benzoyllecgonine, benzodiazepines, barbiturates, methadone, and propoxyphene in human urine.

Comparison to Predicate Device (continued)

Device Characteristics	Subject Device (k133710) LZI Multiple Analyte Set A DAU Calibrators and Controls	Subject Device (k133710) LZI Multiple Analyte Set B DAU Calibrators and Controls	Subject Device (k133710) LZI Multiple Analyte Set C DAU Calibrators and Controls	Predicate Device (k051088) LZI Multiple Analyte DAU Calibrators and Controls
Analyte	benzylecgonine, d-methamphetamine, methadone, 6-acetylmorphine	morphine, oxazepam, secobarbital, MDMA	benzylecgonine d-methamphetamine, methadone, morphine, oxazepam, secobarbital	benzylecgonine, d-methamphetamine, methadone, morphine, oxazepam, secobarbital, phenylethidine, propoxyphene
Materials are applicable for assays with the cutoff concentrations listed (ng/mL)	COC - 150 MAMP or AMP - 500 MTD - 300 6AM - 10	OPI - 2,000 BZO - 200 BARB - 200 MDMA - 500	COC - 150 MAMP or AMP - 500 MTD - 300 OPI - 2,000 BZO - 200 BARB - 200	COC - 150 MAMP or AMP - 500 MTD - 300 OPI - 2,000 BZO - 200 BARB - 200 PCP - 25 PPX - 300
Matrix	Urine	Urine	Urine	Urine
Calibrators Level	5 Levels - See Table 1 (under "Device Description")	5 Levels - See Table 2 (under "Device Description")	5 Levels - See Table 3 (under "Device Description")	5 Levels - See Table 4 (below)
Controls Level	2 Levels - See Table 1 (under "Device Description")	2 Levels - See Table 2 (under "Device Description")	2 Levels - See Table 3 (under "Device Description")	2 Levels - See Table 4 (below)
Storage	2-8 °C until expiration date	2-8 °C until expiration date	2-8 °C until expiration date	2-8 °C until expiration date

Table 4: Predicate Multiple Analyte Calibrators and Controls (k051088)

Material	Low Calibrator ng/mL	Cutoff Calibrator ng/mL	Intermediate Calibrator ng/mL	High Calibrator ng/mL	Control Level 1 ng/mL	Control Level 2 ng/mL
d-Methamphetamine	250	500	750	1000	375	625
Morphine	1000	2000	4000	6000	1500	2500
Phencyclidine	12.5	25	50	100	18	35
Benzoyllecgonine	75	150	300	1000	110	190
Oxazepam	100	200	500	1000	100	300
Secobarbital	100	200	500	1000	100	300
Propoxyphene	150	300	600	1000	225	375
Methadone	150	300	600	1000	225	375

* 5th Calibrator is the Negative Calibrator

Performance Characteristics Summary:

Multiple Analyte Sets A,B, and C Drugs of Abuse Calibrators and Controls

Stability:

Based on the largest % change for any individual concentration or % change in Total Separation (versus Day 0 measurements) at Cold (2-8 °C), Room Temperature, and Accelerated (30 °C) Temperature studies, the LZI Multiple Analyte Set A, B, and C Drugs of Abuse Calibrators and Controls did not vary greater than $\pm 10\%$. In Set A, the only analyte for which this differed was the 6- acetylmorphine (6AM) Calibrator and Control. In Set B and C, the only analyte for which this differed was the Benzodiazepine (Oxazepam) Calibrator and Control. All calibrators and controls held steady when stored at 2 °C to 8 °C (refrigerated) up to Day 366 real time. Based on these data we have capped our predicted shelf-life to the current cold temperature real-time study for the LZI Multiple Analyte Sets A, B, and C Drugs of Abuse Calibrators and Controls to at least 12 months.

Traceability and Value Assignment:

The starting materials for the calibrators and controls were commercially available standard stock solutions in methanol. The standard solution was found to be 99% in purity, analyzed by GC/MS. Secondary stock solutions were prepared using negative urine calibrator matrix and the concentrations were verified gravimetrically using balances calibrated with NIST traceable weights. The secondary stock solutions were then spiked into the negative urine calibrator matrix to the desired concentration and verified by GC/MS to be within $\pm 10\%$ of the target concentration.

Summary:

The information provided in this pre-market notification demonstrates that the LZI Multiple Analyte Sets A, B, and C Drugs of Abuse Calibrators and Controls are substantially equivalent to the legally marketed predicate device for its general intended use. Substantial equivalence was demonstrated through comparison of intended use and physical properties to the commercially available predicate device as confirmed by gas chromatography/mass spectrometry (GC/MS), an independent analytical method. The information supplied in this pre-market notification provides reasonable assurance that the The LZI Multiple Analyte Sets A, B, and C Drugs of Abuse Calibrators and Controls are safe and effective for its stated intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G009
Silver Spring, MD 20993-0002

LIN-ZHI INTERNATIONAL, INC.
BERNICE LIN, PH.D.
VP OPERATIONS
670 ALMANOR AVE
SUNNYVALE CA 94085

February 21, 2014

Re: K133710

Trade/Device Name: Lin-Zhi Multiple Analyte Set A, B and C Urine Drugs of Abuse
Calibrators
Lin-Zhi Multiple Analyte Set A, B and C Urine Drugs of Abuse
Controls

Regulation Number: 21 CFR 862.3200

Regulation Name: Clinical toxicology calibrator

Regulatory Class: II

Product Code: DKB, DIF

Dated: December 2, 2013

Received: December 4, 2013

Dear Dr. Lin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the

electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRI's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Ruth A. Chesler -S

for
Courtney H. Elias, Ph.D.
Director
Division of Chemistry and Toxicology Devices
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
k133710

Device Name

Lin-Zhi Multiple Analyte Set A, B, and C Urine Drugs of Abuse Calibrators; Lin-Zhi Multiple Analyte Set A, B, and C Urine Drugs of Abuse Controls

Indications for Use (Describe)

The Lin-Zhi Multiple Analyte Set A, B, and C Urine Drugs of Abuse Calibrators are intended for in vitro diagnostic use for the calibration of assays for the analytes currently listed in the package insert: Benzoylecgonine, Methamphetamine, Methadone, 6-acetylmorphine, MDMA, Morphine, Oxazepam, and Secobarbital. The calibrators are designed for prescription use with homogeneous enzyme immunoassays on automated clinical chemistry analyzers.

The Lin-Zhi Multiple Analyte Set A, B, and C Urine Drugs of Abuse Controls are intended for in vitro diagnostic use to monitor the performance of assays for the analytes currently listed in the package insert: Benzoylecgonine, Methamphetamine, Methadone, 6-acetylmorphine, MDMA, Morphine, Oxazepam, and Secobarbital. The controls are designed for prescription use with homogeneous enzyme immunoassays on automated clinical chemistry analyzers.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

Avis T. Danishefsky -S

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."