

APR 24 2000



K000432

## 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

---

**MANUFACTURER:** DePuy Orthopaedics, Inc.  
700 Orthopaedic Drive  
Warsaw, Indiana 46581-0988

**510(K) CONTACT:** Arlene C. Saull, RAC  
Sr. Submissions Associate  
DePuy Orthopaedics, Inc.  
700 Orthopaedic Drive  
Warsaw, IN 46581-0988

**TRADE NAME:** Response 2000 Cemented Hip Stem

**COMMON NAME:** Femoral Hip Prosthesis

**CLASSIFICATION:** Class II, per 21 CFR, 888.3350: Hip joint metal/polymer semi-constrained cemented prosthesis.

**DEVICE PRODUCT CODE:** 87 JDI Prosthesis, Hip, Semi-Constrained, Metal/Polymer, Cemented.

**SUBSTANTIALLY EQUIVALENT DEVICES:** Paramount Cemented Hip Stem

**INTENDED USE AND DEVICE DESCRIPTION:**

The subject Response 2000 Cemented Hip Stem is intended for cemented use to replace the femoral portion of the hip joint. The subject device is a smooth-stemmed prosthesis which is very similar to the Paramount Hip Stem design, which has been previously cleared by FDA for cemented use.

**BASIS OF SUBSTANTIAL EQUIVALENCE:**

The subject device is similar to the predicate device in that it uses the same material, has the same stem lengths, same intended use, same indications, same distal tip geometry and uses the same taper neck design.

---

7S&ESummary

000006

2082



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

APR 24 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Arlene C. Saull, RAC  
Senior Submissions Associate  
DePuy Orthopaedics, Inc.  
P.O. Box 988  
700 Orthopaedic Drive  
Warsaw, Indiana 46581-0988

Re: K000432  
Trade Name: Response 2000 Cemented Hip Stem  
Regulatory Class: II  
Product Code : LZO and JDI  
Dated: February 8, 2000  
Received: February 9, 2000

Dear Ms. Saull:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Arlene C. Saull, RAC

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

  
Celia M. Witten, Ph.D., M.D.  
Director  
Division of General and  
Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS

510(k) Number (if known) K000432

Device Name: Response 2000 Cemented Hip Stem

### Indications:

*(The package insert includes intended use/indications for several types of femoral components. The below text, which is extrapolated from the package insert (see Exhibit D), contains only the intended use/indications pertinent to the Response 2000 Cemented Hip Stem.)*

Total hip arthroplasty is intended to provide increased patient mobility and reduce pain by replacing the damaged hip joint articulation in patients where there is evidence of sufficient sound bone to seat and support the components. Total hip replacement is indicated in the following conditions:

1. A severely painful and/or disabled joint from osteoarthritis, traumatic arthritis, rheumatoid arthritis, or congenital hip dysplasia.
2. Avascular necrosis of the femoral head.
3. Acute traumatic fracture of the femoral head or neck.
4. Failed previous hip surgery including joint reconstruction, internal fixation, arthrodesis, hemiarthroplasty, surface replacement arthroplasty, or total hip replacement.
5. Certain cases of ankylosis.

### Cemented Components:

Femoral stem and acetabular cup total hip components labeled "For cemented use only" are indicated only for use with bone cement.

---

Concurrence of CDRH, Office of Device Evaluation:

Denise R. Lochner.  
(Division Sign-Off)  
Division of General Restorative Devices  
510(k) Number K000432

Prescription Use X OR Over-The Counter Use \_\_\_\_ (Per 21 CFR 801.109)

000004