510(K) SUMMARY

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the K012186 requirements of SMDA and 21 CFR §807.92

DEC 0 7 2001

Submitter's Name:

Master & Frank Enterprise Co., Ltd.

Address:

15F-1, No. 57, Sec. 2, Tun Hea S. Rd. Taipei, Taiwan, R.O.C.

Phone:

886-2-2325-5066

Fax:

886-2-2702-6577

Contact:

Mr. Frank Wu (General Manager)

2. Device Name

Trade Name:

Master & Frank Surgical Gowns(Sterile)

Common Name:

Sterile Surgical Gowns

Classification name:

GOWN, SURGICAL

Classification: 3.

Class II

Predicate Device: 4.

Medline Disposable Surgical Drapes & Gowns (K964142)

5. Device Description: Master & Frank Surgical Gowns (Sterile), is manufactured from non-woven

fabric. This surgical Gown is supplied sterile and for single use only.

6. Intended Use: Master & Frank Surgical Gowns (Sterile) is a single use article of surgical apparel worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.

7. Performance Summary: In terms of Physical specification -- ASTM D1424, ASTM D5034 & NFPA Flammability standards---etc, Biological specification ISO 10993 series & Sterilization Specification ISO 17137 & ISO 11607-1, the device are designed to meet applicable standards...

Conclusions:

The Master & Frank Surgical Gowns (Sterile) have the same intended use and similar technological characteristics as the Medline Disposable Surgical Dyapes & Gowns (K964142). Moreover, bench testing contained in this submission demonstrate that any differences in their technological characteristics do not raise any new questions of safety or effectiveness. Thus, the Master & Frank Surgical Gowns (Sterile) is substantially equivalent to the predicate devices.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

DEC 0 7 2001

Master & Frank Enterprise Company Limited C/O Ms. Jennifer Reich Harvest Consulting, Incorporated 3892 South America West Trail Flagstaff, Arizona 86001

Re: K012186

Trade/Device Name: Master & Frank Surgical Gowns (Sterile)

Regulation Number: 878.4040

Regulation Name: Sterile Surgical Gowns

Regulatory Class: II Product Code: FYA

Dated: November 20, 2001 Received: November 23, 2001

Dear Ms. Reich:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours

Timothy A. Ulatowski

Director

Division of Dental, Infection Control and General Hospital Devices Office of Device Evaluation Center for Devices and Radiological Health

		,	
510 (k) NUMBER (IF KN	10WN): <u>K012</u>	186	``.
DEVICE NAME: Mater &	Frank surgical Gov	vns (Sterile)	
INDICATIONS FOR USE:			
Master & Frank Surgical Go intended to be worn by open protect both the surgical pat microorganisms, body fluid	rating room person tient and the operat	nel during surgical procedur ing room personnel from tra	es to help
(PLEASE DO NOT WRITE NEEDED)	E BELOW THIS L	INE-CONTINUE ON ANO	THER PAGE II
Concurrence of	of CDRH, Office of	f Device Evaluation	
Prescription Use(Per 21 CFR 801.109)	OR	Cver-The-Counter (Optional Format)	V

Page __1_ of _1_